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References:

1. Lansdown ABG. Silver 1: its antibacterial properties and mechanism of action. *J Wound Care*. 2002;11(4):125-130.
2. Lansdown ABG. A review of the use of silver in wound care: facts and fallacies. *Br J Nurs*. 2004;13(suppl):S6-S19.
3. Demling RH, DeSanti L. Part 1: effects of silver on wound management. *Wounds*. 2001;13(suppl A):5-15.
4. Bowler PG, Jones SA, Walker M, Parsons D. Microbicidal properties of a silver-containing Hydrofiber dressing against a variety of burn wound pathogens. *J Burn Care Rehabil*. 2004;25(2):192-196.
5. Driver VR. Silver dressings in clinical practice. *Ostomy Wound Manage*. 2004;50(suppl 9A):11S-15S.
6. Bowler P. Progression toward healing: wound infection and the role of an advanced silver-containing Hydrofiber dressing. *Ostomy Wound Manage*. 2003;49(suppl):2-5.
7. Percival SL, Bowler PG, Russell D. Bacterial resistance to silver in wound care. *J Hosp Infect*. 2005;60:1-7.
8. Tosti A, Vincenzi C, Bardazzi F, Mariani R. Allergic contact dermatitis due to povidone-iodine. *Contact Dermatitis*. 1990;23(3):197-198.
9. Hansson C. The effects of cadexomer iodine paste in the treatment of venous leg ulcers compared with hydrocolloid dressing and paraffin gauze dressing. *Int J Dermatol*. 1998;37(5):390-396.
10. Klein DG, Fritsch DE, Amin SG. Wound infection following trauma and burn injuries. *Crit Care Nurse Clin N Am*. 1995;7(4):627-642.
11. Monafu WW, West MA. Current treatment recommendations for topical burn therapy. *Drugs*. 1990;40(3):364-373.
12. Cutting KF. Wound healing, bacteria and topical therapies. *EWMA Journal*. 2003;3(1):17-19.
13. Foster L, Moore P. The application of a cellulose-based fibre dressing in surgical wounds. *J Wound Care*. 1997;6(10):469-473.
14. Walker M, Hobot JA, Newman GR, Bowler PG. Scanning electron microscopic examination of bacterial immobilisation in a carboxymethyl cellulose (AQUACEL®) and alginate dressings. *Biomaterials*. 2003;24:883-890.
15. Bowler PG, Jones SA, Davies BJ, Coyle E. Infection control properties of some wound dressings. *J Wound Care*. 1999;8(10):499-502.
16. Walker M, Cochrane CA. Protease sequestration studies: a comparison between AQUACEL and Promogran in their ability to sequester proteolytic enzymes. WHRI 2494 WA139. May 27, 2003. Data on file, ConvaTec.
17. Jones SA, Bowler PG, Walker M, Parsons D. Controlling wound bioburden with a novel silver-containing Hydrofiber dressing. *Wound Rep Reg*. 2004;12:288-294.
18. Lydon MJ. The development of AQUACEL® Hydrofibre® dressing. In: Krieg T, Harding KG, eds. *AQUACEL® Hydrofibre® dressing: the next step in wound dressing technology*. Proceedings of a satellite symposium at the 6th Congress of the European Academy of Dermatology & Venereology, Dublin, Ireland, 11-15, September 1997. London, UK: Churchill Communications Europe Ltd; 1998:1-3.
19. Robinson BJ. The use of a hydrofibre dressing in wound management. *J Wound Care*. 2000;9(1):32-34.
20. Innes ME, Umraw N, Fish JS, Gomez M, Cartotto RC. The use of silver coated dressings on donor site wounds: a prospective, controlled matched pair study. *Burns*. 2001;27:621-627.
21. Walker M, Cochrane CA, Bowler PG, Parsons D, Bradshaw P. Silver deposition and tissue staining associated with wound dressings containing silver. *Ostomy Wound Manage*. 2006;52(1):42-50.
22. Percival SL, Bowler P, Dolman J. The efficacy of silver dressings against biofilm bacteria. Poster presented at: EWMA Conference; September 15-17, 2005; Stuttgart, Germany. Data on file, ConvaTec.
23. Parsons D, Bowler PG, Myles V, Jones S. Silver antimicrobial dressings in wound management: a comparison of antibacterial, physical, and chemical characteristics. *Wounds*. 2005;17(8):222-232.
24. Jones S. Antimicrobial efficacy of silver-containing dressings competitors report. WHRI2668 MA085. October 20, 2004. Data on file, ConvaTec.
25. Bowler PG. Antimicrobial efficacy of AQUACEL Ag against MRSA. WHRI2474 MA078. January 16, 2003. Data on file, ConvaTec.
26. Jones S, Bowler PG, Walker M. Antimicrobial activity of silver-containing dressings is influenced by dressing conformability with a wound surface. *Wounds*. 2005;17(9):263-270.
27. Vanscheidt W, Lazareth I, Routkovsky-Norval C. Safety evaluation of a new ionic silver dressing in the management of chronic ulcers. *Wounds*. 2003;15(11):371-378.
28. Cains J. A randomised comparative study to assess the performance of AQUACEL Ag in the management of non-ischaemic diabetic foot ulcers. CW-0142-02-U318. June 14, 2005. Data on file, ConvaTec.
29. A phase III, prospective, comparative evaluation of AQUACEL Ag dressing and silver sulfadiazine in the management of partial thickness burns. CW-0142-02-A075. June 9, 2005. Data on file, ConvaTec.
30. Caruso DM, Foster KN, Hermans MHE, Rick C. Aquacel Ag in the management of partial-thickness burns: results of a clinical trial. *J Burn Care Rehabil*. 2004;25:89-97.
31. Lewandowicz G. AQUACEL Ag in the management of open surgical and traumatic wounds. CW-0142-02-U319. June 13, 2005. Data on file, ConvaTec.
32. Cutting KF, White RJ. Maceration of the skin and wound bed 1: its nature and causes. *J Wound Care*. 2002;11(7):275-278.

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