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## Indications for use:

The product is indicated for subcutaneous infusion of medication administered by an external pump.

### PRODUCT DESCRIPTION

The product (drawing 1a, 1b, 1c, page 6) consists of a cannula housing (H) with a 90-degree soft cannula (I) and a tubing (J) with a standard pump connector (luer-lock) (L) combined with an integrated insertion device (A). The product is delivered ready for use as a sterile unit.

### CONTRAINDICATIONS

The product is neither intended nor indicated for intravenous (I.V.) infusion of medication, including blood and blood products.

### CLINICAL BENEFITS

- Continuous subcutaneous infusion of medication to treat or relief disease symptoms per intended use.
- Less invasive administration route associated with fewer adverse events and more locations to infuse medication compared to the intravascular route.
- Less injections (insertions) compared to injection therapy.

### WARNINGS

- Carefully examine the instructions for use before inserting the product. Failure to follow instructions may result in pain or injury.
- Do not use the product if the box and/or blister pack are already opened or damaged. Make sure that the safeguard (B) on top of the insertion device (A) has not been removed. If this occurs please use a new product.
- Do not reposition the product as this may compromise the adhesion of the adhesive tape (E).
- If blood is seen in the cannula housing (H) and/or tubing (J), change to a new product and a new site.
- If the adhesive tape (E) becomes loose, change to a new product and a new site.
- The product integrity may be compromised if exposed to solvents containing substances, e.g. disinfectants, perfumes and deodorants and/or cosmetics and skincare products.
- Do not reuse the infusion set. Reuse of the infusion set may cause infection, site irritation, or damage to the cannula/needle. A damaged cannula/needle may lead to inaccurate medication delivery.
- Never attempt to repair a clogged tubing (J) while the product is connected to the body.

An uncontrolled quantity of medication may accidentally be infused. If clogged, change to at new product and a new site.

- Choice of product and insertion site is closely linked and should always be done by a healthcare professional.
- Never point the insertion device towards any body part where insertion is not desired.

### PRECAUTIONS

- When using the product for the first time, do so in the presence of a healthcare professional.
- The product is a single use item which must be disposed of after use. Do not clean or re-sterilize.
- Wash hands with soap and water when temporarily disconnecting the product. Seal the cannula housing (H) and site connector (K) with the covers (B and M) provided. Consult with a healthcare professional on how to compensate for missed medication when disconnected, and for how long the product can be disconnected from the pump.
- Do not leave air in the tubing (J) when priming. Make sure the medication is kept and handled according to the instructions for the medication.
- Avoid contact with a used soft cannula (I).
- Do not change to a new or unknown product without consulting a healthcare professional for correct handling.

- Change the product in consultation with a healthcare professional after maximum 72 hours of use.
- If using a new insertion device and cannula housing (1b), follow illustrations 1 to 10 in this instruction. Prime and connect the tubing (J) separately referring to illustrations G, H and I.

### STORAGE & DISPOSAL

- Store infusion sets in a cool, dry location at room temperature (15-25°C). Do not store infusion sets in direct sunlight, or in high humidity.
- Dispose of the insertion device in a proper sharps container according to local laws.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

### ADDITIONAL INFORMATION FOR THE PATIENT

- Use an insertion site as recommended by your healthcare professional. Do not use an insertion site which is infected, swollen, red, irritated, bruised or is bleeding. Do not insert in an area with lumps, hard tissue or bumps. Be careful not to insert into bone or muscle as this may cause pain, damage the product and affect medication absorption. Choose an area that does not have constrictive clothing or other objects that can pull on the infusion set.

- Please refer to your healthcare professional for insertion site rotation and general recommendations.
- Check the insertion site frequently to ensure that the soft cannula (I) remains firmly in place. If it is not in place, change to a new product and a new site. The soft cannula (I) must always be completely inserted to ensure delivery of the full amount of medication. If the insertion site becomes red or inflamed, change to a new product and a new site.
- The product is available with various soft cannula (I) and tubing (J) lengths. Follow the recommendations of your healthcare professional regarding appropriate soft cannula (I) and tubing (J) length.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.


### ADDITIONAL INFORMATION FOR THE HEALTHCARE PROFESSIONAL

- Please consider treatment and patient specific factors when choosing an insertion site.
- The choice of soft cannula length should be based on treatment and patient specific factors such as physiology and activity level. A too long soft cannula may result in insertion pain and/or insertion into a bone or muscle and a too short soft cannula in leakage or irritation at the insertion site. In both cases the infusion must be terminated and treatment should not be resumed until an infusion set with another soft cannula length is inserted at a new insertion site.
- Please inform the patient about the correct choice of insertion site in relation to the choice of soft cannula length.


## Preparation and insertion

Go through the instructions for use carefully before inserting the product, as failure to follow instructions may result in inaccurate medication delivery, infection and/or site irritation. For information about pump therapy, including connection and priming procedures, possible errors, and potential risks and their causes, consult your pump manufacturer's instructions for use before using the product. When using the product for the first time, do so in the presence of a healthcare professional.


### THE PRODUCT (1a)

### INSERTION DEVICE & CANNULA HOUSING (1b)

### TUBING (1c)

### 1 CONTENTS

[A] Insertion device

[B] Safeguard (Optional: transparent cover for cannula housing)

[C] Activation button (for insertion of the infusion set)

[D] Indication for tubing direction

[E] Adhesive tape

[F] Paper

[G] Insertion site window

[H] Cannula housing

[I] Soft cannula

[J] Tubing

[K] Site Connector

[L] Pump connector (luer-lock)

[M] White cover for site connector

- 2 Wash your hands with soap and water.
- 3 Use an insertion site as recommended by your healthcare professional. The choice of site depends on treatment and patient specific factors. Check that the insertion site is free of skin irritation such as redness, scar tissue, etc.
- 4 Clean the area of insertion with a disinfectant as recommended by your healthcare professional. Be sure the area is air dried before inserting the product. If needed, remove hair around the insertion site to ensure proper attachment of the adhesive tape (E) to the skin.
- 5 Open the blister pack and take out the contents.
- 6 Pull gently to remove the paper (F) from the adhesive tape (E). Be careful not to touch the adhesive tape (E).
- 7 Remove the safeguard (B) from the insertion device (A) by gently squeezing the sides of the safeguard (B) while pulling it straight out.
- 8 Before placing the insertion device (A) against the skin, note the arrows (D) on the insertion device (A) showing where the tubing (J) will be connected after insertion.
- 9 Gently stretch the skin and ensure that the surface is flat and taut. Position the insertion device (A) on the skin and press the red activation button (C) completely down to insert the infusion set.

The insertion is complete when the insertion device (A) is released from the cannula housing (H) with adhesive tape (E). The product is now in safe mode with the activation button (C) pressed down. The product is now inserted.

- 10 Press the adhesive tape (E) onto the skin. Change to a new product and a new site, if the adhesive tape (E) becomes loose.
- 11 Remove the paper from the tubing (J).
- 12 Connect the tubing (J) to the pump or syringe driver by turning the pump connector (L) clockwise.
- 13 Remove the white cover (M) from the site connector (K).
- 14 Prime\* the tubing (J) according to the instructions provided by the pump manufacturer. Be sure there is no air in the tubing (J).
- 15 Place a finger on the cannula housing (H) while pushing the site connector (K) straight in until you hear a “click”.
- 16 Frequently check the insertion site to ensure that the soft cannula (I) remains in place and that the site is not irritated. Change to a new product and a new site if necessary. Remember to store the covers (B and M) for later use. Dispose of trash according to local recommendations.

### DISCONNECTING THE TUBING

- A Wash your hands with soap and water.
- B Place a finger on the cannula housing (H) and gently squeeze the sides of the site connector (K). Pull the site connector (K) straight out from the cannula housing (H).
- C Place the white cover (M) on the site connector (K).
- D Insert the transparent cover (B) into the cannula housing (H) until you hear a “click.”

### RECONNECTING THE TUBING

- E Wash your hands with soap and water.
- F Remove the transparent cover (B) by placing a finger on the cannula housing (H) and gently squeezing the sides of the cover. Pull the cover straight out of the cannula housing (H).
- G Remove the white cover (M) from the site connector (K).
- H Make sure there is medication flow and no air in the tubing (J).
- I Place a finger on the cannula housing (H) while pushing the site connector (K) straight in until you hear a “click”.

### REMOVAL OF THE INFUSION SET

- J Wash your hands with soap and water.
- K Place a finger on the cannula housing (H) and gently squeeze the sides of the site connector (K). Pull the site connector (K) straight out from the cannula housing (H).
- L Remove the infusion set by carefully loosening the adhesive tape (E) around the cannula housing (H) before pulling the soft cannula (I) out of the skin.
- M Disconnect the tubing (J) from the pump or syringe driver by turning the pump connector (L) counter clock-wise. Dispose of the infusion set according to local recommendations. For available sharps containers, please visit your local pharmacy.

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