

Product Description

Exudate Management – For moderately to heavily exuding wounds

AQUACEL® Ag Hydrofiber® dressing is Sodium Carboxymethylcellulose with 1.2% ionic silver. It offers the unique gelling properties of Hydrofiber® technology with the power of ionic silver. The dressing fibres of AQUACEL® Ag gel on contact with wound fluid by hydrophilic action. The fibres swell as they lock bacterial exudate away from the wound by vertical wicking¹, creating a large fluid absorption capacity².

AQUACEL® Ag dressing can kill a broad spectrum of wound pathogens in the dressing that can cause infection including MRSA, VRE and Pseudomonas aeruginosa and aids in reducing the wound bioburden³.

Bacteria and other harmful wound exudate components are locked within the gelled AQUACEL® Ag dressing fibres¹, where the ionic silver contained within the dressing kills them³.

Application*

Before applying, carefully cleanse the wound with saline and dry surrounding skin. Apply directly onto the wound leaving an adequate overlap onto the surrounding skin of at least 1cm. When using AQUACEL® Ag dressing ribbon in cavity wounds, loosely pack wound to about 80% capacity leaving at least 2.5cm outside the wound for easy retrieval. The AQUACEL® Ag dressing should then be covered with an appropriate secondary wound care dressing.

*Please refer to pack insert for full instructions prior to use.

Specifications

Size	Dressings Per Box	Product Code
5 cm x 5 cm	10	403706
10 cm x 10 cm	10	403708
15 cm x 15 cm	5	403710
20 cm x 30 cm	5	403711
2 cm x 45 cm rope	5	403771

Surgical Size Dressings

Size	Dressings Per Box	Product Code
4 cm x 10 cm	10	403739
4 cm x 20 cm	10	403740
4 cm x 30 cm	10	403741

ARTG Number: 156976

References

1. Bowler PG, Jones SA, Davies BJ, Coyle E. Infection control properties of some wound dressings. J Wound Care. 1999; 8: 499-502.
2. Performance Analysis of AQUACEL® vs Sorbsan® and Sorbsan® plus AH28. Data on file, ConvaTec.
3. Jones SA, Bowler P, Walker M, Parsons D. Controlling wound bioburden with a novel silver containing Hydrofiber® dressing. Wound Repair and Regeneration. May-June 2004, 288-294.
4. Waring M, Parsons D, Clay C. A new wound care. Data on file, ConvaTec.
5. Lyndon MJ. The development of AQUACEL® Hydrofiber® dressing. The development of AQUACEL® Hydrofiber® dressing. AQUACEL® Hydrofiber® dressing: the next step in wound dressing.
6. Harding KG, Price P, Robinson B, Thomas S, Holman D. Cost and dressing evaluation of Hydrofiber® and alginate dressings in the management of community-based patients with chronic leg ulceration. Wounds. 2001; 13:229-236.
7. B.J. Robinson. The use of a hydrofibre dressing in wound management. J Wound Care. 2000; 9 (1): 32-34.
8. Walker M. Demonstration of the role of intimate contact in the killing capacity of silver, poster presented at SAWC, April 2005.



Clinical Benefits

- Supports moist wound healing environment
- Contains ionic silver which can kill a broad spectrum of wound pathogens that can cause infection – including MRSA, VRE and Pseudomonas aeruginosa³
- Absorbs and retains higher levels of exudate than most alginates^{4,5,6} and gauze^{4,5}
- Reduces maceration by locking exudate into its fibres and away from the skin⁷
- Cohesive gel allows easy, atraumatic removal
- Soft and conformable allowing intimate contact with the wound surface. This reduces “dead space” where bacteria may reside⁸

*AQUACEL® Ag dressing is of the same construction as AQUACEL® dressing with the addition of ionic silver

Indications

AQUACEL® Ag dressing is indicated as a primary dressing for moderately to heavily exuding wounds where there is an infection or an increased risk of infection.

- Chronic wounds: leg ulcers, pressure ulcers, diabetic ulcers and fungating lesions
- Acute wounds: post-surgical wounds, abrasions, lacerations and partial thickness burns



Call your local ConvaTec Territory Manager or ConvaTec Support Centre for further information:

Australia: 1800 335 276 or New Zealand: 0800 441 763

www.convatec.com

ConvaTec (Australia) Pty Limited. ABN 70 131 232 570. Unipark Monash, Building 2, Ground Floor, 195 Wellington Road, Clayton VIC 3170 Australia. PO Box 63, Mulgrave, Vic 3170. Phone: (03) 9239 2700 Facsimile: (03) 9239 2743. ConvaTec (New Zealand) Limited. AK2135265 Level 4, 369 Queen Street, Auckland 1010 New Zealand. PO Box 62663, Kalmia Street, Auckland 1544 New Zealand. Phone: (09) 306 8833 Facsimile: (09) 306 8831.

®/™Indicates trademark of ConvaTec Inc. © 2009 ConvaTec Inc. AP-006687-AU. March 2009.

Our world is what we make of it™