

## Product Description

### Exudate Management – As a primary or secondary cover dressing for moderately to heavily exuding wounds

CombiDERM ACD® dressing is an absorbent dressing that is made up of three main parts: an absorbent hydrocolloid adhesive, a central absorbent island pad and gel forming polyacrylate granules. The super absorbent granules act to absorb, gel and hold exudate away from the skin. This helps to create a moist wound healing environment.

CombiDERM ACD® dressing retains exudate even under pressure<sup>1</sup> helping prevent maceration and allowing longer wear time<sup>1</sup>. The outer polyurethane film layer is waterproof and bacterial resistant<sup>2</sup>.

## Application\*

Carefully cleanse the wound with saline, and dry surrounding skin. Choose a dressing to ensure that the central pad size is larger than the wound area. Remove the release paper and place the pad directly over the wound. The absorbent portion of the pad will swell and soften as it absorbs wound exudate.

To remove the dressing, press down gently on the skin, carefully lift the blue tab and continue until all edges are free.

CombiDERM ACD® dressing may be used as a primary dressing or as a secondary dressing in combination with other appropriate dressings such as, but not limited to, AQUACEL®, AQUACEL® Ag, KALTOSTAT® or DuoDERM® Hydroactive® Gel.

\*Please refer to pack insert for full instructions prior to use.

## Specifications

Size	Dressings Per Box	Product Code
10 cm x 10 cm	10	651031
14 cm x 14 cm	10	187725
15 cm x 25 cm	5	651029

### Sacral Size Dressings

Size	Dressings Per Box	Product Code
15 cm x 18 cm	5	651027

### RPBS Codes

Size	RPBS Code	Product Code
10 cm x 10 cm	4692Y	651031
15 cm x 18 cm sacral	4693B	651027

ARTG Number: 158691

## References

1. Chen J. CombiDERM® Absorbent Cover Dressing. Data on file, ConvaTec. (1998).
2. Ameen H, Moore K, Lawrence J.C, Harding K.G. Investigating the bacterial barrier properties of four contemporary wound dressings. J Wound Care. 2000 9 (8) 385-387.
3. Hudak J, Constantine B, Buglino D, Eberhardt D. An in vitro evaluation of the fluid absorption capabilities under compression of a new non-adhesive composite dressing. Data on file, ConvaTec. (1998).



## Clinical Benefits

- Supports moist wound healing environment
- Outer film layer provides waterproof and bacterial barrier\*<sup>2</sup>
- Poly-acrylate granules absorb and lock in exudate reducing the risk of maceration<sup>1</sup>
- Non-adherent wound contact layer allowing for atraumatic removal<sup>1</sup>
- Effective under compression when tested against comparative products\*<sup>3</sup>

\* in vitro

## Indications

CombiDERM ACD® dressing is indicated for the management of moderately to heavily exuding wounds.

- Chronic wounds: leg ulcers, diabetic ulcers, pressure ulcers
- Acute wounds: abrasions, lacerations, surgical wounds

Call your local ConvaTec Territory Manager or ConvaTec Support Centre for further information:

**Australia: 1800 335 276 or New Zealand: 0800 441 763**

**www.convatec.com**

ConvaTec (Australia) Pty Limited. ABN 70 131 232 570. Unipark Monash, Building 2, Ground Floor, 195 Wellington Road, Clayton VIC 3170 Australia. PO Box 63, Mulgrave, Vic 3170. Phone: (03) 9239 2700 Facsimile: (03) 9239 2743. ConvaTec (New Zealand) Limited. AK2135265 Level 4, 369 Queen Street, Auckland 1010 New Zealand. PO Box 62663, Kalmia Street, Auckland 1544 New Zealand. Phone: (09) 306 8833 Facsimile: (09) 306 8831.

®/™Indicates trademark of ConvaTec Inc. © 2009 ConvaTec Inc. AP-006690-AU. March 2009.

Our world is what we make of it™

