
STERILE/ESTERIL  
**DuoDERM**  
**CGF** **BORDER DRESSING**  
CONTROL GEL FORMULA BORDER DRESSING

STERILE/ESTERIL  
**DuoDERM**  
**CGF** **BORDER DRESSING**  
CON FORMULA DE GEL CONTROLADA  
APOSITO CON REBORDE

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
**III** ConvaTec

## PRODUCT DESCRIPTION

DuoDERM® Control Gel Formula Border Dressings interact with wound exudate producing a soft mass that enables removal of the dressing with little or no damage to newly formed tissues. They help isolate the wound against bacterial and other external contamination. These dressings include an adhesive border so that additional taping (picture framing) is not required.

## INDICATIONS

- Dermal Ulcers, including full thickness wounds
  - Pressure Sore Management, Stage I through Stage IV
  - Leg Ulcer Management
- Superficial Wounds — e.g. Minor Abrasions
- Burns — Second Degree
- Donor Sites

## CONTRAINDICATION

- Use on individuals with a known sensitivity to the dressing or its components.

## PRECAUTIONS AND OBSERVATIONS

DuoDERM® CGF® Border Dressings are not recommended for use on third degree burns.

### When Used on Dermal Ulcers:

1. Initial use of this product should be under the direction of a health professional.
2. DuoDERM CGF Border Dressings only provide local management of the wound site. In pressure sore care, other aspects such as repositioning of the patient and nutritional support should not be neglected. In leg ulcer care, lack of adequate rest in patients with vascular (arterial or venous) insufficiency can increase the amount of local edema and hinder potential healing.
3. Increased Wound Size: Deeper tissue damage may have already occurred under an apparent superficial dermal ulcer. When using any occlusive dressing in the presence of necrotic material, the wound may increase in size and depth during the initial phase of management as the necrotic debris is cleaned away. Leg ulcers resulting from vasculitis may rapidly deteriorate during exacerbation of the underlying disorder.

**Odor:** Wounds, particularly those that are large or necrotic, are often accompanied by a disagreeable odor; however, this is not necessarily indicative of infection. The odor should disappear when the wound is cleansed (see infection).

**Infection:** If signs of clinical infection should develop, such as: uncharacteristic odor or change in the color of the exudate, fever or cellulitis (tenderness and erythema in the area of the wound), a bacterial culture of the wound site should be taken. If clinical signs of infection are present, appropriate medical treatment should be initiated. DuoDERM dressings may be continued during the treatment at the discretion of the clinician.

**Granulation:** Excessive granulation tissue may develop in some wounds when using "occlusive" dressings.

## INSTRUCTIONS FOR USE

### Preparing and Cleansing the Wound Site

DuoDERM CGF Border dressings are sterile and should be handled appropriately.

Choose a dressing that allows the DuoDERM CGF mass to extend beyond the wound margin by at least 1/4", covering all unhealthy tissue. Cleanse the wound according to hospital practice. Irrigate with saline and dry the surrounding skin to ensure it is free of any greasy substance. Use of the dressings help facilitate the liquefaction and removal of dead tissue; however, eschar that is particularly thick or fused to the wound margins should be removed prior to application of the dressing.

**NOTE:** DuoDERM® Paste or DuoDERM® Hydroactive® Granules may be used with this dressing. See the instructions for use with each product.

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


## APPLYING THE DRESSING


1. Remove top backing paper *only*.


Note: The triangle shaped dressing can be applied in several directions depending on the location of the ulcer.


In sacral ulcers, the dressing may be folded in half lengthwise to make it easy to apply in the sacral folds.


2. Apply the dressing over the wound. Smooth into place, especially at edges of the center adhesive.

3. Fold back the border and remove the release papers from the border adhesive; press the borders into place.

## REMOVING THE DRESSING

1. Press down on the skin and carefully lift an edge of the dressing. Continue around until all edges are free. Repeat cleansing procedure. It is unnecessary to remove all residual dressing material from the surrounding skin.

Leave the dressing in place (not more than 7 days) unless it is uncomfortable, leaking, or there are clinical signs of infection.

If further information or guidance is needed, please contact ConvaTec Professional Services.

ConvaTec  
A Division of E.R. Squibb & Sons, L.L.C.  
Princeton, NJ 08543  
1-800-422-8811

[www.convatec.com](http://www.convatec.com)

Made in USA

© 2004 E.R. Squibb & Sons, L.L.C.  
®/™ is a registered trademark of E.R. Squibb & Sons, L.L.C.

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