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## Referencias

1. Malbrain ML, De laet IE: Intra-abdominal hypertension: evolving concepts. *Clin Chest Med* 2009, 30(1):45-70, viii.
2. Malbrain ML, Cheatham ML, Kirkpatrick A, Sugrue M, Parr M, De Waele J, Balogh Z, Leppaniemi A, Olvera C, Ivatury R et al: Results from the International Conference of Experts on Intra-abdominal Hypertension and Abdominal Compartment Syndrome. I. Definitions. *Intensive care medicine* 2006, 32(11):1722-1732.
3. Fietsam R, Jr., Villalba M, Glover JL, Clark K: Intra-abdominal compartment syndrome as a complication of ruptured abdominal aortic aneurysm repair. *Am Surg* 1989, 55(6):396-402.
4. Malbrain ML, De laet I, Cheatham M: Consensus conference definitions and recommendations on intra-abdominal hypertension (IAH) and the abdominal compartment syndrome (ACS)—the long road to the final publications, how did we get there? *Acta Clin Belg Suppl* 2007, 62(1):44-59.
5. Malbrain ML, Cheatham ML, Kirkpatrick A, Sugrue M, De Waele J, Ivatury R: Abdominal compartment syndrome: it's time to pay attention! *Intensive care medicine* 2006, 32(11):1912-1914.
6. Malbrain ML: You don't have any excuse, just start measuring abdominal pressure and act upon it! *Minerva Anestesiologica* 2008, 74(2):1-2.
7. Kirkpatrick AW, Brennehan FD, McLean RF, Rapanos T, Boulanger BR: Is clinical examination an accurate indicator of raised intra-abdominal pressure in critically injured patients? *CanuSurg* 2000, 43(3):207-211.
8. Sugrue M, Bauman A, Jones F, Bishop G, Flabouris A, Parr M, Stewart A, Hillman K, Deane SA: Clinical examination is an inaccurate predictor of intraabdominal pressure. *World JSurg* 2002, 26(12):1428-1431.
9. Malbrain ML, De laet I, Van Regenmortel N, Schoonheydt K, Dits H: Can the abdominal perimeter be used as an accurate estimation of intra-abdominal pressure? *Crit Care Med* 2009, 37(1):316-319.
10. Malbrain ML: Different techniques to measure intra-abdominal pressure (IAP): time for a critical re-appraisal. *Intensive care medicine* 2004, 30(3):357-371.
11. Malbrain M, Jones F: Intra-abdominal pressure measurement techniques. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 19-68.
12. De Potter TJ, Dits H, Malbrain ML: Intra- and interobserver variability during in vitro validation of two novel methods for intra-abdominal pressure monitoring. *Intensive care medicine* 2005, 31(5):747-751.
13. Malbrain ML: The assumed problem of air bubbles in the tubing during intra-abdominal pressure measurement - author reply. *Intensive care medicine* 2004, 30(8):1693.
14. Cheatham ML, Malbrain ML, Kirkpatrick A, Sugrue M, Parr M, De Waele J, Balogh Z, Leppaniemi A, Olvera C, Ivatury R et al: Results from the International Conference of Experts on Intra-abdominal Hypertension and Abdominal Compartment Syndrome. II. Recommendations. *Intensive care medicine* 2007, 33(6):951-962.
15. Saggi B, Ivatury R, Sugerman HJ: Surgical critical care issues: Abdominal compartment syndrome. In: *Surgical Treatment Evidence-Based and Problem-Oriented*. Edited by Holzheimer RG, Mannick JA. München: W. Zuckschwerdt Verlag München; 2001.
16. Davis PJ, Koottayi S, Taylor A, Butt WW: Comparison of indirect methods of measuring intra-abdominal pressure in children. *Intensive care medicine* 2005, 31(3):471-475.
17. Suominen PK, Pakarinen MP, Rautiainen P, Mattila J, Sairanen H: Comparison of direct and intravesical measurement of intraabdominal pressure in children. *J Pediatr Surg* 2006, 41(8):1381-1385.
18. Sanchez NC, Tenofsky PL, Dort JM, Shen LY, Helmer SD, Smith RS: What is normal intra-abdominal pressure? *Am Surg* 2001, 67(3):243-248.
19. Sugerman H, Windsor A, Bessos M, Wolfe L: Intra-abdominal pressure, sagittal abdominal diameter and obesity comorbidity. *UlnterMed* 1997, 241(1):71-79.
20. Sugerman HJ: Effects of increased intra-abdominal pressure in severe obesity. *SurgClinNorth Am* 2001, 81(5):1063-1075, vi.
21. De Keulenaer BL, De Backer A, Schepens DR, Daelmans R, Wilmer A, Malbrain ML: Abdominal compartment syndrome related to noninvasive ventilation. *Intensive Care Med* 2003, 29(7):1177-1181.
22. De Keulenaer BL, De Waele JJ, Powell B, Malbrain ML: What is normal intra-abdominal pressure and how is it affected by positioning, body mass and positive end-expiratory pressure? *Intensive care medicine* 2009, 35(6):969-976.
23. Hering R, Vorwerk R, Wrigge H, Zinserling J, Schroder S, von Spiegel T, Hoeff A, Putensen C: Prone positioning, systemic hemodynamics, hepatic indocyanine green kinetics, and gastric intramucosal energy balance in patients with acute lung injury. *Intensive Care Med* 2002, 28(1):53-58.
24. Hering R, Wrigge H, Vorwerk R, Brensing KA, Schroder S, Zinserling J, Hoeff A, Spiegel TV, Putensen C: The effects of prone positioning on intraabdominal pressure and cardiovascular and renal function in patients with acute lung injury. *AnesthAnalg* 2001, 92(5):1226-1231.
25. Citerio G, Berra L: Central nervous system. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 144-156.
26. Citerio G, Vascotto E, Villa F, Celotti S, Pesenti A: Induced abdominal compartment syndrome increases intracranial pressure in neurotrauma patients: a prospective study. *Crit Care Med* 2001, 29(7):1466-1471.
27. Deeren D, Dits H, Malbrain MLNG: Correlation between intra-abdominal and intracranial pressure in nontraumatic brain injury. *Intensive care medicine* 2005, 31(1):1577-1581.
28. Bloomfield G, Saggi B, Blocher C, Sugerman H: Physiologic effects of externally applied continuous negative abdominal pressure for intra-abdominal hypertension. *The Journal of trauma* 1999, 46(6):1009-1014; discussion 1014-1006.
29. Bloomfield GL, Dalton JM, Sugerman HJ, Ridings PC, DeMaria EJ, Bullock R: Treatment of increasing intracranial pressure secondary to the acute abdominal compartment syndrome in a patient with combined abdominal and head trauma. *The Journal of trauma* 1995, 39(6):1168-1170.
30. Bloomfield GL, Ridings PC, Blocher CR, Marmarou A, Sugerman HJ: Effects of increased intra-abdominal pressure upon intracranial and cerebral perfusion pressure before and after volume expansion. *The Journal of trauma* 1996, 40(6):936-941; discussion 941-933.
31. Bloomfield GL, Ridings PC, Blocher CR, Marmarou A, Sugerman HJ: A proposed relationship between increased intra-abdominal, intrathoracic, and intracranial pressure. *Crit Care Med* 1997, 25(3):496-503.
32. Saggi BH, Bloomfield GL, Sugerman HJ, Blocher CR, Hull JP, Marmarou AP, Bullock MR: Treatment of intracranial hypertension using nonsurgical abdominal decompression. *The Journal of trauma* 1999, 46(4):646-651.
33. Irgau I, Koyfman Y, Tikellis JI: Elective intraoperative intracranial pressure monitoring during laparoscopic cholecystectomy. *Arch Surg* 1995, 130(9):1011-1013.
34. Joseph DK, Dutton RP, Aarabi B, Scalea TM: Decompressive laparotomy to treat intractable intracranial hypertension after traumatic brain injury. *The Journal of trauma* 2004, 57(4):687-693; discussion 693-685.
35. Josephs LG, Este-McDonald JR, Birkett DH, Hirsch EF: Diagnostic laparoscopy increases intracranial pressure. *The Journal of trauma* 1994, 36(6):815-818; discussion 818-819.
36. Sugerman HJ, DeMaria EJ, Felton WL, III, Nakatsuka M, Sismanis A: Increased intra-abdominal pressure and cardiac filling pressures in obesity-associated pseudotumor cerebri. *Neurology* 1997, 49(2):507-511.
37. Sugerman HJ, Felton IW, III, Sismanis A, Saggi BH, Doty JM, Blocher C, Marmarou A, Makhoul RG: Continuous negative abdominal pressure device to treat pseudotumor cerebri. *IntJObesRelat Metab Disord* 2001, 25(4):486-490.
38. Sugerman H, Windsor A, Bessos M, Kellum J, Reines H, DeMaria E: Effects of surgically induced weight loss on urinary bladder pressure, sagittal abdominal diameter and obesity co-morbidity. *IntJObesRelat Metab Disord* 1998, 22(3):230-235.
39. Kashtan J, Green JF, Parsons EQ, Holcroft JW: Hemodynamic effect of increased abdominal pressure. *The Journal of surgical research* 1981, 30(3):249-255.
40. Ridings PC, Bloomfield GL, Blocher CR, Sugerman HJ: Cardiopulmonary effects of raised intra-abdominal pressure before and after intravascular volume expansion. *The Journal of trauma* 1995, 39(6):1071-1075.
41. Richardson JD, Trinkle JK: Hemodynamic and respiratory alterations with increased intra-abdominal pressure. *The Journal of surgical research* 1976, 20(5):401-404.
42. Malbrain ML, Cheatham ML: Cardiovascular effects and optimal preload markers in intra-abdominal hypertension. In: *Yearbook of Intensive Care and Emergency Medicine*. Edited by Vincent J-L. Berlin: Springer-Verlag; 2004: 519-543.
43. Cheatham M, Malbrain M: Abdominal perfusion pressure. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 69-81.
44. Cheatham M, Malbrain M: Cardiovascular implications of elevated intra-abdominal pressure. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 89-104.
45. Cheatham ML, Block EF, Nelson LD, Safcsak K: Superior predictor of the hemodynamic response to fluid challenge in critically ill patients. *Chest* 1998, 114(4):1226-1227.
46. Cheatham ML, Nelson LD, Chang MC, Safcsak K: Right ventricular end-diastolic volume index as a predictor of preload status in patients on positive end-expiratory pressure. *Crit Care Med* 1998, 26(11):1801-1806.
47. Schachtrupp A, Graf J, Tons C, Hoer J, Fackeldey V, Schumpelick V: Intravascular volume depletion in a 24-hour porcine model of intra-abdominal hypertension. *JTrauma* 2003, 55(4):734-740.
48. Michard F, Alaya S, Zarka V, Bahloul M, Richard C, Teboul JL: Global end-diastolic volume as an indicator of cardiac preload in patients with septic shock. *Chest* 2003, 124(5):1900-1908.
49. Michard F, Teboul JL: Predicting fluid responsiveness in ICU patients: a critical analysis of the evidence. *Chest* 2002, 121(6):2000-2008.
50. Malbrain ML, de Laet I: Functional hemodynamics and increased intra-abdominal pressure: same thresholds for different conditions ...? *Crit Care Med* 2009, 37(2):781-783.
51. Malbrain ML, Reuter DA: Assessing fluid responsiveness with the passive leg raising maneuver in patients with increased intra-abdominal pressure: Be aware that not all blood returns! *Crit Care Med* 2010, 38(9):1912-1915.
52. Simon RJ, Friedlander MH, Ivatury RR, DiRaimo R, Machiedo GW: Hemorrhage lowers the threshold for intra-abdominal hypertension-induced pulmonary dysfunction. *The Journal of trauma* 1997, 42(3):398-403; discussion 404-395.
53. Burchard KW, Ciombor DM, McLeod MK, Slothman GJ, Gann DS: Positive end expiratory pressure with increased intra-abdominal pressure. *SurgGynecolObstet* 1985, 161(4):313-318.
54. Pelosi P, Ravagnan I, Giurati G, Panigada M, Bottino N, Tredici S, Eccher G, Gattinoni L: Positive end-expiratory pressure improves respiratory function in obese but not in normal subjects during anesthesia and paralysis. *Anesthesiology* 1999, 91(5):1221-1231.
55. Sugrue M, D'Amours S: The problems with positive end expiratory pressure (PEEP) in association with abdominal compartment syndrome (ACS). *JTrauma* 2001, 51(2):419-420.
56. Sussman AM, Boyd CR, Williams JS, DiBenedetto RJ: Effect of positive end-expiratory pressure on intra-abdominal pressure. *SouthMedJ* 1991, 84(6):697-700.
57. Cheatham ML, White MW, Sagraves SG, Johnson JL, Block EF: Abdominal perfusion pressure: a superior parameter in the assessment of intra-abdominal hypertension. *The Journal of trauma* 2000, 49(4):621-626; discussion 626-627.
58. Malbrain ML: Abdominal perfusion pressure as a prognostic marker in intra-abdominal hypertension. In: *Yearbook of Intensive Care and Emergency Medicine*. Edited by Vincent JL. Berlin: Springer-Verlag; 2002: 792-814.
59. Mertens zur Borg IR, Verbrugge SJ, Olvera C: Pathophysiology: Respiratory. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 105-118.
60. Mutoh T, Lamm WJ, Embree LJ, Hildebrandt J, Albert RK: Abdominal distension alters regional pleural pressures and chest wall mechanics in pigs in vivo. *JApplPhysiol* 1991, 70(6):2611-2618.
61. Mutoh T, Lamm WJ, Embree LJ, Hildebrandt J, Albert RK: Volume infusion produces abdominal distension, lung compression, and chest wall stiffening in pigs. *J Appl Physiol* 1992, 72(2):575-582.
62. Quintel M, Pelosi P, Caironi P, Meinhardt JP, Luecke T, Herrmann P, Taccone P, Rylander C, Valenza F, Carlesso E et al: An increase of abdominal pressure increases pulmonary edema in oleic acid-induced lung injury. *Am J Respir Crit Care Med* 2004, 169(4):534-541.
63. Pelosi P, Quintel M, Malbrain ML: Effect of intra-abdominal pressure on respiratory mechanics. *Acta Clin Belg Suppl* 2007, 62(1):78-88.
64. Diebel LN, Wilson RF, Dulchavsky SA, Saxe J: Effect of increased intra-abdominal pressure on hepatic arterial, portal venous, and hepatic microcirculatory blood flow. *The Journal of trauma* 1992, 33(2):279-282; discussion 282-273.
65. Wendon J, Biancofiore G, Auzinger G: Intra-abdominal hypertension and the liver. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 138-143.
66. Biancofiore G, Bindi ML, Boldrini A, Consani G, Bisa M, Esposito M, Urbani L, Catalano G, Filippini F, Mosca F: Intraabdominal pressure in liver transplant recipients: incidence and clinical significance. *Transplant Proc* 2004, 36(3):547-549.
67. Biancofiore G, Bindi ML, Romanelli AM, Boldrini A, Consani G, Bisa M, Filippini F, Vagelli A, Mosca F: Intra-abdominal pressure monitoring in liver transplant recipients: a prospective study. *Intensive Care Med* 2003, 29(1):30-36.
68. Michelet P, Roch A, Gannier M, Sainty JM, Auffray JP, Papazian L: Influence of support on intra-abdominal pressure, hepatic kinetics of indocyanine green and extravascular lung water during prone positioning in patients with ARDS: a randomized crossover study. *Critical care (London, England)* 2005, 9(3):R251-257.
69. Schein M: Abdominal compartment syndrome: historical background. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 1-7.
70. Biancofiore G, Bindi ML, Romanelli AM, Bisa M, Boldrini A, Consani G, Filippini F, Mosca F: Postoperative intra-abdominal pressure and renal function after liver transplantation. *ArchSurg* 2003, 138(7):703-706.
71. Sugrue M, Hallal A, D'Amours S: Intra-abdominal pressure hypertension and the kidney. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 119-128.
72. Sugrue M, Buist MD, Hourihan F, Deane S, Bauman A, Hillman K: Prospective study of intra-abdominal hypertension and renal function after laparotomy. *BrJSurg* 1995, 82(2):235-238.
73. Sugrue M, Jones F, Deane SA, Bishop G, Bauman A, Hillman K: Intra-abdominal hypertension is an independent cause of postoperative renal impairment. *ArchSurg* 1999, 134(10):1082-1085.
74. Bradley SE, Mudge GH, Blake WD, Alphonse P: The effect of increased intra-abdominal pressure on the renal excretion of water and electrolytes in normal human subjects and in patients with diabetes insipidus. *Acta Clin Belg* 1955, 10(3):209-223.
75. Harman PK, Kron IL, McLachlan HD, Freedlander AE, Nolan SP: Elevated intra-abdominal pressure and renal function. *AnnSurg* 1982, 196(5):594-597.
76. Ivatury R, Diebel L: Intra-abdominal hypertension and the splanchnic bed. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 129-137.
77. Sugrue M, Jones F, Lee A, Buist MD, Deane S, Bauman A, Hillman K: Intraabdominal pressure and gastric intramucosal pH: is there an association? *World J Surg* 1996, 20(8):988-991.
78. Ivatury RR, Porter JM, Simon RJ, Islam S, John R, Stahl WM: Intra-abdominal hypertension after life-threatening penetrating abdominal trauma: prophylaxis, incidence, and clinical relevance to gastric mucosal pH and abdominal compartment syndrome. *The Journal of trauma* 1998, 44(6):1016-1021; discussion 1021-1013.
79. Balogh Z, McKinley BA, Cocanour CS, Kozar RA, Valdivia A, Sailors RM, Moore FA: Supranormal trauma resuscitation causes more cases of abdominal compartment syndrome. *Arch Surg* 2003, 138(6):637-642; discussion 642-633.
80. Balogh Z, McKinley BA, Cox Jr CS, Allen SJ, Cocanour CS, Kozar RA, Moore EE, Miller IC, Weisbrodt NW, Moore FA: Abdominal compartment syndrome: the cause or effect of postinjury multiple organ failure. *Shock* 2003, 20(6):483-492.
81. Moore FA: The role of the gastrointestinal tract in postinjury multiple organ failure. *Am J Surg* 1999, 178(6):449-453.
82. Eleftheriadis E, Kotzampassi K, Papanotas K, Heliadis N, Sarris K: Gut ischemia, oxidative stress, and bacterial translocation in elevated abdominal pressure in rats. *World J Surg* 1996, 20(1):11-16.
83. Ivy ME: Secondary abdominal compartment syndrome in burns. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 178-186.
84. Demling RH, Crawford G, Lind L, Read T: Restrictive pulmonary dysfunction caused by the grafted chest and abdominal burn. *Crit Care Med* 1988, 16(8):743-747.
85. Greenhalgh DG, Warden GD: The importance of intra-abdominal pressure measurements in burned children. *JTrauma* 1994, 36(5):685-690.
86. Hobson KG, Young KM, Ciraulo A, Palmieri TL, Greenhalgh DG: Release of abdominal compartment syndrome improves survival in patients with burn injury. *JTrauma* 2002, 53(6):1129-1133.
87. Ivy ME, Atweh NA, Palmer J, Possenti PP, Pineau M, D'Aluio M: Intra-abdominal hypertension and abdominal compartment syndrome in burn patients. *JTrauma* 2000, 49(3):387-391.
88. Ivy ME, Possenti PP, Kepros J, Atweh NA, D'Aluio M, Palmer J, Pineau M, Burns GA, Caushaj PF: Abdominal compartment syndrome in patients with burns. *JBurn Care Rehabil* 1999, 20(5):351-353.
89. Latenser BA, Kowal-Vern A, Kimball D, Chakrin A, Dujovny N: A pilot study comparing percutaneous decompression with decompressive laparotomy for acute abdominal compartment syndrome in thermal injury. *J Burn Care Rehabil* 2002, 23 190-195.
90. Ziakas PD, Voulgarelis M, Felekouras E, Anagnostou D, Tzelepis GE: Myelofibrosis-associated massive splenomegaly: a cause of increased intra-abdominal pressure, pulmonary hypertension, and positional dyspnea. *Am J Hematol* 2005, 80(2):128-132.
91. Hamad GG, Peitzman AB: Morbid obesity and chronic intra-abdominal hypertension. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 187-194.
92. Sugerman HJ: Increased intra-abdominal pressure in obesity. *IntJObesRelat Metab Disord* 1998, 22(11):1138.
93. Malbrain ML, De Laet IE, Willems A, Van Regenmortel N, Schoonheydt K, Dits H: Localised abdominal compartment syndrome: bladder-over-gastric pressure ratio (B/G ratio) as a clue to diagnosis. *Acta clinica Belgica* 2010, 65(2):98-106.
94. Balogh ZJ, Butcher NE: Compartment syndromes from head to toe. *Crit Care Med* 2010, 38(9 Suppl):S445-451.
95. Scalea TM, Bochicchio GV, Habashi N, McGunn M, Shih D, McQuillan K, Aarabi B: Increased intra-abdominal, intrathoracic, and intracranial pressure after severe brain injury: multiple compartment syndrome. *The Journal of trauma* 2007, 62(3):647-656.
96. De laet I, Citerio G, Malbrain ML: The influence of intraabdominal hypertension on the central nervous system: current insights and clinical recommendations, is it all in the head? *Acta Clin Belg Suppl* 2007, 62(1):89-97.
97. Malbrain ML, Wilmer A: The polycompartment syndrome: towards an understanding of the interactions between different compartments! *Intensive care medicine* 2007, 33(11):1869-1872.
98. Malbrain MLNG, De laet I: A new concept: the polycompartment syndrome - Part 2. *Int J Intensive Care* 2009, Spring 2009:19-25.
99. Malbrain MLNG, De laet I: A new concept: the polycompartment syndrome - Part 1. *Int J Intensive Care* 2008, Autumn 2008:19-24.
100. Reintam A, Parm P, Kitus R, Kern H, Starkopf J: Primary and secondary intra-abdominal hypertension-different impact on ICU outcome. *Intensive care medicine* 2008.
101. Vidal M, Ruiz Weissner E, Gonzalez F, Toro M, Loudet C, Balasini C, Canales H, Reina R, Estenssoro E: Incidence and clinical effects of intraabdominal hypertension in critically ill patients. *Crit Care Med* 2008, in press.
102. Malbrain ML, De Laet I: AIDS is coming to your ICU: be prepared for acute bowel injury and acute intestinal distress syndrome. *Intensive care medicine* 2008.
103. Kimball EJ: Intra-abdominal hypertension and the abdominal compartment syndrome: "ARDS" of the gut. *Int J Intensive Care* 2006, Spring:1-7.
104. De Laet IE, Hoste EA, De Waele JJ: Survey on the perception and management of the abdominal compartment syndrome among Belgian surgeons. *Acta Chir Belg* 2007, 107(6):648-652.
105. Mayberry JC: Prevention of abdominal compartment syndrome. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 221-229.
106. Parr M, Olvera C: Medical management of abdominal compartment syndrome. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 230-237.
107. Mertens zur Borg IR, Verbrugge SJ, Kolkman KA: Anesthetic considerations in abdominal compartment syndrome. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 252-263.
108. Balogh Z, Moore FA: Postinjury secondary abdominal compartment syndrome. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 170-177.
109. Balogh Z, Moore FA, Goettler CE, Rotondo MF, Schwab CW, Kaplan MJ: Management of abdominal compartment syndrome. In: *Abdominal compartment syndrome*. Edited by Ivatury R, Cheatham M, Malbrain M, Sugrue M. Georgetown: Landes Bioscience; 2006: 264-294.
110. De Laet IE, Ravvys M, Vidts W, Valk J, De Waele JJ, Malbrain ML: Current insights in intra-abdominal hypertension and abdominal compartment syndrome: open the abdomen and keep it open! *Langenbecks Arch Surg* 2008.
111. Cheatham ML, Safcsak K, Llerena LE, Morrow CE, Jr., Block EF: Long-term physical, mental, and functional consequences of abdominal decompression. *The Journal of trauma* 2004, 56(2):237-241; discussion 241-232.
112. Cheatham ML, Safcsak K: Is the evolving management of intra-abdominal hypertension and abdominal compartment syndrome improving survival? *Crit Care Med* 2010, 38(2):402-407.
113. Malbrain ML: Is it wise not to think about intraabdominal hypertension in the ICU? *Curr Opin Crit Care* 2004, 10(2):132-145.
114. Malbrain ML, Chiumello D, Pelosi P, Bihari D, Innes R, Ranieri VM, Del Turco M, Wilmer A, Brienza N, Malcangi V et al: Incidence and prognosis of intraabdominal hypertension in a mixed population of critically ill patients: a multiple-center epidemiological study. *Crit Care Med* 2005, 33(2):315-322.
115. Burch JM, Moore EE, Moore FA, Franciose R: The abdominal compartment syndrome. *Surg Clin North Am* 1996, 76(4):833-842.
116. Ivatury RR, Sugerman HJ, Peitzman AB: Abdominal compartment syndrome: recognition and management. *AdvSurg* 2001, 35:251-269.

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Adaptado de la Figura 1. El algoritmo de evaluación intraabdominal (pg. 954) según lo originalmente publicado en el artículo, resulta de la Conferencia Internacional de Expertos sobre Hipertensión Intraabdominal y el Síndrome Compartimental Abdominal. II. Recomendaciones; Intensive Care Medicine. Algoritmo reproducido con la autorización de la Sociedad Mundial del Síndrome Compartimental Abdominal y Springer Healthcare. ©Springer-Verlag 2007. Reservados todos los derechos.