



Application Card

How to Use ConvaTec Little Ones® One-Piece Extra Small Ostomy System*

Begin by cleansing the skin around your child's stoma with a mild, non-oily soap and warm water. If you prefer to use commercial wet wipes on the skin around the stoma, choose a brand that does not contain oils or lotions, since they make it difficult for the pouching system to stick. The skin should be free of any solvent or greasy substances. Pat dry thoroughly.

Skin Barrier/Pouch Preparation and Application



1 **Measure** the stoma size from the base of the stoma with the measuring guide provided.



2 **Trace** the stoma outline you have selected onto the release paper on the skin barrier. The small starter hole should be in the center of your tracing. (Gloving may not be necessary; consult with your health-care professional).



3 **Important:** Before cutting, separate the front and back panels of the pouch to avoid cutting through the pouch. Pouch should contain a small amount of air. Cut along the tracing with curved, short-blade scissors.



4 **Cut** the center hole the same size as the stoma, using guidelines on the release paper. Do not cut beyond the last line.



5 Prior to applying the pouch, it is recommended that you close the InvisiClose® clipless tail closure system by folding upwards and pressing the interlocking closures together until you feel them lock into place. A "click" indicates a secure closure.

Attachment of Pouch



6 Remove the release paper from the skin barrier.



7 Place the pouching system gently in position over your child's stoma.



8 Press the skin barrier against the skin, paying particular attention to the area closest to the stoma.



9 During application, maintain gentle pressure on the skin barrier for 30 seconds to improve adhesion.

Pouch Removal



10 To remove the old pouch and skin barrier, have a wet paper towel or soft washcloth ready to clean the skin. Hold the edge of the pouch with one hand and push down gently on the skin with your other hand to release the skin from the pouch and the skin barrier. Go slowly and gently.

Discard as necessary.

How to Use InvisiClose® Tail Closure System*

These instructions apply to the Little Ones® Extra Small Ostomy System

Directions to Close the Pouch



The InvisiClose® Tail Closure System is comprised of a **security flap (A)**, **interlocking closures (B)**, and **outlet end-strips (C)**.



Hold the security flap against the front of the pouch revealing the interlocking closures.



Fold the tail end up towards you four times and ensure the interlocking closures overlap.



Press the interlocking closures together until you feel them lock (**click**) into place along their entire length.



Wrap the security flap around to the back of the pouch (**body side**) and press the interlocking closures together until you feel them lock (**click**) into place.

Directions for Draining the Pouch



Turn the back of the pouch (body side) towards you. Carefully lift the security flap up towards you to unlock.



While holding the security flap up against the front of the pouch, unlock the interlocking closure by lifting at the corner and unraveling the folds four times.



To facilitate draining and/or cleaning, you may open the tail end by pushing in both ends of the outlet end-strips with your thumb and forefinger. If desired, pouch tail may be cleaned with gentle soap and water. Pat dry.

Tips:

- Remember to empty your child's pouch when it is one-third full.
- The use of soaps or cleansers in the pouch is not recommended. If desired, rinse only with tepid water.

Little Ones and InvisiClose are registered trademarks of ConvaTec Inc.
Realize the possibilities is a trademark of ConvaTec Inc.

*See package insert for complete Instructions for Use.
©2009 ConvaTec Inc. May 2009 AP-007505-US

For more information or if you have questions, call

1-800-422-8811

or

visit our website at www.ConvaTec.com