

Product Description

Exudate Management – For moderately to heavily exuding wounds

Versiva[®] XC[®] Gelling Foam dressing consists of a waterproof top polyurethane foam/film layer, which protects the wound from external contaminants and provides a bacterial/viral barrier¹. It also manages the moisture vapour transmission of the exudate absorbed by the dressing¹. An absorptive non-woven fibrous layer (with Hydrofiber[®] Technology) absorbs and retains exudate by forming a cohesive gel.

Versiva[®] XC[®] is available as an adhesive dressing or non-adhesive dressing. The adhesive dressing has a thin adhesive border that secures the dressing in place while allowing gentle removal. The non-adhesive dressing has a thin non-adhesive wound contact layer which allows gentle dressing removal.

The dressing absorbs and retains exudate creating a moist wound environment which supports the body's healing process and aids in the removal of non-viable tissue from the wound (autolytic debridement) without damaging new tissue¹.

Application*

Carefully cleanse the wound with saline and dry surrounding skin. Choose a dressing and size, ensuring the pad area is larger than the wound.

Adhesive Dressing - Remove release paper from the dressing, avoiding finger contact with the wound contact surface and the adhesive border. Place the pad directly over the wound, ensuring the wound is in the centre of the dressing pad, and smooth down the adhesive border.

Non-Adhesive Dressing - Remove dressing from the pack being careful to minimise finger contact with the wound contact surface. Place the pad directly over the wound ensuring the wound is in the centre of the dressing pad. An appropriate retention bandage or tape should be used to secure the dressing in place.

Versiva[®] XC[®] dressing may be used as a primary dressing or as a secondary dressing in combination with appropriate dressings such as, but not limited to, AQUACEL[®] or AQUACEL[®] Ag dressing. Change dressing when clinically indicated. Maximum recommended wear time is up to 7 days.

To remove, press down gently on the skin and carefully lift one corner of the dressing until it is no longer adhered to the skin. Continue until all edges are free. Carefully lift away and discard.

*Please refer to pack insert for full instructions, prior to use.

Specifications

Description	Pack Size	Product Code
<i>Non-adhesive</i>		
7.5cm x 7.5cm	10	410606
11cm x 11cm	10	410607
15cm x 15cm	5	410608
20cm x 20cm	5	410614
<i>Adhesive</i>		
10cm x 10cm	10	410609
14cm x 14cm	10	410610
19cm x 19cm	5	410615
22cm x 22cm	5	410611
<i>Adhesive Heel</i>		
18.5cm x 20.5cm	5	410612
<i>Adhesive Sacral</i>		
21cm x 25cm	5	410613

ARTG Number: 158691



Clinical Benefits

- Provides a moist wound environment
- Protects periwound skin and reduces the risk of maceration²
- Comforts patients over time whilst the dressing is *in situ* and upon removal¹
- Provides superior cushioning³
- Easy to apply and remove²
- Provides a bacterial/viral barrier¹

Indications

Versiva[®] XC[®] Adhesive Gelling Foam dressing is indicated for moderately to heavily exuding wounds.

- Chronic wounds: Leg ulcers, pressure ulcers, diabetic foot ulcers
- Acute wounds: Minor cuts, abrasions, lacerations and scalds, surgical wounds, traumatic wounds, second degree burns

References

1. Versiva XC Gelling Foam Dressing [package insert], Skillman, NJ: ConvaTec; 2005.
2. Munter KC, Vanscheidt W, Klovekorn. A prospective study on the use of a non-adhesive gelling foam dressing on exuding leg ulcers. *J Wound Care*. 2007; 16 (6): 261-265.
3. Griffiths et al. The cushioning effect of a new gelling foam dressing compared with a number of other foam dressings in an in vitro study. Poster presented at European Wound Management Association meeting; May 2007; Glasgow.



Call your local ConvaTec Territory Manager or ConvaTec Support Centre for further information:

Australia: 1800 335 276 or New Zealand: 0800 441 763

www.convatec.com

ConvaTec (Australia) Pty Limited. ABN 70 131 232 570. Unipark Monash, Building 2, Ground Floor, 195 Wellington Road, Clayton VIC 3170 Australia. PO Box 63, Mulgrave, Vic 3170. Phone: (03) 9239 2700 Facsimile: (03) 9239 2743. ConvaTec (New Zealand) Limited. AK2135265 Level 4, 369 Queen Street, Auckland 1010 New Zealand. PO Box 62663, Kalmia Street, Auckland 1544 New Zealand. Phone: (09) 306 8833 Facsimile: (09) 306 8831.

®/™Indicates trademark of ConvaTec Inc. © 2009 ConvaTec Inc. AP-006701-AU. March 2009.

Our world is what we make of it™

