Indications for use

The insulin™ is indicated for subcutaneous infusion of medication administered by injection. The infwelling catheter is indicated for single use.

CONTRAINDICATIONS

insulin™ is indicated for subcutaneous use only, insulin™ is not to be used in combination with an external pump, and is not to be used for IV infusion. insulin™ may not be used for the infusion of blood or blood products.

WARNINGS

- Carefully examine the instructions for use before inserting insulin™. Failure to follow instructions may result in pain or injury.
- insulin™ infwelling catheter is only sterile and non-pyrogenic if the packaging is undamaged and intact. Do not use insulin™ if the packaging is already opened or damaged.
- When using the insulin™ for the first time, please do so in the presence of a healthcare provider.
- Do not use a needle larger than 27 gauge when infusing

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OPSLAG EN AFVOER

PRECAUTIONS

- insulin™ is for single use and must be disposed of after use. Do not clean or re-sterilize.
- Be sure to cleanse the insertion site carefully prior to inserting the insulin™. Do not put disinfectants, perfumes, deodorants or other products containing alcohol or disinfectants in contact with the insulin™.
- Replace insulin™ if the adhesive tape becomes loose or is displaced from its original site. Use a new site until the first site has healed.
- When using insulin, insulin™ should be removed every third day. Activity levels, amounts of injections and medication type may vary how often insulin™ needs to be changed. Please consult with a healthcare provider.
- Be sure to rotate insertion sites. If self-administering medication, be sure to place the insulin™ in an area where it is easily accessible.
- If the insertion site becomes irritated or infected, replace the insulin™ and use a new site until the first site has healed.
- Be sure to change rapid acting insulin, check blood glucose level one to two hours after inserting a new insulin™. Due to varying lengths of duration of action and peak time of

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

PRECAUTIONS

mixed insulin, regular insulin and long acting insulin, please check with your healthcare provider regarding checking blood glucose levels after inserting a new insulin™. Be sure to check the insertion site several times a day and monitor blood glucose regularly.

- If infusing insulin and blood glucose levels become abnormally high, remove insulin™ and insert a new one at a new site. The soft cannula in insulin™ could be dislodged or partially collapsed.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

STORAGE & DISPOSAL

- Store insulin™ in a cool, dry location at room temperature. Do not store insulin™ in direct sunlight, or in high humidity or in a easily accessible.
- Dispose of the introducer needle in a proper sharps container according to local laws, to avoid the risk of needle stick.
- Dispose the use infwelling catheter according to local regulations for biohazardous waste.

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