- (ru) Эндотрахеальная трубка
- (da) Trakealtube
- (sv) Trakealtub
- (no) Tracheal Tube
- (fi) Intubaatioputki
- (de) Trachealtubus
- (fr) Sonde Trachéale
- (es) Tubo Traqueal
- (pt) Tubo Traqueal
- (it) Tubo Tracheale
- (nl) Tracheaal Tube
- (pl) Rurka intubacyjna
- (cs) Tracheální rourka
- (ro) Sonda traheala
- (hu) Trachea Tubus
- (gr) Τραχειακός Σωλήνας
- (ar) أنابيب القصبة الهوائية
- (tw) 氣管內管
- (th) ท่อช่วยหายใจ

## (en) Tracheal Tube

### This instruction for use applies to:

- Magill shaped Tracheal Tubes
- Preformed Oral and Nasal Tracheal Tubes cuffed and uncuffed (plain), with and without Murphy eye

### INTENDED USE / INDICATIONS

Tracheal Tubes are used in general anesthesia, intensive care and emergency medicine for airway management and mechanical ventilation. The tube is inserted into the patient's trachea through the nose or mouth in order to maintain the airway.

The Tracheal Tubes described in this document have the following features:

- Magill curve design, suitable for oral and nasal intubation;
- Preformed curved design, specifically indicated for Oral or Nasal intubation;
- Clear thermosensitive material;
- Soft rounded beveled tip;
- Variants with cuff: low pressure cuff, pilot balloon and inflation line;
- Optional Murphy eye on both cuffed and uncuffed versions;
- Optional preloaded Intubation Stylet on specific variants;
- Depth metric markers to indicate the distance to the distal tip. Cuffed versions feature ring markings above the cuff, to determine the position of cuff and distal tip after passing the glottis;
- Standard 15mm ISO connector, color-coded to indicate the size of the tube and to suggest the most suitable suction catheter (ISO 8836);
- Radiopaque line on the full length of the tube;
- Sterilized by ethylene oxide (EO) and packed in peelable pouch.

### CONTRAINDICATIONS

- It is contraindicated to use Tracheal Tubes in procedures that involve the use of a laser or an electrosurgical active electrode

in the immediate area of the device. Especially in presence of oxygen-enriched mixtures, this could result in a rapid combustion of the Tracheal Tube with harmful thermic effects and emission of corrosive and toxic products, including hydrochloric acid (HCl).

- Tracheal Tubes should not be used for patients suffering from throat edema/inflammation, hemorrhage or neck vertebra trauma.

### WARNINGS / PRECAUTIONS (General)

- Do not use if package has been previously opened, damaged or under inadequate storage.
- Avoid exposure to elevated temperatures, ultraviolet and fluorescent light.
- It is not recommended to cut Tracheal Tubes. If a Tracheal Tube is cut, this is at the healthcare professional's discretion and own clinical judgment.
- If patient positioning or Tracheal Tube placement are altered after intubation it is essential to verify if the tube is still in the correct position. Any tube displacement must be corrected immediately.
- Should extreme flexing of the head (chin-to-chest) or repositioning of the patient (e.g. to a lateral or prone position) be anticipated after intubation, the use of a Reinforced Tracheal Tube should be considered.
- Non-standard dimensions of some connectors on ventilators or anesthesia equipment may affect secure mating with the Tracheal Tube 15mm connector. Use only with equipment featuring standard 15mm connectors.
- When choosing the appropriate size of Tracheal Tube, healthcare professionals should pay attention to the anatomical variations of each patient, especially the length of the airway. Depth metric markers on the tube should never substitute an expert clinical judgment.

- Care should be taken when using a lubricant in conjunction with the Tracheal Tube. An excessive amount of lubricant can dry on the inner surface of the tube, forming a plug or a clear film, which can obstruct the airway.
- Do not use after the expiry date.
- For single patient use only.
- Do not resterilize. Resterilization or reuse can result in harm to the patient.
- Dispose of the device in accordance with the applicable national regulations for biological hazardous waste.

#### **WARNINGS / PRECAUTIONS (Cuff related)**

- Prior to use, test the Tracheal Tube's cuff, pilot balloon and valve by inflation. If a dysfunction is detected in any part of the inflation system, the tube must not be used.
- The use of lidocaine topical aerosol has been associated with the formation of pinholes in PVC cuffs (Jayasuriya, K.D., Watson, W.F., "P.V.C. Cuffs and Lidocaine-based aerosol", *Brit. J. Ann.* 53:1368, 1981). Expert clinical judgment must be used when prescribing treatment involving the use of this substance, to help prevent situations of cuff leaks due to pinholes. The same authors report that lidocaine hydrochloride solution does not have this effect.
- Various anatomical structures on the intubation route (e.g. teeth, turbinates) or any intubation tool with sharp surface might represent a threat for cuff integrity. Care must be taken to avoid damaging the thin-walled cuff during intubation, as the patient would be subject to the trauma of extubation and re-intubation with a new device. If the cuff is damaged, the Tracheal Tube must not be used.
- Diffusion of nitrous oxide ( $N_2O$ ) mixture, oxygen or air from a surrounding gas mixture may either increase or decrease

cuff volume and pressure. In order to limit such effect it is recommended to inflate the cuff with the same gas mixture, which will contact its external surface.

- Do not overinflate the cuff. Overinflation can result in the rupture of the cuff with subsequent deflation or cuff distortion, which might lead to airway blockage or prevent the microvascularization of the tissue around the cuff. For indications on the recommended cuff pressure, refer to the actual scientific literature.
- It is recommended to assess cuff pressure using a cuff-pressure measuring device, in conjunction with minimal occluding volume or minimum leak technique. Cuff pressure should be monitored and any deviation from the selected seal pressure should be investigated and corrected immediately.
- Prior to repositioning the tube, deflate the cuff. Moving the Tracheal Tube with inflated cuff could result in patient injury or damage to the cuff, requiring a new device. When the cuff is completely deflated, the syringe will show a definite vacuum and the Tracheal Tube pilot balloon will be collapsed. After repositioning the tube, verify its correct placement once again.
- Syringes, three-way stopcocks or other devices should not be left inserted in the inflation valve for an extended period. The valve housing could crack due to pressure and allow the cuff to deflate.

#### **ADVERSE REACTIONS**

Standard textbooks and scientific literature should be consulted for any specific adverse reaction information.

#### **DIRECTIONS FOR USE (General)**

- Use aseptic technique.
- Intubation and extubation should be performed following currently accepted

medical techniques.

- Expert clinical judgment should be exercised in the selection of the appropriate size of Tracheal Tube for each individual patient.
- If the Tracheal Tube is preloaded with Intubation Stylet, prior to intubation use the stylet to shape the tube. Remove the stylet when intubation is complete, before connecting the tube to the ventilation/anesthesia equipment.
- When intubation is completed, auscultate both lung fields. If breath sounds are reduced over one lung field or absent over one or both fields the Tracheal Tube has to be repositioned.
- If required, Tracheal Tube placement can be checked by means of a chest radiography.

#### DIRECTIONS FOR USE (Cuff related)

- Prior to intubation test the integrity of the cuff: inflate the cuff with a Luer-tip syringe, and then deflate it completely after test inflation.
- When intubation is completed, inflate the cuff using the minimum volume of air to provide an effective seal.
- After cuff inflation, remove the Luer-tip syringe from the valve to effect closing of the one-way valve.
- Prior to extubation, deflate the cuff completely. Extubate following currently accepted medical techniques.

#### CAUTION

- Single use only
- Sterile if package is unopened and undamaged
- Store in a cool, dry place.

#### TRACHEO-BRONCHIAL SUCTIONING

It is recommended to perform tracheal suctioning by utilizing a suction catheter/suction system of appropriate size, matching the color of the 15mm connector

present on the Tracheal Tubes.

Tracheal Tube	Recommended Suction Catheter size		Colour Coding
	I.D. (mm)	O.D. (mm)	
3.0 - 3.5	1.3	4	red
4.0 - 4.5	2.0	6	green
5.0 - 6.0	2.7	8	blue
6.5 - 7.0	3.3	10	black
7.5 - 8.0	4.0	12	white
8.5 - 9.0	4.7	14	green
9.5 - 10.0	5.3	16	orange

#### CUFF SEALING PERFORMANCE

Test was performed according to ISO 5361:2012. Data were collected using a bench test that is intended to provide a comparison of the sealing characteristics of Tracheal Tube cuffs only in a laboratory setting. The bench test is not configured or intended to predict performance in the clinical setting.

Tracheal Tube size I.D. (mm)	Minimum trachea diameter (mm)	Maximum trachea diameter (mm)	Leakage rate range (mL/h)			
			50th percentile	90th percentile	50th percentile	90th percentile
<b>HVLP cuff</b>						
3.0/3.5	9	11	54.9	884.1	0.0	191.3
4.0/4.5	11	13	723.2	1174.1	0.0	146.0
5.0/5.5	13	16	269.1	948.6	0.0	6.6
6.0/6.5	15	19	359.0	647.0	0.0	1.3
7.0/7.5	19	24	376.7	656.5	0.0	2.5
8.0/8.5	20	25	570.6	2082.4	0.0	0.0
9.0/9.5/10.0	21	27	625.9	1777.2	0.0	0.0
<b>Polyurethane cuff, without Murphy-eye</b>						
3.0	8	10	146.8	184.9	0.0	0.0
3.5/4.0	9	11	98.5	102.9	0.0	0.0
4.5/5.0	11	13	69.9	78.2	0.0	0.0

\*Cuff pressure hPa (cm H<sub>2</sub>O) = 27

	Consult instructions for use		Does not contain natural rubber latex
	Do not use if package is damaged		Does not contain phthalate
	Sterilized using ethylene oxide		Catalogue number
	Do not re-sterilize		Batch code
	Do not re-use		Manufacturer
	Keep dry		Use-by date

Authorized representative in the European Community

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Well Lead Medical Co.,Ltd.
Corp GmbH
ConvaTec Limited
Sponsored in Australia by ConvaTec (Australia) Pty Limited