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
Antimicrobial activity of AQUACEL™ Ag+ Extra™, AQUACEL™ Ag Extra™, Acticoat 7, Silvercel® Non Adherent and Urgotul® Ag/Silver dressings against CA-MRSA; illustration of kill and prevention of regrowth over 9 days.

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
Quantity of total viable microorganisms recovered following 24-hour dressing exposure to a 72-hour pre-formed multi-species biofilm.

Dressing: A = AQUACEL™ Ag+ Extra™, B = UrgoClean™ Ag, C = Exufiber™ Ag+, D = Maxorb™ Extra Ag+, E = Kerracel™ Ag

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
Wound status at baseline (light blue) and after evaluation (dark blue).

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
1. Sen CK, Gordillo GM, Roy S, et al. Human skin wounds: a major and snowballing threat to public health and the economy. *Wound Repair Regen.* 2009;17(6):763-71. 2. Ortman JM, Velkoff VA, Hogan H. An Aging Nation: The Older Population in the United States. U.S. Department of Commerce; 2014 [Available from: <https://www.census.gov/prod/2014pubs/p25-1140.pdf> (assessed 16 January 2018)]. 3. Rosenthal RJ, Morton J, Brethauer S, Mattar S, De Maria E, Benz JK, et al. Obesity in America. *Surg Obes Relat Dis.* 2017;13(10):1643-50. 4. Purwins S, Herberger K, Debus S, et al. Cost-of-illness of chronic leg ulcers in Germany. *Int Wound J.* 2010;7(2):97-102. 5. Hjort A, Gottrup F. Cost of wound treatment to increase significantly in Denmark over the next decade. *J Wound Care.* 2010;19(5):173-4, 176, 178, 180, 182, 184. 6. Posnett J, Gottrup F, Lundgren H, Saal G. The resource impact of wounds on health-care providers in Europe. *J Wound Care.* 2009;18(4):154-161. 7. Guest JF, Ayoub N, McLlraith T, et al. Health economic burden that wounds impose on the National Health Service in the UK. *BMJ Open.* 2015;5:e009283. doi:10.1136/bmjopen-2015-009283. 8. Tammelin, A, Lindholm, C, Hambraeus, A. Chronic ulcers and antibiotic treatment. *J Wound Care.* 1998; 7:435-437. 9. Gürgen, M. Excess use of antibiotics in patients with non-healing ulcers. *EWMA Journal.* 2014; 14:17-22. 10. Bowler PG. Antibiotic resistance and biofilm tolerance: A combined threat in the treatment of chronic infections. *J Wound Care.* 2018; 27:273-277. 11. Prompers L, Huijberts M, Schaper N, Apelqvist J et al. Resource utilisation and costs associated with the treatment of diabetic foot ulcers. Prospective data from the Eurodiale Study. *Diabetologia.* 2008; 51:1826-1834. 12. Malone M et al. 2017. The prevalence of biofilm in chronic wounds: a systematic review and meta-analysis of published data. *JWC.* 20-25. 13. Hurlow, Blanz, Gaddy. Clinical investigation of biofilm in non-healing wounds by high resolution microscopy techniques. *J Wound Care WUWHS Suppl* 2016; 25: S11-S22. 14. Wolcott R, Sanford N, Gabriliska R et al. Microbiota is a primary cause of pathogenesis of chronic wounds. *J Wound Care WUWHS Suppl* 2016;25(10): S33-S43. 15. Percival SL, Bowler PG. 2004. Biofilms and their potential role in wound healing. *WOUNDS.* 16: 234-240. 16. National Institutes of Health, 2002. Research on microbial biofilms. <https://grants.nih.gov/grants/guide/pa-files/PA-03-047.html>. Accessed October 2017. 17. Donlan RM, Costerton JW. Biofilms: survival mechanisms of clinically relevant microorganisms. *Clin Micro Rev.* 2002; 15:167-193. 18. Wolcott RD et al. Biofilm maturity studies indicate sharp debridement opens a time dependent therapeutic window. *J Wound Care.* 2010; 19:320-328. 19. Bowler. Antibiotic resistance and biofilm tolerance: a combined threat in the treatment of chronic infections. *JWC Vol 27*;No 5; 2018. 20. Hurlow, J, Couch, K, Laforet, K, Bolton, L, Metcalf, D. et al. (2015). Clinical Biofilms: A Challenging Frontier in Wound Care. *Advances in Wound Care.* 4(5), 295-301. 21. Metcalf, Bowler. Biofilm delays wound healing: a review of the evidence. *Burns Trauma* 2013; 1: 5-12. 22. Gurjala AN et al. Development of a novel, highly quantitative in vivo model for the study of biofilm-impaired cutaneous wound healing. *Wound Rep Reg* (2011) 19 400-410. 23. Hall-Stoodley L et al. Towards diagnostic guidelines for biofilm-associated infections. *FEMS Immunol Med Microbiol.* 2012; 65:127-145. 24. Costerton JW, Stewart PS, Greenberg EP. Bacterial Biofilms: A Common Cause of Persistent Infections. 1999; 284 Science. 25. Saye D.E. Recurring and antimicrobial-resistant infections: Considering the potential role of biofilms in clinical practice. *Ostomy Wound Management* 2007 53:4 (46-62). 26. Bowler PG, Parsons, D. Combatting wound biofilm and recalcitrance with a novel anti-biofilm Hydrofi ber® wound dressing. *Wound Medicine* 14 (2016) 6-11. 27. N Høiby et al. ESCMID guideline for the diagnosis and treatment of biofilm infections 2014. *Clin Microbiol Infect* 2015; 21: S1-S25. 28. Hurlow, Bowler. Potential implications of biofilm in chronic wounds: a case series. *J Wound Care* 2012; 21: 109-119. 29. Parsons D, Meredith K, Rowlands VJ, Short D, Metcalf DG, Bowler PG. Enhanced Performance and Mode of Action of a Novel Antibiofilm Hydrofi ber® Wound Dressing. *Biomed Res Int.* 2016;2016:7616471. 30. Said J, Walker M, Parsons D, Stapleton P, Beezer AE, Gaisford S. An in vitro test of the efficacy of an anti-biofilm wound dressing. *Int J Pharmaceutics.* 2014; 474: 177-181. DOI: 10.1016/j.ijpharm.2014.08.034. 31. Composition comprising antimicrobial metal ions and a quaternary cationic surfactant WO12136968 Parsons World patent application 11th October 2012. 32. Banin E, Brady K.M. & Greenberg E.P. (2006). Chelator Induced Dispersal and Killing of *Pseudomonas aeruginosa* Cells in Biofilm. *Appl. Environ. Microbiol.* 72. 2064-2069. 33. Chen X, Stewart PS, 2000. Biofilm removal caused by chemical treatments. *Wat. Res.* 34: 4229-4233. 34. Seth AK, Zhong A, Nguyen KT, Hong S.J, Leung KP, Galiano RD, Mustoe TA. Impact of a novel, antimicrobial dressing on in vivo, *Pseudomonas aeruginosa* wound biofilm: quantitative comparative analysis using a rabbit ear model. *Wound Repair Regen.* 2014; 22: 712-719. DOI: 10.1111/wrrr.12232. 35. Hobot JA, Walker M, Newman GN, Bowler PG, 2008. Effect of Hydrofi ber® wound dressings on bacterial ultrastructure. *J Electr Micro;* 57: 67-75. 36. T. J. Beveridge, W. S. Fyfe. Metal fixation by bacterial cell walls. *Canadian Journal of Earth Sciences.* 1985, 22(12): 1893-1898. <https://doi.org/10.1139/e85-204>. 37. WHR15860 MA322-In Vitro Antimicrobial Activity of AQUACEL Ag + Extra and AQUACEL Ag Extra against Acticoat 7, Silvercel NA and Urgotul Ag -V1. ConvaTec Data on File. 27 November 2018. 38. Assessment of the Anti-biofilm Activity of Silver-containing gelling fibre wound dressings, using a multi-species CDC reactor model. Hannah Thomas, Perfectus Biomed. Presented at Wounds UK November 2018. <https://epostersonline.com/wounds2018/node/171> Accessed November 2018. 39. Metcalf DG, Parsons D, Bowler PG. Clinical safety and effectiveness evaluation of a new antimicrobial wound dressing designed to manage exudate, infection and biofilm. *Int Wound J* 2017; 14: 203-213.