

Soft convexity... for when extra comfort is needed



ESTEEM™+ Soft Convex



ESTEEM™+ Soft
Convex Closed



ESTEEM™+ Soft
Convex Drainable



ESTEEM™+ Soft
Convex Urostomy

Esteem™+ Soft Convex is the latest addition to our range of one-piece convexity solutions, enhancing your choices when selecting the best convexity option for your patient. The flexibility of Esteem™+ Soft Convex combined with a shallow convexity is designed to be gentle on the skin and deliver comfort to your patient.

Suitable for all patients, particularly those who are post-op, elderly or have a low BMI.

The portfolio includes closed, drainable and urostomy pouches that are either pre-cut or cut-to-fit with different cutting ranges. Helping to ensure suitability for various stoma shapes and sizes.

Key benefits include:

- **Comfort and flexibility**
Three convexity plateau widths provide the support required in the peristomal area to help both improve the seal and help minimise leaks
- **Peace of mind**
Split fabric, easy-view window for accurate positioning and observation
- **Trust and reassurance**
The filter on the closed and drainable pouches helps minimise odours and ballooning
- **Simple secure fastening**
The Safe Seal clipless closure on drainable options helps make cleaning and emptying easy
- **Support and security**
Integral soft flexible belt loops for added security

Three plateau widths

The Esteem™+ Soft Convex system is available in 3 plateau widths to address a range of stoma management issues.

Three plateau widths



V1 50mm



V2 40mm



V3 30mm

Esteem™+ Soft Convex **Closed**

Pouch size 145 x 205mm

Cut-to-fit

Product Code	Plateau Size	Cutting Range	Quantity
422360	V1	20-47mm	30
422361	V2	15-40mm	30
422362	V3	10-28mm	30

Pre-cut

Product Code	Plateau Size	Cutting Range	Quantity
422527	V1	25mm	30
422528	V1	30mm	30
422529	V1	35mm	30
422530	V1	40mm	30
422531	V1	45mm	30
422533	V2	25mm	30
422534	V2	30mm	30
422535	V2	35mm	30
422538	V3	25mm	30

Esteem™+ Soft Convex **Drainable**

Pouch size 150 x 220mm

Cut-to-fit

Product Code	Plateau Size	Cutting Range	Quantity
422363	V1	20-47mm	10
422364	V2	15-40mm	10
422365	V3	10-28mm	10

Pre-cut

Product Code	Plateau Size	Cutting Range	Quantity
422539	V1	25mm	10
422540	V1	30mm	10
422541	V1	35mm	10
422542	V1	40mm	10
422543	V1	45mm	10
422545	V2	25mm	10
422546	V2	30mm	10
422547	V2	35mm	10
422550	V3	25mm	10

Esteem™+ Soft Convex **Urostomy**

Pouch size 145 x 205mm

Cut-to-fit


Product Code	Plateau Size	Cutting Range	Quantity
422551	V1	20-47mm	10
422557	V2	15-40mm	10
422562	V3	10-28mm	10

Pre-cut

Product Code	Plateau Size	Cutting Range	Quantity
422552	V1	25mm	10
422553	V1	30mm	10
422554	V1	35mm	10
422555	V1	40mm	10
422556	V1	45mm	10
422559	V2	25mm	10
422560	V2	30mm	10
422561	V2	35mm	10
422564	V3	20mm	10
422565	V3	25mm	10

Urostomy Night Drainage Solution

Choose a night drainage bag to suit your needs, the Esteem™+ Urostomy pouches comes with connectors to fit a variety of night drainage bags.



Support and advice every step of the way

We aim to help people be much more than their ostomy.

We continue to work hard to innovate and bring new advances in product design, all focused on improving the life of every person living with a stoma. This commitment is providing people with new choices and allowing them the freedom to live the life they want.



Nurse Advisor Line
UK: 0800 085 2516
ROI: 1800 818 988



me+ Advisors
UK: 0800 467 866
ROI: 1800 721 721



Product Specialists
UK: 0800 282 254
ROI: 1800 721 721



Homecare Group
UK: 0800 885 050