



Cure Medical Product Catalog

[CLICK to request samples](#)




Intermittent Catheter Options

Product	Page
Cure Dextra® Closed System	4
Cure Catheter® Closed System	5
Cure Ultra® and Cure Ultra®Plus™	6
Cure Twist®	7
Cure Catheter® (uncoated)	8
Cure Catheter® (hydrophilic)	9
Cure Medical® Pocket Catheter	10
Cure Medical® Catheter Insertion Kits	11
Cure Medical® Extension Tube	11

[CLICK for instructions for use](#)

Size Options

To aid in product selection, the approximate outer dimension of the catheter French (FR) size is shown below. Additionally, the color of the funnel end corresponds to FR sizes for fast, easy identification.



Find Your Ideal Solution

The Cure Medical line of intermittent catheters includes options for men, women, and children. Traditional and innovative styles featuring advanced technologies are designed to provide a safe, clean, fast, and convenient cathing experience.


Cure Medical products offer multiple benefits:



COMFORT
Smooth, fire polished eyelets

Manufacturing processes and material choices for Cure Medical intermittent catheters are intended to offer an assurance of safety and comfort. For example, instead of using a less expensive eyelet 'punch' technique, all eyelets on all Cure Medical catheters are fire polished. This not only helps aid in smooth insertion, it is designed to help minimize friction that can lead to infection.

A photomacrograph image shows an eyelet on a Cure Medical catheter (middle) compared to an eyelet on other catheter brands.



The difference in quality between the Cure Medical catheter eyelet and the other brands is clearly seen – and it can be felt!



RESEARCH
We support medical research

Financial support for medical research has been a core differentiator since the founding of Cure Medical.

INNOVATIVE OPTIONS. BENEFICIAL SOLUTIONS.

Request free samples • www.curemedical.com • 800.422.8811



SAFETY

Not made with DEHP/DINP, BPA, or Natural Rubber Latex

By investing in advanced technology and using high quality manufacturing materials, Cure Medical helps catheter users experience safety, comfort, and peace of mind. Not every catheter manufacturer can make this statement.

Some other catheter manufacturers use Di(2-ethylhexyl) Phthalate (DEHP) or Diisononyl Phthalate (DINP) as an economical way to make catheters flexible – even though these phthalates are included in a published list* of chemicals known to "...cause cancer, birth defects, and reproductive harm".

Additionally, Natural Rubber (NR) Latex may be used by other catheter manufacturers – even though it can cause serious allergic reactions.

In contrast, Cure Medical catheter users can be assured that the product does not put them at risk of exposure to DEHP/DINP, BPA, or NR-Latex.



[CLICK to learn about DEHP](#)



VALUE

'Upgrade Quality' and 'Reimbursement Friendly'

Fire polished eyelets, high quality materials, and innovative product designs are just a few examples of features that you might expect on a premium priced product. However, these 'Upgrade Quality' Cure Medical catheter features are uniquely combined with 'Reimbursement Friendly' availability.**

*<https://www.p65warnings.ca.gov>

**US only. Respects most reimbursement models. Contact your local healthcare product supplier for exact pricing.



USER FRIENDLY

Patient-centric packaging, products, and support

With a foundation of catheter user input, Cure Medical product features are designed to facilitate ease of use and help increase independence.



Evidence of this pursuit can be seen in features such as: a blue Control Stripe on catheters with a coude tip to help guide insertion orientation; packaging with dual openings or reinforced Finger Loops; Grippers and Sleeves to enable clean, controlled 'No Touch' use, and proprietary product features that support individuals with limited dexterity.



The same patient-centric perspective that drives product and packaging innovations also influences support strategies for improving overall well-being.

A unique, social platform called Cure Nation is designed to provide relevant information and enable catheter users to share life enhancing articles online and through a free eNewsletter.



[CLICK to join the Cure Nation](#)

Resources such as product guides, instructional videos, and book downloads are offered online to help support catheter users. Additionally, a backpack with custom flashcards and coloring book are offered to help teach a child 'How To Cath'. Crayons and a trinket add some fun to the educational resource.




INNOVATIVE CATH-IN-A-BAG

Request free samples • www.curemedical-dextra.com • 800.422.8811

Cure Dextra® Closed System

The sterile, single-use Cure Dextra is an innovative closed system catheter featuring smooth, fire polished eyelets on a pre-lubricated catheter integrated with a 1000 ml collection bag. It enables 'No Touch' use which helps reduce the risk of exposure to bacteria that can lead to infection. Proprietary product features are designed to help make the Cure Dextra easy to use – even for individuals with limited hand dexterity.



- Reinforced finger loops on the package help enable easy opening.
- A Support Band on the back of the collection bag helps enable controlled use.
- The Pull Ring helps enable easy removal of the cap on the pre-lubricated introducer tip.
- Tip Advancing Technology with Gripper Arrow helps facilitate ease of use. The catheter advances into the body with each full forward stroke. It does not retract into the bag on the backstroke.
- Perforations in the collection bag help create an opening for easy draining.

Unisex, Single

Offered as a case of 90 units

Part No.	Description	FR Size
DEX12*	Closed System Single, Straight Tip	12
DEX14*	Closed System Single, Straight Tip	14
DEX16*	Closed System Single, Straight Tip	16

[CLICK to go to product webpage with video](#)

*Not available in Canada.

TRADITIONAL CATH-IN-A-BAG

Request free samples • www.curemedical.com • 800.422.8811

Cure Catheter® Closed System

The sterile, single-use Cure Catheter Closed System features smooth, fire polished eyelets on a pre-lubricated catheter with integrated with a 1500 ml collection bag. It is offered with options for straight or coude tip design. The Cure Catheter Closed System is offered as a Single and a Kit with insertion supplies including a nonstaining BZK wipe, povidone-iodine swabsticks, ambidextrous gloves, gauze wipe, and underpad.



Unisex, Single

Offered as a case of 100 units

Part No.	Description	FR Size
CB8	Closed System Single, Straight Tip	8
CB10	Closed System Single, Straight Tip	10
CB12	Closed System Single, Straight Tip	12
CB14	Closed System Single, Straight Tip	14
CB16	Closed System Single, Straight Tip	16

Unisex, Kit

Offered as a case of 100 units

Part No.	Description	FR Size
CS8*	Closed System Kit, Straight Tip	8
CS10*	Closed System Kit, Straight Tip	10
CS12*	Closed System Kit, Straight Tip	12
CS14*	Closed System Kit, Straight Tip	14
CS14C*	Closed System Kit, Coude Tip	14
CS16*	Closed System Kit, Straight Tip	16

*Not available in Canada.

Cure Ultra® and Cure Ultra® Plus™

The sterile, single-use Cure Ultra line of pre-lubricated, ready-to-use catheters has expanded to include an 8" length for women – the Cure Ultra Plus. This conveniently portable product family offers many beneficial features:

- Smooth, fire polished eyelets are designed for comfort and to help reduce irritation that can lead to infection.
- The flexible, attractive packaging and kink-resistant catheter help enable discreet portability – even in small spaces.
- The unique, ergonomic funnel shape offers a 'No Roll' benefit and helps enable controlled handling.
- To help ensure comfortable use, a CoverAll™ feature provides even distribution of gel lubricant over the catheter as it is removed from the package – for ready-to-use convenience.
- Use of minimum product and packaging material means less waste and a more eco-conscious choice compared to many other options.



Male, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
ULTRA M8	Male, 16", Straight Tip	8
ULTRA M10	Male, 16", Straight Tip	10
ULTRA M10C	Male, 16", Coude Tip	10
ULTRA M12	Male, 16", Straight Tip	12
ULTRA M12C	Male, 16", Coude Tip	12
ULTRA M14	Male, 16", Straight Tip	14
ULTRA M14C	Male, 16", Coude Tip	14
ULTRA M16	Male, 16", Straight Tip	16
ULTRA M16C	Male, 16", Coude Tip	16
ULTRA M18	Male, 16", Straight Tip	18
ULTRA M18C	Male, 16", Coude Tip	18

Female, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
ULTRA8	Female, Straight Tip	8
ULTRA10	Female, Straight Tip	10
ULTRA12	Female, Straight Tip	12
ULTRA14	Female, Straight Tip	14
ULTRA16	Female, Straight Tip	16

Female, Single, 8" length

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
ULTRAPlus8*	Female, 8" Straight Tip	8
ULTRAPlus10*	Female, 8", Straight Tip	10
ULTRAPlus12*	Female, 8", Straight Tip	12
ULTRAPlus14*	Female, 8", Straight Tip	14
ULTRAPlus16*	Female, 8", Straight Tip	16

CLICK to go to product webpage with video


*Not available in Canada.

READY-TO-USE

Request free samples • www.curemedical.com • 800.422.8811

Cure Twist®

The small, attractive Cure Twist is similar to the design of a cosmetic item – so it is an ideal choice for women who value portability and discretion. It features an easy, twist open top and smooth, fire polished eyelets on a straight tip. After use, the Cure Twist catheter can go back into the container for clean, secure, and discreet disposal.



The Cure Twist offers ready-to-use convenience and discreet portability. Water soluble gel helps enable comfortable use with a 'No Drip/No Mess' benefit. If needed, the universal funnel can easily connect to standard drainage options.

Female, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
T8	Female, Straight Tip	8
T10	Female, Straight Tip	10
T12	Female, Straight Tip	12
T14	Female, Straight Tip	14
T16	Female, Straight Tip	16



The Cure Twist Kit with insertion supplies includes a BZK wipe, ambidextrous gloves, underpad, and 1500 ml collection bag with universal connector.

Female, Kit

Offered as a case of 90 units (3 boxes of 30)

Part No.	Description	FR Size
T8K*	Female, Straight Tip w/Insert. Kit	8
T10K*	Female, Straight Tip w/Insert. Kit	10
T12K*	Female, Straight Tip w/Insert. Kit	12
T14K*	Female, Straight Tip w/Insert. Kit	14
T16K*	Female, Straight Tip w/Insert. Kit	16

[CLICK to go to product webpage with video](#)

*Not available in Canada.

UNCOATED CATHETERS

Request free samples • www.curemedical.com • 800.422.8811

Cure Catheter®

The sterile, single-use, uncoated Cure Catheter features smooth, fire polished eyelets for comfort. It is offered in male, female, and pediatric lengths. This product family offers several beneficial features:

- All Cure Catheter packages feature unique, dual openings for increased ease of use.
- The Cure Catheter with coude tip for men features a Control Stripe to help guide insertion orientation.
- Resources to help teach a child 'How To Cath' are available for pediatric Cure Catheter users.



Male, Straight Tip, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
M8	Male, 16", Straight Tip	8
M10	Male, 16", Straight Tip	10
M12	Male, 16", Straight Tip	12
M14	Male, 16", Straight Tip	14
M16	Male, 16", Straight Tip	16
M18	Male, 16", Straight Tip	18

Male, Coude Tip, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
M10C	Male, 16", Coude Tip	10
M12C	Male, 16", Coude Tip	12
M14C	Male, 16", Coude Tip	14
M16C	Male, 16", Coude Tip	16
M18C	Male, 16", Coude Tip	18

Female, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
F8	Female, 6", Straight Tip	8
F10	Female, 6", Straight Tip	10
F12	Female, 6", Straight Tip	12
F14	Female, 6", Straight Tip	14
F14NC	Female, 6", Straight Tip, No Connector	14
F16	Female, 6", Straight Tip	16
F18	Female, 6", Straight Tip	18

Pediatric, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
P8	Pediatric, 10", Straight Tip	8
P10	Pediatric, 10", Straight Tip	10
P12	Pediatric, 10", Straight Tip	12
P14	Pediatric, 10", Straight Tip	14

HYDROPHILIC CATHETERS

Request free samples • www.curemedical.com • 800.422.8811

Hydrophilic Cure Catheter®

The sterile, single-use, Hydrophilic Cure Catheter features smooth, fire polished eyelets for comfort. It has a smooth, slippery surface when wet. It is easy to hydrate the Hydrophilic Cure Catheter in the package:

- Compress the water packet to release purified water into the sealed catheter package.
- Open the package just enough to be able to reach the catheter funnel. For convenience, the package can be temporarily adhered to a vertical surface to keep the purified water contained within the package.



Male, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
HM12	Hydrophilic, Male, 16", Straight Tip	12
HM12C	Hydrophilic, Male, 16", Coude Tip	12
HM14	Hydrophilic, Male, 16", Straight Tip	14
HM14C	Hydrophilic, Male, 16", Coude Tip	14
HM16	Hydrophilic, Male, 16", Straight Tip	16
HM16C	Hydrophilic, Male, 16", Coude Tip	16
HM18C	Hydrophilic, Male, 16", Coude Tip	18

Male, Kit

Offered as a case of 90 units (3 boxes of 30)

The Hydrophilic Cure Catheter Kit with insertion supplies includes BZK wipe, ambidextrous gloves, underpad, and 1500 ml collection bag for men.

Part No.	Description	FR Size
HM12UK*	Hydro., Male, 16", Straight Tip, Kit	12
HM14UK*	Hydro., Male, 16", Straight Tip, Kit	14
HM16UK*	Hydro., Male, 16", Straight Tip, Kit	16



Pediatric, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
HP8	Hydrophilic, Pediatric, 10", Straight Tip	8
HP10	Hydrophilic, Pediatric, 10", Straight Tip	10
HP12	Hydrophilic, Pediatric, 10", Straight Tip	12
HP14	Hydrophilic, Pediatric, 10", Straight Tip	14

*Not available in Canada.

[CLICK to go to product webpage with video](#)

DISCREET PORTABILITY

Request free samples • www.curemedical.com • 800.422.8811

Cure Medical® Pocket Catheter

The sterile, single-use, U-shaped Cure Medical Pocket Catheter fits in most pant pockets and small spaces for discreet portability. The package features large flaps to help facilitate easy opening. The Cure Medical Pocket Catheter features smooth, fire polished eyelets for comfort. It is offered in multiple configurations to meet a variety of needs and preferences:

- (U) Single catheter
- (UL) Single catheter with lubricant
- (XL) Extra long, single catheter with lubricant and gripper sleeve
- (UK) Kit



Male, Single

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
M12U	Male, 16", Straight Tip	12
M12ULC*	Male, 16", Coude Tip with Lubricant	12
M14U	Male, 16", Straight Tip	14
M14XL*	Male, 25", Straight Tip with Lubricant and Gripper Sleeve	14
M14UL*	Male, 16", Straight Tip with Lubricant	14
M14ULC*	Male, 16", Coude Tip with Lubricant	14
M16U	Male, 16", Straight Tip	16
M16ULC*	Male, 16", Coude Tip with Lubricant	16

Male, Kit

Offered as a case of 90 units (3 boxes of 30)

The Cure Medical Pocket Catheter Kit with insertion supplies includes nonstaining BZK wipe, packet of lubricant, ambidextrous gloves, underpad, and 1500 ml collection bag.

Part No.	Description	FR Size
M14UK*	Male, 16", Straight Tip, Kit	14

[CLICK to go to product webpage with video](#)

*Not available in Canada.

*Not available in Canada.

CONVENIENT OPTIONS

Request free samples • www.curemedical.com • 800.422.8811

Cure Medical® Catheter Insertion Kit

When combined with a catheter (not included), the sterile, single-use, unisex Cure Medical Catheter Insertion Kit is an alternative to using a complete closed system product. Insertion Kits include a nonstaining BZK wipe, packet of lubricant, ambidextrous gloves, underpad, and one of two collection bag options.



Offered as a case of 100 units

Part No.	Description
K1	Insertion Kit supplies plus resealable disposal bag



Offered as a case of 100 units

Part No.	Description
K2-100	Insertion Kit supplies plus 1500 ml collection bag with universal connector




Offered as a case of 90 units

Part No.	Description
K2-90*	Insertion Kit supplies plus 1500 ml collection bag with universal connector

Cure Medical® Extension Tube

The non-sterile, flexible Cure Medical Extension Tube easily connects to the funnel of all Cure Medical catheters to help facilitate draining.






Offered as a case of 30 units.

Part No.	Description
ET1*	29" tube

*Not available in Canada.



CLICK to read testimonials



Cure Medical
1434 E. Lake Mead Parkway, Ste. 210
Henderson, NV 89015

phone 800.422.8811 or 702.564.2901
fax 702.564.2903
www.curemedical.com

Follow us on   

Cure Medical, the Cure logo, Cure Catheter, Cure Dextra, Cure Twist, Cure Ultra, and Cure Ultra Plus are registered trademarks of Cure Medical LLC.

April 2023 • 900,045-T-Digital • AP-62423-USA-ENGU-v1