



Convatec

**Convatec  
Mouldable  
Technology™ and  
Esteem synergy®**



# There's only one ConvaTec Mouldable Technology™ And it's changing lives!

## How It Works:

- The only skin barrier with patented mouldable adhesive for a custom fit every time.
- It's easy! Simply mould it with your fingers for a personalised fit – no scissors needed!
- Hydrocolloid flexible fabric collar provides an extra level of security and improves flexibility and comfort by gently adhering to sensitive peristomal skin.



Hydrocolloid flexible fabric collar

- Built in Rebounding Memory Technology™ helps the adhesive gently hug the stoma for a secure seal.
- Durahesive® Technology absorbs effluent and swells to gently turtleneck and hug the stoma – for even greater leakage protection.
- Flexible adhesive moves with the patient to maintain a secure seal, minimise leakage and maintain healthy skin!



## Multi-Centre Evaluation:<sup>1</sup> High Scores Across the Board!

Both Patients and Healthcare Professionals rated ConvaTec Mouldable Technology™ as “Excellent” or “Very Good” across 10 criteria which included:

- Ease of creating a custom fit
- Adherence
- Overall satisfaction
- Effective skin protection
- Overall comfort
- Easy to teach

### Overall Satisfaction Rated Excellent or Very Good

	Patient	HCP
Colostomy	76%	83%
Ileostomy	72%	93%
Urostomy	88%	83%

**87% of clinicians felt teaching  
the usage of ConvaTec Mouldable  
Technology™ was Easy.<sup>1</sup>**

# Your effective 2-piece solution **just got better!**



## The Low Pressure You Want:

- Because there are no plastic rings to snap together – **application is low pressure.**
- Fingertips fit under “landing zone” to **minimise removal pressure.**
- Ability to gently separate pouch from “landing zone” to view the stoma (or release gas).



## The Secure Seal You Trust:

- Innovative foam adhesive (Adhesive Coupling Technology™) **seals the pouch to the clear skin barrier landing zone.**
- Lab Tested Security – Seal maintained adhesion when tested at the equivalent of 3.5 times the capacity of the standard 30 cm (12 inch) pouch.<sup>†,2,3</sup>

<sup>†</sup> Standard Pouch

## The Low Profile & Flexibility That Builds Patient Confidence:

- Low profile design for patient comfort and discretion.
- A system and skin barrier that moves as the patient moves, this flexible ostomy system bends and curves to fit any abdomen size or shape for proven performance.

# Meeting Patient's Needs at Each Stage of Care

	NEEDS	BENEFITS
Immediately Post-surgery	<ul style="list-style-type: none"> <li>■ A low-pressure coupling system that is gentle on tender post-surgery abdomens</li> <li>■ A versatile product to meet the needs of a variety of patient types</li> <li>■ A system that is easy to use and teach</li> <li>■ A unique, patented skin barrier designed to help prevent leakage</li> </ul>	<ul style="list-style-type: none"> <li>■ Adhesive Coupling Technology™ places less stress on abdomen during application</li> <li>■ ConvaTec Mouldable Technology™ skin barriers can be adjusted to fit a variety of stoma shapes and sizes</li> <li>■ The ConvaTec Mouldable Technology™ skin barrier simply rolls back to create an opening for the stoma, without the need for scissors or templates</li> <li>■ Mouldable adhesive on the ConvaTec Mouldable Technology™ skin barrier reduces gaps between the skin barrier and stoma, reducing the likelihood of leakage</li> </ul>
Discharge/ Recovery	<ul style="list-style-type: none"> <li>■ A system that is easy to learn and self-manage</li> <li>■ A product that will transition from hospital to home</li> <li>■ A secure seal around the stoma</li> <li>■ A discreet, low-profile system</li> </ul>	<ul style="list-style-type: none"> <li>■ The Esteem synergy® system with ConvaTec Mouldable Technology™ is easy to learn and to teach</li> <li>■ Flexible adhesive on the ConvaTec Mouldable Technology™ skin barrier adjusts to the changing size and shape of post-surgery stomas</li> <li>■ Mouldable adhesive on the ConvaTec Mouldable Technology™ skin barrier simply rolls open to fit desired stoma size and shape; no cutting or templates required</li> <li>■ Forgiving adhesive reduces likelihood of gaps by conforming to unique stoma contours, minimising risk of leakage</li> <li>■ Esteem synergy® system is designed to be slim under clothing</li> </ul>
Ongoing Lifestyle	<ul style="list-style-type: none"> <li>■ Helping to prevent leakage that leads to peristomal skin breakdown</li> <li>■ Comfort and peace of mind</li> </ul>	<ul style="list-style-type: none"> <li>■ The ConvaTec Mouldable Technology™ skin barrier adhesive design provides multiple layers of durable protection against caustic effluent</li> <li>■ The Adhesive Coupling Technology™ creates an extremely secure seal</li> <li>■ Adhesive moves with the ever-changing stoma as it functions throughout the day and night</li> </ul>





**ConvaTec**

**Because Your  
Patients Deserve  
Flexibility, Security  
& Adaptability**



# Ordering Information

Esteem synergy® mouldable wafers

## Esteem synergy® Flat Mouldable

Stoma Range	Color Code	Skin Barrier Ordering Number	Pack Size	Pouch Size
13-22 mm	●	411657	10	48 mm
22-33 mm	●	411658	10	48 mm
33-45 mm	●	411659	10	61 mm

## Esteem synergy® Convex Mouldable

Stoma Range	Color Code	Skin Barrier Ordering Number	Pack Size	Pouch Size
13-22 mm	●	409268	10	35 mm
22-33 mm	●	409269	10	48 mm
33-45 mm	●	409270	10	61 mm

Works with a range of  
Esteem synergy® pouch options

**More information and product orders:**

**Australia: 1800 006 609**  
**connection.au@convatec.com**

**New Zealand: 0800 441 763**  
**connection.nz@convatec.com**

#### References:

1. Hoeflok, J Guy, D, Allen S, St-Cyr D. A prospective multicenter evaluation of a Mouldable stoma skin barrier. *Ostomy Wound Manage.* 2009;55(5):62-69.
2. To determine shear resistance under load and ease of removal of Esteem synergy® coupling system. A2-0087. March 12, 2002. Data on file, ConvaTec.
3. Determine capacities for Synergy pouches with medium collar. A2-0108. March 26, 2002. Data on file, ConvaTec.

ConvaTec (Australia) Pty Limited. ABN 70 131 232 570. Unipark Monash, Building 2, Ground Floor, 195 Wellington Road, Clayton VIC 3168 Australia. PO Box 63, Mulgrave, Vic 3170. Phone: (03) 9239 2700 Facsimile: (03) 9239 2743.

ConvaTec (New Zealand) Limited. AK2135265 Level 4, 369 Queen Street, Auckland 1010 New Zealand. PO Box 62663, Kalmia Street, Auckland 1544 New Zealand. Phone: (09) 306 8833 Facsimile: (09) 306 8831.

®/™Indicates trademarks of ConvaTec Inc.

© 2010 ConvaTec Inc. AP-008240-MM January 2010 O255.