

# Post-Operative Kit Portfolio Upgrades

Creating Healthy Bonds with Complete Ostomy Products, Care and Support

- + Simple post-operative solutions
- + Quickly and correctly get your patient into the right system after surgery
- + Reduce material waste
- + Contains our most advanced 2-piece technologies



## Post-Operative Kit Selector Quick Guide

Flange Type	Cutting Type	Closure	Cut/Mold Range (mm)	New SKU	Description	Collar
ACCORDION LOW PRESSURE	Cut-to-Fit	InvisiClose (drainable)	13-45	423548	Natura™ POKit LPA Invisi 70mm CTF 1X5	Acrylic
	Moldable		33-45	423549	Natura™ POKit LPA Invisi 70mm CMT 1X5	Hydrocolloid
	Cut-to-Fit	Accuseal (urostomy)	13-45	423550	Natura™ POKit LPA UroAccu 70mm CTF 1X5	Acrylic
			13-33	423551	Natura™ POKit LPA UroAccu 57mm CTF 1X5	
	Moldable		33-45	423552	Natura™ POKit LPA UroAccu 70mm CMT 1X5	Hydrocolloid
			22-33	423553	Natura™ POKit LPA UroAccu 57mm CMT 1X5	
FLAT	Cut-to-Fit	InvisiClose (drainable)	13-89	423554	Natura™ POKit Invisi 100mm CTF 1X5	No Collar
			13-57	423555	Natura™ POKit Invisi 70mm CTF 1X5	
			13-45	423556	Natura™ POKit Invisi 57mm CTF 1X5	
	Moldable	InvisiClose (drainable)	45-56	423557	Natura™ POKit Invisi 70mm CMT 1X5	Hydrocolloid
			33-45	423558	Natura™ POKit Invisi 57mm CMT 1X5	
			22-33	423559	Natura™ POKit Invisi 45mm CMT 1X5	




# The right barrier and pouching solutions help support a healthy bond

The most common ostomy self-care challenges are related to leakage, skin issues, frequency of pouch changes and difficulty caring for an ostomy.<sup>1</sup> Our post-operative kit options have been simplified to include only our best technologies, improving patient outcomes and generating less waste.

## Convatec Durahesive®

All barriers use Durahesive®, featuring initial tack, strong adhesion and extended wear time.



**30% of patients** exhibit peristomal skin complications within 90 days of surgery.<sup>1</sup>




**95% of people** who started on Convatec Moldable Technology™ kept their skin healthy<sup>2</sup>

## Convatec Moldable Technology™ Skin Barrier

Patented Rebounding Memory Technology™ – Customizable so that it molds to the stoma for a perfect fit every time

## Accordion

Our Accordion flange is a post-operative option designed to prevent pressure or trauma. The accordion flange lifts up to provide generous finger room to make coupling easier and more comfortable.



**2x the lifting height** of any other lifting flange in the market<sup>3</sup>



Created utilizing feedback from **1,000 WOC Nurses<sup>4</sup>** and **400 patients<sup>5</sup>**

\*Only for drainable pouches

## Pouches

Convatec pouches are designed for comfort, ease and security.

**InvisiClose™ Clipless Tail Closure\***: intuitive, easy to open without compromising security

**Wider Outlet\***: enables easier drainage

**Accuseal™ Tap**: a faucet type tap designed for easy manipulation and emptying



# TRANSITION GUIDE

Use the chart below to note the changes in our new, simplified Post-Operative Kits. All kits include 5 packs per box.



## DRAINABLE

### Natura™ Drainable **Cut-to-Fit**

Previous SKU	New SKU	Change to Durahesive®	Larger Cut Range (cuts to 45mm)	No Loop Rod	No Low Pressure Adapter
416919	423556	●	●		
416923	423556	●			●
416927	423555	●		●	●
416929	423554	●		●	●

### Natura™ Drainable **Moldable**

Previous SKU	New SKU	Change to Durahesive®	No Low Pressure Adapter
416917	423559		●
416918	423559	●	●
416921	423558		●
416922	423558	●	●
416925	423557		●
416926	423557	●	●

### Natura™ Drainable **Moldable Accordion**

Previous SKU	New SKU	Change to Durahesive®	New Molding Range (now 33-45mm)	Flange Size Change
421042	423549	●	●	
421043	423549	●		
421044	423549	●		
421045	423549	●		
421046	423549	●	●	
421047	423549	●	●	●

# UROSTOMY

## Natura™ Urostomy **Cut-to-Fit**

Previous SKU	New SKU	Change to Durahesive®	Accordion Flange for Low-pressure Support	Flange Size Change
401923	423551		●	●
401924	423550	●	●	

## Natura™ Urostomy **Moldable**

Previous SKU	New SKU	Change to Durahesive®	Molding Range (now 33-45mm)	No Low Pressure Adapter	Accordion Flange for Low-pressure Support	Flange Size Change
416932	423553			●	●	●
416934	423553	●		●	●	●
416940	423552			●	●	●
416942	423552	●		●	●	●
416948	423552		●	●	●	
416950	423552	●	●	●	●	

## ALSO AVAILABLE

### Convatec Loop Ostomy Rod

Designed for external use in loop ostomy surgery. For use with SUR-FIT Natura™ Ostomy and Surgical Post-operative Systems.

SKU	Flange Size	Size	Quantity	Sterile
022355	2 3/4" (70mm)	65mm	10	●
022356	4" (100mm)	90mm	10	●



For more information, or to enroll your patient in our me+™ program, scan this code or call 1-800-422-8811 or visit [convatec.com/ostomy](https://convatec.com/ostomy)