

Soft convexity... for when extra comfort is needed



ESTEEM™+ Soft Convex



ESTEEM™+ Soft
Convex Closed



ESTEEM™+ Soft
Convex Drainable

Esteem™+ Soft Convex is the latest addition to our range of one-piece convexity solutions, enhancing your choices when selecting the best convexity option for your patient. The flexibility of Esteem™+ Soft Convex combined with a shallow convexity is designed to be gentle on the skin and deliver comfort to your patient.

Suitable for all patients, particularly those who are post-op, elderly or have a low BMI.

The portfolio includes closed and drainable pouches with different cutting ranges suitable for various stoma shapes and sizes.


Key benefits include:

- **Comfort and flexibility**
Three convexity plateau widths provide the amount of pressure required in the peristomal area to help both improve the seal and help minimise leaks
- **Peace of mind**
Split fabric, easy-view window for accurate positioning and observation
- **Trust and reassurance**
The filter on the closed and drainable pouches helps minimise odours and ballooning
- **Simple secure fastening**
The Safe Seal clipless closure on drainable options helps make cleaning and emptying easy
- **Support and security**
Integral soft flexible belt loops for added security



Three plateau widths

The Esteem™+ Soft Convex system is available in 3 plateau widths to address a range of stoma management issues.



Support and advice every step of the way

We aim to help people be much more than their ostomy.

We continue to work hard to innovate and bring new advances in product design, all focused on improving the life of every person living with a stoma.

This commitment is providing people with new choices and allowing them the freedom to live the life they want.

To find out more about our ostomy solutions or if you would like to talk to us, please get in touch on **XXXX XXX XXX**.

Esteem™+ Soft Convex **Closed**

Pouch size 145 x 205mm

Cut-to-fit

Product Code	Plateau Size	Cutting Range	Quantity
422360	V1	20-47mm	30
422361	V2	15-40mm	30
422362	V3	10-28mm	30

Esteem™+ Soft Convex **Drainable**

Pouch size 150 x 220mm

Cut-to-fit

Product Code	Plateau Size	Cutting Range	Quantity
422363	V1	20-47mm	10
422364	V2	15-40mm	10
422365	V3	10-28mm	10