

UnoVac™

Wound drainage system

Intended use:

The UnoVac™ wound drainage system is a low-vacuum wound drainage system designed for post-operative wound drainage procedure. The UnoVac™ system is connected to a drain which is placed directly into the surgical wound. The product is intended for single use only.



List of product range:

UnoVac™	Length of drain mm	Length of connecting tube mm	Non-return valve (NRV) in bag	CH 10	CH 12	CH 14	CH 16	CH 18
				Inner/outer diameter of drain, mm				
				2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3	3.5/6.0
				Length of perforation, mm				
				78	98	116	112	112
REF*	REF*	REF*	REF*	REF*				
UnoVac™ with Ulmer drain and solid trocar, with universal drain connector 	500	1050	No	2601418x	2601618x	2601818x	2602018x	2602218x
	750		No	2601518x	2601718x	2601918x	2602118x	2602318x
	750		Yes	2605818x		2605918x		
UnoVac™ with Ulmer drain, without trocar, with universal drain connector 	750	1050	No	2605318x	2605418x	2605518x	2605618x	2605718x
UnoVac™ with Ulmer drain and solid trocar, with cut to fit connector 	750	1050	No	2604718x	2604818x	2604918x		2605218x
UnoVac™ without Ulmer drain, without trocar, with universal drain connector		1050	No				2601118x	
			Yes				2606118x	
UnoVac™ wound drainage bag			No			2601318x		
UnoVac™ wound drainage bag			Yes			2606018x		

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

- Contains latex? No
- Contains DEHP? Yes
- Universal drain connector Polyethylene, Styrolux,
Nitril butadien rubber
- Cut to fit connector Thermoplastic elastomeric
- Connecting tube Polyvinyl chloride
- Front part Polyurethane
- Reverse part Styrolux
- Non-return valves Silicone rubber
- Bag and drain Polyvinyl chloride
- Trocar Steel

Certification:

- CE-marking? Yes (CE 0086)
- Product classification Class IIA (EU) (Systems with drain)
Class I (EU) (Systems without drain)
- Quality system certified to ISO 13485
- Global Medical Device Nomenclature (GMDN). 35824

Sterilisation:

- Sterile? Yes
- Sterilized by Ethylene oxide
- Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
2601118x	*	*	10	391 x 304 x 243	2.4	0.031	120
2601318x	40	294 x 125 x 360	120	392 x 306 x 375	4.0	0.045	120
2601418x	*	*	10	391 x 304 x 243	2.6	0.031	120
2601518x	*	*	10	391 x 304 x 243	2.6	0.031	120
2601618x	*	*	10	391 x 304 x 243	2.6	0.031	120
2601718x	*	*	10	391 x 304 x 243	3.0	0.031	120
2601818x	*	*	10	391 x 304 x 243	2.8	0.031	120
2601918x	*	*	10	391 x 304 x 243	2.8	0.031	120
2602018x	*	*	10	391 x 304 x 243	2.9	0.031	120
2602118x	*	*	10	391 x 304 x 243	2.9	0.031	120
2602218x	*	*	10	391 x 304 x 243	3.0	0.031	120
2602318x	*	*	10	391 x 304 x 243	3.0	0.031	120
2604718x	*	*	10	391 x 304 x 243	2.6	0.031	120
2604818x	*	*	10	391 x 304 x 243	2.9	0.031	200
2604918x	*	*	10	391 x 304 x 243	2.8	0.031	200
2605218x	*	*	10	391 x 304 x 243	3.0	0.031	200
2605318x	*	*	10	391 x 304 x 243	2.2	0.031	120
2605418x	*	*	10	391 x 304 x 243	2.3	0.031	120
2605518x	*	*	10	391 x 304 x 243	2.4	0.031	120
2605618x	*	*	10	391 x 304 x 243	2.5	0.031	120
2605718x	*	*	10	391 x 304 x 243	2.6	0.031	120
2605818x	*	*	10	391 x 304 x 243	2.6	0.031	120
2605918x	*	*	10	391 x 304 x 243	2.8	0.031	120
2606018x	40	294 x 125 x 360	120	392 x 306 x 375	4.0	0.045	120
2606118x	*	*	10	391 x 304 x 243	2.6	0.031	120

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.