


Suction Catheters

ProFlo™ and ProFlo™ Plain with Funnel

Intended use:

Suction catheters are intended for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via tube before being connected to vacuum source.
The product is intended for single use only.



List of product range:

	Standard Material (Shore A 78) REF**	Item description	CH	O.D. mm	Length* mm	Color coding
ProFlo™, straight tip with Funnel, with metric cm-marking on catheter						
	2531010101xx	ProFlo™, straight, Extended Length, Metric	10	3.3	570	Black 
	2531011101xx	ProFlo™, straight, Extended Length, Metric	12	4.0	570	White 
	2531012101xx	ProFlo™, straight, Extended Length, Metric	14	4.7	570	Green 
	2531013101xx	ProFlo™, straight, Extended Length, Metric	16	5.3	570	Orange 
ProFlo™ Plain, straight tip with Funnel, with metric cm-marking on catheter						
	2531001101xx	ProFlo™ Plain, straight, Extended Length, Metric	10	3.3	570	Black 
	2531002101xx	ProFlo™ Plain, straight, Extended Length, Metric	12	4.0	570	White 
	2531003101xx	ProFlo™ Plain, straight, Extended Length, Metric	14	4.7	570	Green 
	2531004101xx	ProFlo™ Plain, straight, Extended Length, Metric	16	5.3	570	Orange 
ProFlo™, coudé (bent) tip with Funnel						
	2531006101xx	ProFlo™ Coudé	10	3.3	570	Black 
	2531007101xx	ProFlo™ Coudé	12	4.0	570	White 
	2531008101xx	ProFlo™ Coudé	14	4.7	570	Green 
	2531009101xx	ProFlo™ Coudé	16	5.3	570	Orange 

* Effective length

** Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Catheter material. Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 10749

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Suction Catheters

ProFlo™ and ProFlo™ Plain with Funnel

Packaging:

Suction Catheters ProFlo™ and ProFlo™ Plain with Funnel

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m³	
ProFlo™, straight tip with Funnel							
2531010101xx	100	145 x 74 x 734	1,000	391 x 304 x 743	9.4	0.088	2,400
2531011101xx	100	145 x 74 x 734	1,000	391 x 304 x 743	9.4	0.088	2,400
2531012101xx	100	145 x 74 x 734	1,000	391 x 304 x 743	10.5	0.088	2,400
2531013101xx	100	145 x 74 x 734	1,000	391 x 304 x 743	10.0	0.088	2,400
ProFlo™ Plain, straight tip with Funnel							
2531001101xx	100	145 x 74 x 734	1,000	391 x 304 x 743	13.8	0.088	2,000
2531002101xx	100	145 x 93 x 734	800	391 x 304 x 743	12.7	0.088	1,600
2531003101xx	100	145 x 74 x 734	1,000	391 x 304 x 743	13.8	0.088	2,000
2531004101xx	100	145 x 93 x 734	800	391 x 304 x 743	13.0	0.088	4,800
ProFlo™, coudé (bent) tip with Funnel							
2531006101xx	100	145 x 93 x 734	800	391 x 304 x 743	9.4	0.088	2,400
2531007101xx	100	145 x 93 x 734	800	391 x 304 x 743	9.4	0.088	2,400
2531008101xx	100	145 x 93 x 734	800	391 x 304 x 743	10.5	0.088	2,400
2531009101xx	100	145 x 93 x 734	800	391 x 304 x 743	10.0	0.088	2,400

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

