


# Oxygen tubing

**Intended use:**

Oxygen tubes are used to connect the oxygen source with the oxygenation device during administration of oxygen. The product is intended for single use only.



**List of product range:**

Oxygen tubing (Inner diameter 0.394 cm)		Tubing		
REF*	Description	Length, m	Type	Connector
400MM	Oxygen tubing with standard lumen	2.1	Standard lumen	Standard
3400MM	Oxygen tubing with no-crush lumen	2.1	Standard lumen	Fits-all
402MM	Oxygen tubing with no-crush lumen	2.1	No-crush lumen	Standard
3402MM	Oxygen tubing with no-crush lumen	2.1	No-crush lumen	Fits-all
404MM	Oxygen tubing with no-crush lumen	4.2	No-crush lumen	Standard
407MM	Oxygen tubing with no-crush lumen	7.6	No-crush lumen	Standard
3407MM	Oxygen tubing with no-crush lumen	7.6	No-crush lumen	Fits-all
408MM	Oxygen tubing with no-crush lumen	15.2	No-crush lumen	Standard

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

Contains latex? . . . . . No  
 Contains DEHP? . . . . . Yes  
 Tube . . . . . Polyvinyl chloride  
 Connectors . . . . . Polyvinyl chloride

**Certification:**

CE-marking? . . . . . Yes (CE 0120)  
 Product classification . . . . . Class IIA (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . . 12875

**Sterilisation:**

Sterile? . . . . . No  
 Shelf life . . . . . 5 years

# Oxygen tubing

## Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
400MM	*	*	50	381 x 254 x 203	2.27	0.020	1,250
3400MM	*	*	50	381 x 254 x 203	2.52	0.020	1,250
402MM	*	*	50	406 x 355 x 254	2.36	0.037	1,250
3402MM	*	*	50	381 x 254 x 203	2.52	0.020	1,250
404MM	*	*	50	406 x 355 x 254	4.40	0.037	1,250
407MM	*	*	25	355 x 406 x 305	4.04	0.044	500
3407MM	*	*	25	406 x 355 x 254	3.83	0.037	625
408MM	*	*	25	355 x 406 x 305	6.92	0.044	500

\* = Single packed - no shelf box

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.

