

# Oxygen Mask





## High Concentration

**Intended use:**

High concentration masks are intended for administration of oxygen. They are equipped with a reservoir bag allowing an adequate amount of oxygen to be available to meet unpredictable breathing patterns and tidal volumes. Oxygen concentration delivered is 80%-100% depending on the patient's breathing rate and tidal volume. The product is intended for single use only.



**Product range:**

|   | REF*  | Item                | Size      | Description   | Side valves | Tube length |
|---|---|---------------------|-----------|---|-------------|-------------|
|   | 93-105MM  | Rebreather Mask     | Adult     | Adult High Concentration Oxygen Mask<br>Rebreathing mask (elongated) with nose clip,<br>No side valves, No tubing   | 0           | -           |
|   | 93-104MM<br>93-104MMV<br>(Vanilla scented mask) | Rebreather Mask     | Adult     | Adult High Concentration Oxygen Mask,<br>Rebreathing mask (elongated) with nose clip,<br>No side valves, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector        | 0           | 2.1m (7')   |
|   | 93-225MM  | Rebreather Mask     | Pediatric | Pediatric High Concentration Oxygen Mask,<br>Rebreathing mask (elongated) with nose clip,<br>No side valves, No tubing  | 0           | -           |
|   | 93-224MM<br>93-224MMV<br>(Vanilla scented mask) | Rebreather Mask     | Pediatric | Pediatric High Concentration Oxygen Mask,<br>Rebreathing mask (elongated) with nose clip,<br>No side valves, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector    | 0           | 2.1m (7')   |
|  | 93-265MM  | Rebreather Mask     | Infant    | Infant High Concentration Oxygen Mask,<br>Rebreathing mask with nose clip,<br>No side valves, No tubing   | 0           | -           |
|   | 93-264MM  | Rebreather Mask     | Infant    | Infant High Concentration Oxygen Mask,<br>Rebreathing mask with nose clip,<br>No side valves, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector                   | 0           | 2.1m (7')   |
|  | 93-106MM  | Non-Rebreather Mask | Adult     | Adult High Concentration Oxygen Mask,<br>Non-rebreathing mask (elongated) with nose clip,<br>1 side valve, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector      | 1           | 2.1m (7')   |
|   | 93-226MM  | Non-Rebreather Mask | Pediatric | Pediatric High Concentration Oxygen Mask,<br>Non-rebreathing mask (elongated) with nose clip,<br>1 side valve, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector  | 1           | 2.1m (7')   |
|   | 93-268MM  | Non-Rebreather Mask | Infant    | Infant High Concentration Oxygen Mask,<br>Non-rebreathing mask with nose clip,<br>1 side valve, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector                 | 1           | 2.1m (7')   |
|  | 93-109MM  | Non-Rebreather Mask | Adult     | Adult High Concentration Oxygen Mask,<br>Non-rebreathing mask (elongated) with nose clip,<br>2 side valves, No tubing   | 2           | -           |
|   | 93-108MM<br>93-108MMV<br>(Vanilla scented mask) | Non-Rebreather Mask | Adult     | Adult High Concentration Oxygen Mask,<br>Non-rebreathing mask (elongated) with nose clip,<br>2 side valves, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector     | 2           | 2.1m (7')   |
|   | 93-230MM<br>93-230MMV<br>(Vanilla scented mask) | Non-Rebreather Mask | Pediatric | Pediatric High Concentration Oxygen Mask,<br>Non-rebreathing mask (elongated) with nose clip,<br>2 side valves, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector | 2           | 2.1m (7')   |
|   | 93-267MM  | Non-Rebreather Mask | Infant    | Infant High Concentration Oxygen Mask,<br>Non-rebreathing mask with nose clip,<br>2 side valves, No tubing  | 2           | -           |
|   | 93-266MM  | Non-Rebreather Mask | Infant    | Infant High Concentration Oxygen Mask,<br>Non-rebreathing mask with nose clip,<br>2 side valves, 2.1m (7') No-crush oxygen tube,<br>Fits-All connector                | 2           | 2.1m (7')   |

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

# Oxygen Mask

## High Concentration



### Material specifications:

**Contains latex?** . . . . . No  
**Contains phthalates?** . . . . . No  
**Contains bisphenols?** . . . . . No  
**Mask** . . . . . Polyvinyl chloride  
**Elastic strap** . . . . . Polyester  
**Nose clip** . . . . . Aluminum  
**Oxygen nipple connector** . . . . . High density polyethylene  
**Oxygen tube** . . . . . Polyvinyl chloride  
**Tube connector** . . . . . Polyvinyl chloride  
**Reservoir bag** . . . . . Polyvinyl chloride  
**Check valve** . . . . . Silicone  
**Side valves** . . . . . Silicone

### Certifications:

**CE-marking?** . . . . . Yes (CE 0120)  
**Product classification** . . . . . Class IIA (EU)  
**Quality system certified to** . . . . . ISO 13485

### Sterilization:

**Sterile?** . . . . . No  
**Shelf life** . . . . . 5 years

### Packaging specifications:

| REF       | Primary package | Secondary package/Polybag | Tertiary package/Shipper |                        |               | Pallet   |
|-----------|-----------------|---------------------------|--------------------------|------------------------|---------------|----------|
|           | Qty, pcs        | Qty, pcs                  | Qty, pcs                 | Dimensions, mm (LxWxD) | Gross Wt., kg | Qty, pcs |
| 93-105MM  | 1               | 10                        | 50                       | 400 x 300 x 240        | 2.09          | 1,600    |
| 93-104MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 4.27          | 900      |
| 93-104MMV | 1               | 10                        | 50                       | 400 x 380 x 260        | 4.27          | 900      |
| 93-225MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 2.30          | 1,600    |
| 93-224MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 3.51          | 900      |
| 93-224MMV | 1               | 10                        | 50                       | 400 x 380 x 260        | 1.45          | 900      |
| 93-265MM  | 1               | 10                        | 50                       | 350 x 300 x 200        | 2.09          | 1,600    |
| 93-264MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 3.30          | 900      |
| 93-106MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 4.27          | 900      |
| 93-226MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 3.51          | 900      |
| 93-268MM  | 1               | 10                        | 50                       | 400 x 300 x 240        | 3.30          | 1,600    |
| 93-109MM  | 1               | 10                        | 50                       | 400 x 300 x 240        | 2.09          | 1,600    |
| 93-108MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 4.27          | 900      |
| 93-108MMV | 1               | 10                        | 50                       | 400 x 380 x 260        | 4.27          | 900      |
| 93-230MM  | 1               | 10                        | 50                       | 400 x 380 x 260        | 3.51          | 900      |
| 93-230MMV | 1               | 10                        | 50                       | 400 x 380 x 260        | 3.51          | 900      |
| 93-267MM  | 1               | 10                        | 50                       | 350 x 300 x 200        | 1.45          | 1,600    |
| 93-266MM  | 1               | 10                        | 50                       | 400 x 300 x 240        | 3.30          | 1,600    |



# Oxygen Mask

## High Concentration







### EAN codes:

| REF       | Primary package    | Secondary package  | Tertiary package   |
|-----------|--------------------|--------------------|--------------------|
| 93-105MM  | (01)06945099301116 | (01)26945099301110 | (01)46945099301114 |
| 93-104MM  | (01)06945099301147 | (01)26945099301141 | (01)46945099301145 |
| 93-104MMV | (01)06945099301369 | (01)26945099301363 | (01)46945099301367 |
| 93-225MM  | (01)06945099301154 | (01)26945099301158 | (01)46945099301152 |
| 93-224MM  | (01)06945099301161 | (01)26945099301165 | (01)46945099301169 |
| 93-224MMV | (01)06945099301376 | (01)26945099301370 | (01)46945099301374 |
| 93-265MM  | (01)06945099301185 | (01)26945099301189 | (01)46945099301183 |
| 93-264MM  | (01)06945099301208 | (01)26945099301202 | (01)46945099301206 |
| 93-106MM  | (01)06945099306050 | (01)26945099306054 | (01)46945099306058 |
| 93-226MM  | (01)06945099306067 | (01)26945099306061 | (01)46945099306065 |
| 93-268MM  | (01)06945099306074 | (01)26945099306078 | (01)46945099306072 |
| 93-109MM  | (01)06945099301123 | (01)26945099301127 | (01)46945099301121 |
| 93-108MM  | (01)06945099301130 | (01)26945099301134 | (01)46945099301138 |
| 93-108MMV | (01)06945099301352 | (01)26945099301356 | (01)46945099301350 |
| 93-230MM  | (01)06945099301178 | (01)26945099301172 | (01)46945099301176 |
| 93-230MMV | (01)06945099301383 | (01)26945099301387 | (01)46945099301381 |
| 93-267MM  | (01)06945099301192 | (01)26945099301196 | (01)46945099301190 |
| 93-266MM  | (01)06945099301215 | (01)26945099301219 | (01)46945099301213 |

### Packaging material:

|                   |  |
|-------------------|--|
| Primary package   | Poly bag - Low density polyethylene / Insert - 80g paper     |
| Secondary package | Poly bag - Low density polyethylene / Sticker - Coated paper |
| Tertiary package  | Carton - Corrugated board / Sticker - Coated paper           |

### Storage and handling:

|   |   |
|---|---|
|  | <b>FRAGILE</b> - Contents of the transport package are fragile and shall therefore be handled with care |
|  | <b>KEEP DRY</b>   |
|  | <b>THIS WAY UP</b> - Indicates correct upright position of the transport package                        |
|  | <b>DO NOT USE HAND HOOKS</b> - Hooks are prohibited for handling the transport package                  |



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.  
Hangzhou Bay New Zone, 315336, Ningbo City  
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg  
Germany

Distributed by



ConvaTec Limited  
First Avenue, Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

[www.convatec.com](http://www.convatec.com)

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2019 ConvaTec Inc. AP-019546-MM, ver. FEB 2019