





Aerosol Tubing and Accessories

Intended use:

An aerosol tube is used to connect an oxygen source to a nebulization device, such as an aerosol mask or face tent in order to administer medication via respiration. Aerosol mouthpieces and T-pieces are interfaces to connect to nebulization devices to administer aerosolized medication or for delivering oxygen. The products are intended for single use only.



Product range:

| | REF* | Item | Description | Tube length |
|---|-----------|--------------------|---|-------------|
|  | 93-600MM | Aerosol Tube | Aerosol Tube 30.5m, Cuttable every 15cm | 30.5m |
| | 93-626MM | Aerosol Tube | Aerosol Tube 1.8m, Cuttable every 15cm | 1.8m |
|  | 93-865MM | Aerosol Mouthpiece | - | - |
|  | 93-1379MM | Aerosol T-piece | - | - |

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Contains bisphenols? No
Aerosol tube. Ethylene vinyl acetate copolymer,
 Low density polyethylene
Aerosol mouthpiece High density polyethylene
Aerosol T-piece High density polyethylene

Certifications:

CE-marking? Yes (CE 0120)
Product classification IIA
Quality system certified to EN ISO 13485:2012

Sterilization:

Sterile? No
Shelf life 5 years

Packaging specifications:

| REF | Primary package | Secondary package/Polybag | Tertiary package/Shipper | | | Pallet |
|-----------|-----------------|---------------------------|--------------------------|------------------------|---------------|----------|
| | Qty, pcs | Qty, pcs | Qty, pcs | Dimensions, mm (LxWxD) | Gross Wt., kg | Qty, pcs |
| 93-600MM | 1 | - | 1 | 400 x 380 x 260 | 2.03 | 18 |
| 93-626MM | 1 | 5 | 10 | 400 x 400 x 300 | 1.41 | 180 |
| 93-865MM | 1 | 10 | 50 | 300 x 200 x 200 | 0.82 | 3,200 |
| 93-1379MM | 1 | 10 | 50 | 300 x 200 x 200 | 1.02 | 3,200 |





EAN codes:

| REF | Primary package | Secondary package | Tertiary package |
|-----------|--------------------|--------------------|--------------------|
| 93-600MM | (01)06945099313010 | - | (01)46945099313018 |
| 93-626MM | (01)06945099313027 | (01)26945099313021 | (01)46945099313025 |
| 93-865MM | (01)06945099313041 | (01)26945099313045 | (01)46945099313049 |
| 93-1379MM | (01)06945099313034 | (01)26945099313038 | (01)46945099313032 |

Packaging material:

| | |
|-------------------|--|
| Primary package | Poly Bag - Low density polyethylene / Insert - 80g Paper |
| Secondary package | Poly Bag - Low density polyethylene / Sticker - Coated Paper |
| Tertiary package | Sticker - Coated Paper / Carton - Corrugated Board |

Storage and handling:

| | |
|--|---|
|  | FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care |
|  | KEEP AWAY FROM RAIN - Transport package shall be kept away from rain |
|  | THIS WAY UP - Indicates correct upright position of the transport package |
|  | DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package |



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017332-MM, ver. FEB 2017