


# Oxygen Catheters

**Intended use:**

Oxygen catheters are intended for direct administration of oxygen via the nasopharyngeal route. The product is intended for single use only.



**List of product range:**

REF*	Type	Tube Colour	Size	Length cm	Lateral Eyes	Connector	Tube Shore A Hardness
1900518x	Plain	Green	CH 08	40	8	Male	74
1900618x	Plain	Green	CH 10	40	8	Male	74
1900718x	Plain	Green	CH 12	40	8	Male	74
1900818x	Plain	Green	CH 14	40	8	Male	74
1900918x	Plain	Green	CH 16	40	8	Male	74
1902018x	Plain	Green	CH 10	40	8	Male	74
1902818x	Plain	White	CH 10	50	8	Funnel Black	74
1903418x	Plain	Blue	CH 08	40	8	Male	74
1903518x	Plain	Blue	CH 10	40	8	Male	74
1903718x	Plain	Blue	CH 12	40	8	Male	74
1903818x	Plain	Blue	CH 14	40	8	Male	74
1901018x	Moist, Fixed Compress	Clear	CH 24	28	-	Funnel Grey	64
1901118x	Moist, Fixed Compress	Clear	CH 20	28	-	Funnel Yellow	64
1901318x	Fixed Compress	White	CH 10	40	-	Male	74
1901418x	Fixed Compress	White	CH 12	40	-	Male	74
1901518x	Fixed Compress	White	CH 14	40	-	Male	74
1901718x	Fixed Compress	Green	CH 10	40	-	Male	74
1901818x	Fixed Compress	Green	CH 12	40	-	Male	74
1901918x	Fixed Compress	Green	CH 14	40	-	Male	74
1903918x	Fixed Compress	Blue	CH 08	40	-	Male	74
1904018x	Fixed Compress	Blue	CH 10	40	-	Male	74
1904118x	Fixed Compress	Blue	CH 12	40	-	Male	74
1904218x	Fixed Compress	Blue	CH 14	40	-	Male	74
1902918x	Adjustable Compress	Blue	CH 08	40	8	Male	74
1902418x	Adjustable Compress	Green	CH 10	40	8	Male	74
1903018x	Adjustable Compress	Blue	CH 10	40	8	Male	74
1903118x	Adjustable Compress	Blue	CH 12	40	8	Male	74
1903218x	Adjustable Compress	Blue	CH 14	40	8	Male	74

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

**Contains latex?** ..... No  
**Contains DEHP?** ..... Yes  
**Tube** ..... Polyvinyl chloride  
**Funnel** ..... Polyvinyl chloride  
**Male connector** ..... Polyvinyl chloride  
**Compress** ..... Polyurethane foam

**Certification:**

**CE-marking?** ..... Yes (CE 0086)  
**Product classification** ..... IIA rule 5 (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device Nomenclature (GMDN)** ..... 35203

**Sterilisation:**

**Sterile?** ..... Yes  
**Sterilized by** ..... Ethylene oxide  
**Shelf life** ..... 5 years

# Oxygen Catheters

## Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
1900518x	100	145 x 74 x 446	1,000	391 x 304 x 455	6.7	0.056	3,000
1900618x	100	145 x 74 x 446	1,000	391 x 304 x 455	7.5	0.056	3,000
1900718x	100	145 x 74 x 446	1,000	391 x 304 x 455	8.3	0.056	3,000
1900818x	100	145 x 74 x 446	1,000	391 x 304 x 455	9.1	0.056	3,000
1902018x	100	145 x 74 x 446	1,000	391 x 304 x 455	7.5	0.056	3,000
1903418x	100	145 x 74 x 446	1,000	391 x 304 x 455	5.0	0.056	3,000
1903518x	100	145 x 74 x 446	1,000	391 x 304 x 455	5.3	0.056	3,000
1903718x	100	145 x 74 x 446	1,000	391 x 304 x 455	5.7	0.056	3,000
1903818x	100	145 x 74 x 446	1,000	391 x 304 x 455	6.4	0.056	3,000
1901018x	100	145 x 124 x 446	600	391 x 304 x 455	8.9	0.056	1,800
1901118x	100	145 x 124 x 446	600	391 x 304 x 455	7.2	0.056	1,800
1901318x	100	145 x 124 x 446	600	391 x 304 x 455	5.3	0.056	1,800
1901418x	100	145 x 124 x 446	600	391 x 304 x 455	5.7	0.056	1,800
1901518x	100	145 x 124 x 446	600	391 x 304 x 455	6.4	0.056	1,800
1901718x	100	145 x 124 x 446	600	391 x 304 x 455	5.3	0.056	1,800
1901818x	100	145 x 124 x 446	600	391 x 304 x 455	5.7	0.056	1,800
1901918x	100	145 x 124 x 446	600	391 x 304 x 455	6.4	0.056	1,800
1903918x	100	145 x 124 x 446	600	391 x 304 x 455	5.0	0.056	1,800
1904018x	100	145 x 124 x 446	600	391 x 304 x 455	5.3	0.056	1,800
1904118x	100	145 x 124 x 446	600	391 x 304 x 455	5.7	0.056	1,800
1904218x	100	145 x 124 x 446	600	391 x 304 x 455	6.4	0.056	1,800
1902918x	100	145 x 124 x 446	600	391 x 304 x 455	5.3	0.056	1,800
1902418x	100	145 x 124 x 446	600	391 x 304 x 455	5.3	0.056	1,800
1903018x	100	145 x 124 x 446	600	391 x 304 x 455	4.6	0.056	1,800
1903118x	100	145 x 124 x 446	600	391 x 304 x 455	6.4	0.056	1,800
1903218x	100	145 x 124 x 446	600	391 x 304 x 455	6.4	0.056	1,800
1900918x	100	145 x 93 x 446	800	391 x 304 x 455	8.7	0.056	2,400
1902818x	100	145 x 74 x 576	1,000	391 x 304 x 585	9.2	0.072	3,000

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.

