



AT A GLANCE

Helps minimise pressure on the abdomen

Accordion Flange™ skin barriers lift a generous finger's width away from the body so your patients can attach their pouch with minimal pressure on their abdomen.



- Intuitive design so patients can quickly get back to their everyday life
- **Lift 'UP'** flange to easily attach the pouch away from the body
- Designed to be suitable for people with compromised dexterity or limited eyesight



- Low profile in the 'DOWN' position so it's not visible under clothing
- Skin barrier moves as your patient's body flexes

Comprehensive portfolio to suit a wide range of patients

Tailored solutions to suit patients' lifestyles and a variety of abdomens and stomas, ensuring continuum of care from hospital to home and throughout their life with an ostomy.

Designed for intuitive ease of use

Smart features to help patients quickly adapt to managing an ostomy, including the lift-UP flange for easy coupling away from the body and a click lock for security and confidence.



- Designed to help minimise abdominal pressure and pain when attaching the pouch
- **Audible 'click' and physical 'snap'** to confirm the pouch is firmly attached



NATURA™ Accordion Flange™ Flat ConvaTec Moldable Technology™

CODE	COLOUR CODE	DESCRIPTION	STOMA SIZE	FLANGE SIZE	QTY
421030	●	Moldable Stomahesive™ Wafer 2 1/4" (57mm)	1/2" - 7/8" (13-22mm)	2 1/4" (57mm)	5
421031	●	Moldable Stomahesive™ Wafer 2 1/4" (57mm)	7/8" - 1 1/4" (22-33mm)	2 1/4" (57mm)	5
421032	●	Moldable Stomahesive™ Wafer 2 3/4" (70mm)	1 1/4" - 1 3/4" (33-45mm)	2 3/4" (70mm)	5
421033	●	Moldable Stomahesive™ Wafer 2 1/4" (57mm)	1/2" - 7/8" (13-22mm)	2 1/4" (57mm)	10
421034	●	Moldable Stomahesive™ Wafer 2 1/4" (57mm)	7/8" - 1 1/4" (22-33mm)	2 1/4" (57mm)	10
421035	●	Moldable Stomahesive™ Wafer 2 3/4" (70mm)	1 1/4" - 1 3/4" (33-45mm)	2 3/4" (70mm)	10
421036	●	Moldable Duraheasive™ Wafer 2 1/4" (57mm)	1/2" - 7/8" (13-22mm)	2 1/4" (57mm)	5
421037	●	Moldable Duraheasive™ Wafer 2 1/4" (57mm)	7/8" - 1 1/4" (22-33mm)	2 1/4" (57mm)	5
421038	●	Moldable Duraheasive™ Wafer 2 3/4" (70mm)	1 1/4" - 1 3/4" (33-45mm)	2 3/4" (70mm)	5
421039	●	Moldable Duraheasive™ Wafer 2 1/4" (57mm)	1/2" - 7/8" (13-22mm)	2 1/4" (57mm)	10
421040	●	Moldable Duraheasive™ Wafer 2 1/4" (57mm)	7/8" - 1 1/4" (22-33mm)	2 1/4" (57mm)	10
421041	●	Moldable Duraheasive™ Wafer 2 3/4" (70mm)	1 1/4" - 1 3/4" (33-45mm)	2 3/4" (70mm)	10



NATURA™ Accordion Flange™ Flat Cut-To-Fit

CODE	COLOUR CODE	DESCRIPTION	STOMA SIZE	FLANGE SIZE	QTY
421453	●	Duraheasive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	5
421454	●	Duraheasive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	10
421455	●	Stomahesive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	5
421456	●	Stomahesive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	10
421457	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	5
421458	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	10
421459	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	5
421460	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	10
421461	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	5
421462	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	10
421463	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	5
421464	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	10



NATURA™ Accordion Flange™ Convex Cut-To-Fit

CODE	COLOUR CODE	DESCRIPTION	STOMA SIZE	FLANGE SIZE	QTY
421632	●	Duraheasive™ Skin Barrier with Acrylic Collar (45mm)	1/2" - 7/8" (13-22mm)	1 3/4" (45mm)	5
421633	●	Duraheasive™ Skin Barrier with Hydrocolloid Collar (45mm)	1/2" - 7/8" (13-22mm)	1 3/4" (45mm)	5
421634	●	Duraheasive™ Skin Barrier with Acrylic Collar (57mm)	1/2" - 1 1/4" (13-35mm)	2 1/4" (57mm)	5
421635	●	Duraheasive™ Skin Barrier with Hydrocolloid Collar (57mm)	1/2" - 1 1/4" (13-35mm)	2 1/4" (57mm)	5
421636	●	Duraheasive™ Skin Barrier with Acrylic Collar (70mm)	1/2" - 1 3/4" (13-48mm)	2 3/4" (70mm)	5
421637	●	Duraheasive™ Skin Barrier with Hydrocolloid Collar (70mm)	1/2" - 1 3/4" (13-48mm)	2 3/4" (70mm)	5
421638	●	Duraheasive™ Skin Barrier with Acrylic Collar (45mm)	1/2" - 7/8" (13-22mm)	1 3/4" (45mm)	10
421639	●	Duraheasive™ Skin Barrier with Hydrocolloid Collar (45mm)	1/2" - 7/8" (13-22mm)	1 3/4" (45mm)	10
421640	●	Duraheasive™ Skin Barrier with Acrylic Collar (57mm)	1/2" - 1 1/4" (13-35mm)	2 1/4" (57mm)	10
421641	●	Duraheasive™ Skin Barrier with Hydrocolloid Collar (57mm)	1/2" - 1 1/4" (13-35mm)	2 1/4" (57mm)	10
421642	●	Duraheasive™ Skin Barrier with Acrylic Collar (70mm)	1/2" - 1 3/4" (13-48mm)	2 3/4" (70mm)	10
421643	●	Duraheasive™ Skin Barrier with Hydrocolloid Collar (70mm)	1/2" - 1 3/4" (13-48mm)	2 3/4" (70mm)	10

SUPPORT AND ADVICE EVERY STEP OF THE WAY

To find out more about our ostomy solutions or if you'd like to talk to us, please get in touch on

XXXX XXX XXX

©2017 ConvaTec Inc.™ Indicates a trademark of a ConvaTec Company AP-017586-MM



ConvaTec NATURA™ ACCORDION FLANGE™

SIMPLY UPLIFTING

Jane

Enjoying Thursday night yoga again

More space for coupling

Reassuring 'click-lock'


Low profile in 'DOWN' position

The Natura™ Accordion Flange™

HELPS MINIMISE PRESSURE ON ABDOMEN

COMPREHENSIVE RANGE TO SUIT ALL PATIENTS

INTUITIVE AND EASY-TO-USE DESIGN



NATURA™ ACCORDION FLANGE™

Allows people living with ostomy to get on with living life to their fullest.

Three types of Accordion Flange™ skin barriers, available with a wide range of Natura+™ comfort pouch variants.

Designed to make pouch application much easier for people living with an ostomy.

Available in Flat Cut-To-Fit, Convex Cut-to-Fit and Flat ConvaTec Moldable Technology™.




1. Natura™ Accordion Flange™ Flat Cut-To-Fit

The Natura™ Accordion Flange™ Flat Cut-To-Fit skin barrier is designed for regular-shaped abdomens and protruding stomas.

The flange lifts gently away from the body to securely click and lock with our Natura+™ pouch. Patients can simply cut the skin barrier to fit their stoma shape and size.

Flat Cut-To-Fit is available in Stomahesive™ and Durahesive™ with Acrylic Collar



- The flange lifts a generous finger's width away from the body
- Helping patients transition smoothly from hospital to their home

2. Natura™ Accordion Flange™ Convex Cut-To-Fit

The Natura™ Accordion Flange™ Convex Cut-To-Fit skin barrier is designed for flush, retracted, loop or telescopic stomas.

The skin barrier pushes gently down on peristomal skin to help protrude a flush or retracted stoma. Additionally, the Accordion Flange™ lifts a finger's width away from the abdomen to make coupling as easy and comfortable as possible.

Convex Cut-To-Fit is available in Durahesive™ with Hydrocolloid and Acrylic Collar



Convexity shown to improve quality of life³



- Easier to manage a variety of stoma and abdomen sizes and shapes¹
- Presses skin evenly to improve the protrusion of the stoma²
- Designed to help minimise leaks for healthier skin³

1. Rolstad BS. Principles and techniques in the use of convexity. Ostomy Wound Manage. 1996; 42(1):24-32.
2. Hanley J. Convex stoma appliances: are we getting it right? British Journal of Nursing (Stoma Supplement). 2013;12-16
3. Cronin E. A guide to the appropriate use of convex stoma care products. Gastrointestinal nursing.2008;6(2)Accessed 29 Feb.2016.

3. Natura™ Accordion Flange™ Flat Moldable

The Natura™ Accordion Flange™ is available with our unique ConvaTec Moldable Technology™ which gently molds to your patient's stoma. Clinically proven skin protection¹ for regular-shaped abdomens and protruding stomas.

Flat ConvaTec Moldable Technology™ is available in Stomahesive™ and Durahesive™ with Hydrocolloid Collar



Over 95% of people who started with Moldable Technology™ kept healthy skin⁴



- Soft, pliable elastic-like seal that rolls open to fit the size and shape of your patient's stoma.
- Flexible but strong enough to hold firm, comfortably 'hugging' the stoma to help minimise gaps and leaks
- Designed to provide a secure and snug fit that helps minimise leaks for healthier skin⁴

References 4. Clinical Summary: Osmose Study. Multinational Evaluation of the Peristomal Condition in New Ostomates Using Moldable Skin Barriers. Data on file. ConvaTec Inc.