






# Precision Handles

**Intended use:**

Precision Handles for the removal of fluids during surgery.  
The product is intended for single use only.



**List of product range:**

Precision handle	REF*	Bending	Connector	Length mm	Inner/Outer diameter of suction handle, mm	CH/FG (handle)	Vacuum control	Number of relieving eyes
	3300118x	Single	Clear	120	1.3/2.0	6	No	0
	3300218x	Double	Clear	120	1.3/2.0	6	No	0
	3300418x	Single	Funnel	115	1.3/2.0	6	No	0
	3300518x	Double	Funnel	115	1.3/2.0	6	No	0

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

Contains latex? . . . . . No  
 Contains DEHP? . . . . . Yes  
 Handle . . . . . Polyvinyl chloride  
 Clear connector . . . . . Polyvinyl chloride  
 Funnel connector . . . . . Polyvinyl chloride

**Certification:**

CE-marking? . . . . . Yes (CE 0086)  
 Product classification . . . . . Class IIA (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . . 35917

**Sterilisation:**

Sterile? . . . . . Yes  
 Sterilized by . . . . . Ethylene oxide  
 Shelf life . . . . . 5 years

# Precision Handles

## Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume m <sup>3</sup>	
3300118x	100	145 x 74 x 234	1,000	391 x 304 x 243	2.9	0.029	3,000
3300218x	100	145 x 74 x 234	1,000	391 x 304 x 243	2.9	0.029	3,000
3300418x	100	145 x 74 x 234	1,000	391 x 304 x 243	2.9	0.029	3,000
3300518x	100	145 x 74 x 234	1,000	391 x 304 x 243	2.9	0.029	3,000

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
First Avenue  
Address: Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

