


Universal™ drainage sets

Wound drainage system

Intended use:

Low-vacuum wound drainage system for continuous postoperative drainage. The product is intended for single use only.



List of product range:

Universal™ drainage system, includes suction bellows, a 1200 mm connecting tube and clamp	Length of drain mm	Volume of suction bellows, ml	CH 10	CH 12	CH 14	CH 16
			Inner/outer diameter of drain, mm			
			2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3
			Length of perforation, mm			
			80	100	120	120
Connecting tube: inner/outer diameter 3.5/6.0 mm, Shore A 79 Drains: Shore A 69			REF*	REF*	REF*	REF*
Universal™ drainage set with a Ulmer drain, solid trocar and Y-connector	1,000	400		1304518x		
Universal™ drainage set with an Ulmer drain, solid trocar and silicone cut-to-fit connector	500	400		1304418x		
Universal™ drainage set with an Ulmer drain, solid trocar and cut-to-fit connector	500	100	1305318x			
	500	400	1302718x	1302818x	1302918x	1303018x
Universal™ suction bellows only	-	100	1300518x			
	-	400	1300618x			
Accessories	Length of drain mm		CH 10	CH 12	CH 14	
Ulmer drain with connector and clamp for Universal™ bellows Drain: Shore A 69			Inner/outer diameter of drain, mm			
			2.0/3.33	2.2/4.0	2.4/4.7	
			Length of perforation, mm			
			80	100	116	
	REF*	REF*	REF*			
	500		1509518x	1509618x	1509718x	
			REF*			
Connector for connecting tube and 2 drains, CH 08-16			8406618x			

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Suction bellows Polyethylene
 Drains and connecting tube Polyvinyl chloride
 Trocar Stainless steel
 Connector for 2 drains Silicone
 Guide clamp Acrylonitrile-butadiene-styrene
 Closure ring Polypropylene

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 35824

Sterilisation:

Sterile? Yes
 Sterilised by Ethylene oxide
 Shelf life 5 years

Unoversal™ drainage sets

Wound drainage system

Packaging:

REF*	Shelf box		Shipping carton				MOQ. pcs.
	Qty., pcs	Dimensions. mm	Qty., pcs	Dimensions. mm	Gross Wt., kg	Volume m ³	
1300518x	25	297 x 190 x 366	100	391 x 304 x 743	4.1	0.092	200
1300618x	10	297 x 190 x 366	40	391 x 304 x 743	4.1	0.092	120
1302718x	10	297 x 190 x 366	40	391 x 304 x 743	6.3	0.092	40
1302818x	10	297 x 190 x 366	40	391 x 304 x 743	6.3	0.092	120
1302918x	10	297 x 190 x 366	40	391 x 304 x 743	6.3	0.092	40
1303018x	10	297 x 190 x 366	40	391 x 304 x 743	6.3	0.092	40
1304418x	10	297 x 190 x 366	40	391 x 304 x 743	7.0	0.092	120
1304518x	10	297 x 190 x 366	40	391 x 304 x 743	7.0	0.092	120
1305318x	10	297 x 190 x 366	40	391 x 304 x 743	6.3	0.092	120
1509518x	100	297 x 190 x 366	400	391 x 304 x 743	8.1	0.092	800
1509618x	100	297 x 190 x 366	400	391 x 304 x 743	9.0	0.092	1,200
1509718x	100	297 x 190 x 366	400	391 x 304 x 743	10.1	0.092	1,200
8406618x	100	145 x 93 x 234	800	391 x 304 x 243	3.8	0.031	800

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

