

Neria™Guard

Intended use

The product is intended for the subcutaneous infusion of immunoglobulins for the treatment of primary immune deficiency, apomorphine and foslevodopa/foscarbidopa for parkinson's disease, morphine (hydromorphone, morphine sulfate and morphine chloride) for pain management therapy, and deferoxamine mesylate for thalassemia. The product has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion. Clinical evaluation supports infusion set use from age 2 years and up.



Product specifications

Product code	Description	Cannula length	Tubing length	GTIN (Inner box)	Pcs per inner box	Pcs per outer box
704012-5226	Neria™Guard infusion set	6 mm	12 cm	05705244021242	10	80
704030-5226	Neria™Guard infusion set	6 mm	30 cm	05705244021303	10	80
704060-5226	Neria™Guard infusion set	6 mm	60 cm	05705244021365	10	80
704080-5226	Neria™Guard infusion set	6 mm	80 cm	05705244021426	10	80
704110-5226	Neria™Guard infusion set	6 mm	110 cm	05705244021488	10	80
704012-5229	Neria™Guard infusion set	9 mm	12 cm	05705244021273	10	80
704030-5229	Neria™Guard infusion set	9 mm	30 cm	05705244021334	10	80
704060-5229	Neria™Guard infusion set	9 mm	60 cm	05705244021396	10	80
704080-5229	Neria™Guard infusion set	9 mm	80 cm	05705244021457	10	80
704110-5229	Neria™Guard infusion set	9 mm	110 cm	05705244021518	10	80

For further information on these products, please contact your local Convatec sales representative.

Neria™ Guard

Product specifications:

Insertion angle	90 degrees
Soft cannula length	6 mm and 9 mm
Soft cannula outer diameter	~ 0,68 mm
Introducer needle	27 gauge
Tubing length in cm	12, 30, 60, 80 and 110
Tubing inner diameter	0,385 mm

Material specifications:

PVC free	Yes
Made with natural rubber latex	No
Raw materials or derivation of animal origin	No
The product contains phthalates	No

Certifications:

CE-marking	Yes
Product classification	Class IIB (EU)
Quality system certified to	ISO 13485
Tariff code	HS9018390000
GMDN	35833





Sterilization:

Sterile	Yes
Sterilized by	Ethylene oxide
Shelf life	3 years
Usage time (drug depended)	Maximum 72 hours

Drug compatibility/stability tested with¹:

- Immunoglobulins
- Apomorphine
- Foslevodopa/foscarbidopa,
- Morphine
((hydromorphone, morphine sulphate & morphine chloride)
- Deferoxamine mesylate

Storage and handling

	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	DO NOT RE-USE
	DO NOT USE IF PACKAGE IS DAMAGED
	CAUTION - See Instruction for Use



Manufacturer:
Unomedical A/S - A Convatec company
 Aaholmvej 1-3, Østød,
 4320 Iejre, Denmark
www.convatec.com

For further information on these products, please contact your local Convatec sales representative.
 Please refer to the instructions for use for additional information, warnings, and precautions.
 With reservations for misprints and alterations in product specifications. Printed version not controlled.

AP- AP-68087-GBL-ENG-v1 (v1.3)

© 2024 Convatec

Neria is a trademark of Unomedical A/S. The Convatec name, Hug icon and strapline are trademarks of Convatec.

Reference: 1. Convatec IC, data on file, D012-002912, 2023.