


Corrugated drains

Wound drains

Intended use:

Corrugated drains are intended for open passive drainage, used in situations where wounds must be kept open. They are placed directly into the wounds left to heal by secondary intent. The product is intended for single use only.



List of product range:

Corrugated drainage sheet with safety pin

| | REF* | Drain length, mm | Drain width, mm |
|--|----------|------------------|-----------------|
| | 1508518x | 250 | 25 |
| | 1508618x | 400 | 25 |

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Corrugated drains Polyvinyl chloride
 Safety pin Hard steel with Nickel coating

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 11305

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Packaging:

| REF* | Shipping carton | | | | MOQ, pcs. |
|----------|-----------------|-----------------|---------------|------------------------|-----------|
| | Qty., pcs. | Dimensions, mm | Gross Wt., kg | Volume, m ³ | |
| 1508518x | 600 | 391 x 304 x 455 | 10.0 | 0.056 | 600 |
| 1508618x | 600 | 391 x 304 x 585 | 11.6 | 0.072 | 600 |

Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.