

Intended use

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency.

Neria™ Multi has been tested to be compatible with immunoglobulins approved for subcutaneous infusion.



Product specifications

Product code	Description	Needle length	Tubing length	Furcated	GTIN (Inner box)	Pcs per inner box	Pcs per outer box
721060-2738	Neria™ Multi infusion set	8 mm	60 cm + 30 cm	2	05705244002074	10	120
721060-2731	Neria™ Multi infusion set	10 mm	60 cm + 30 cm	2	05705244002098	10	120
721060-3738	Neria™ Multi infusion set	8 mm	60 cm + 30 cm	3	05705244002111	10	120
721060-3731	Neria™ Multi infusion set	10 mm	60 cm + 30 cm	3	05705244002135	10	120
721060-3737	Neria™ Multi infusion set	12 mm	60 cm + 30 cm	3	05705244002159	10	120
721060-4738	Neria™ Multi infusion set	8 mm	60 cm + 30 cm	4	05705244002173	10	120
721060-4731	Neria™ Multi infusion set	10 mm	60 cm + 30 cm	4	05705244002197	10	120
721060-4737	Neria™ Multi infusion set	12 mm	60 cm + 30 cm	4	05705244002210	10	120

For further information on these products, please contact your local Convatec sales representative.

Neria™Multi

Product specifications:

Insertion angle	90 degrees
Needle length	8 mm, 10 mm and 12 mm
Needle gauge	27 gauge
Tubing length in cm	Total of 90 cm
Sets available as	bi-, tri- and quad-furcated sets
Tubing inner diameter (after furcater)	0,385 mm
Tubing inner diameter (supply tube)	0.650 mm
Tubing outer diameter	1,50 mm

Material specifications:

PVC free	Yes
Made with natural rubber latex	No
Raw materials or derivation of animal origin	No
The product contains phthalates	No

Certifications:

CE-marking	Yes
Product classification	Class IIB (EU)
Quality system certified to	ISO 13485
Tariff code	HS9018390000
GMDN	35833

Drug compatibility/stability tested with¹:

- Immunoglobulins

Storage and handling

	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	DO NOT RE-USE
	DO NOT USE IF PACKAGE IS DAMAGED
	CAUTION - See Instruction for Use



Unomedical A/S - A Convatec company
 Aaholmvej 1-3, Osted,
 4320 lejre, Denmark
www.convatec.com

Made in Mexico

For further information on these products, please contact your local Convatec sales representative. Please refer to the instructions for use for additional information, warnings, and precautions.
 With reservations for misprints and alterations in product specifications. Printed version not controlled.

AP-68086-GBL-ENG-v1

© 2024 Convatec

Neria is a trademark of Unomedical A/S. The Convatec name, Hug icon and strapline are trademarks of Convatec.

Reference: 1. Neria™Multi instructions for use, INS-809-3 2019-11