

# neria™ multi

## Intended use:

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency.

neria™ multi has been tested to be compatible with immunoglobulins approved for subcutaneous infusion.



## Product range:

Product code	Description	Needle length	Tubing length	Furcated	GTIN (Inner box)	Pcs per inner box	Pcs per outer box
<b>neria™</b>							
721060-2731	neria™ multi infusion set	10 mm	60 cm + 30 cm	2	05705244002098	10	120
721060-2738	neria™ multi infusion set	8 mm	60 cm + 30 cm	2	05705244002074	10	120
721060-3737	neria™ multi infusion set	12 mm	60 cm + 30 cm	3	05705244002159	10	120
721060-3731	neria™ multi infusion set	10 mm	60 cm + 30 cm	3	05705244002135	10	120
721060-3738	neria™ multi infusion set	8 mm	60 cm + 30 cm	3	05705244002111	10	120
721060-4737	neria™ multi infusion set	12 mm	60 cm + 30 cm	4	05705244002210	10	120
721060-4731	neria™ multi infusion set	10 mm	60 cm + 30 cm	4	05705244002197	10	120
721060-4738	neria™ multi infusion set	8 mm	60 cm + 30 cm	4	05705244002173	10	120

## Product specifications:

**Insertion angle** . . . . . 90 degrees  
**Needle length** . . . . . 8 mm, 10 mm and 12 mm  
**Needle gauge** . . . . . 27 gauge  
**Tubing length in cm** . . . . . Total of 90 cm  
**Sets available as** . . . . . bi-, tri- and quad-furcated sets  
**Tubing inner diameter** . . . . . 0,385 mm  
**Tubing outer diameter** . . . . . 1,50 mm

## Material specifications:

**PVC free** . . . . . Yes  
**Made with natural rubber latex** . . . . . No  
**Raw materials or derivation of animal origin** . . . . . No  
**The product contains phthalates** . . . . . No  
**Introducer & connector needles may contain** . . ≤ 10% nickel

## Certifications:

**CE-marking** . . . . . Yes  
**Product classification** . . . . . Class IIB (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Tariff code** . . . . . HS9018390000  
**GMDN** . . . . . 35833






## Sterilization:

**Sterile** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years  
**Usage time** . . . . . Maximum 12 hours

## Drug compatibility/stability tested with<sup>1</sup>:

Immunoglobulins

## Storage and handling:

	<b>AVOID DIRECT SUNLIGHT</b> - Keep the product away from direct sunlight and any humid environment
	<b>KEEP DRY</b> - Store in a dry place at room temperature
	<b>DO NOT RE-USE</b>
	<b>DO NOT USE</b> - If package is damaged
	<b>CAUTION</b> - See Instruction for Use