


# Ulmer drains

## Wound drains





### Intended use:

Ulmer drain is placed directly into surgical wounds by means of a trocar or other surgical instrument to perform wound drainage.

Ulmer drains have graduated hole sizes for equal suction effect along full drain length. Must be connected to a collection device. All drains are radiopaque and have metric markings at 20, 40 and 60 mm from the end of the perforation. The products are intended for single use only.



### Product range:

	Drain length, mm	Perforation length, mm	Inner/outer diameter of drain, mm 						
			1.6/2.7	2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3	3.5/6.0	
			CH 08	CH 10	CH 12	CH 14	CH 16	CH 18	
			REF*	REF*	REF*	REF*	REF*	REF*	
<b>Ulmer drain without trocar</b>  	500	50		1501818x	1501918x				
		60	1502018x						
		80		1502118x					
		100			1502418x				
		112						1502818x	1503018x
		116					1502618x		
	750	100			1502518x				
		112					1502918x	1503118x	
		116				1502718x			
<b>Ulmer drain with trocar</b>  	500	80		1505618x					
		100			1505818x				
		116				1506018x			
	750	80		1505718x					
		100			1505918x				
		112					1506318x	1506518x	
		116				1506118x			
<b>Ulmer drain with connector and clamp for Unoversal drainage system (without trocar)</b>  	500	80		1509518x					
		100			1509618x				
		116				1509718x			

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Material specifications:

**Contains latex?** ..... No  
**Contains DEHP?** ..... No  
**Drain** ..... Polyvinyl chloride  
**Trocar** ..... Stainless steel  
**Connector and clamp for Unoversal drainage system** ..... Polypropylene

### Certifications:

**CE-marking?** ..... Yes (CE 2797)  
**Product classification** ..... Class IIA (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device Nomenclature (GMDN)** ..... 11305

### Sterilization:

**Sterile?** ..... Yes  
**Sterilized by** ..... Ethylene oxide  
**Shelf life** ..... 5 years

# Ulmer drains

Wound drains



## Packaging specifications:

REF	Shelf box		Shipping carton			
	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m <sup>3</sup>
1501818x	100	297 x 92 x 239	1,200	391 x 304 x 743	11.0	0.088
1501918x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.088
1502018x	100	297 x 92 x 239	1,200	391 x 304 x 743	9.1	0.088
1502118x	100	297 x 92 x 239	1200	391 x 304 x 743	11.0	0.088
1502418x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.088
1502518x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.088
1502618x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.088
1502718x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.088
1502818x	50	297 x 92 x 239	600	391 x 304 x 743	8.5	0.088
1502918x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.088
1503018x	50	297 x 92 x 239	600	391 x 304 x 743	9.1	0.088
1503118x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.088
1505618x	100	297 x 190 x 366	200	392 x 306 x 375	5.8	0.045
1505718x	100	297 x 190 x 366	200	392 x 306 x 375	5.8	0.045
1505818x	100	297 x 190 x 366	200	392 x 306 x 375	6.4	0.045
1505918x	100	297 x 190 x 366	200	392 x 306 x 375	6.4	0.045
1506018x	100	297 x 190 x 366	200	392 x 306 x 375	8.0	0.045
1506118x	100	297 x 190 x 366	200	392 x 306 x 375	8.0	0.045
1506318x	75	297 x 190 x 366	150	392 x 306 x 375	6.4	0.045
1506518x	75	297 x 190 x 366	150	392 x 306 x 375	6.9	0.045
1509518x	100	297 x 190 x 366	400	391 x 304 x 743	8.1	0.092
1509618x	100	297 x 190 x 366	400	391 x 304 x 743	9.0	0.092
1509718x	100	297 x 190 x 366	400	391 x 304 x 743	10.1	0.092

## Storage and handling:



**ConvaTec Limited**

First Avenue  
Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

Made in Belarus



**Unomedical A/S**

Aaholmvej 1-3, Osted, 4320 Lejre  
Denmark



**ConvaTec**

ConvaTec Limited  
First Avenue, Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

[www.convatec.com](http://www.convatec.com)

**For further information on these products, please contact your local ConvaTec sales representative.**

*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

*All trade marks and logos are the property of a ConvaTec group company.*

© 2020 ConvaTec Inc. AP-030921-EM-NP-SP, ver. DEC 2020