

We listen to your clinical needs

For inspired improvements
in faecal management

Sampling port

- Provides access to the catheter for safe and easy stool collection
- Designed to help decrease risk of exposure to infectious stool



Medication delivery

- Equipped to deliver prescribed medication rectally via the irrigation port
- Simple delivery process without complicating the use of the device



Extended catheter


- Adapts to a range of bed sizes
- Accommodates a variety of patient heights and weights

Adhesive labels

- 5 adhesive labels with date printing format, available in kit pack

Featuring the unique SIGNAL™ indicator

- Aids in assuring a comfortable and customised seal
- Indicates maximum fill capacity of soft retention balloon, at any volume ≤ 45 mL, based on size of the rectal vault*



*Maximum fill volume is ≤ 45 mL. Signal indicator could pop before 45 mL if the rectal space is smaller than maximum balloon size. Filling should stop when the indicator pops out and stays out. Never add more than 45 mL total.



Flexi-Seal™ FMS

Used by more nurses and hospitals than any other faecal management system¹

Engineered by a company with a heritage of over 30 years' experience in wound, ostomy, and skin care



Flexi-Seal™ SIGNAL™ FMS Product Information

Flexi-Seal™ SIGNAL™ FMS Kit (1 kit/box, 3 bags)	418000
Charcoal Filter Collection Bag (10/box)	411102



Contact your ConvaTec representative today to learn more about **Flexi-Seal™ SIGNAL™ FMS** improvements inspired by you
ConvaTec Customer Service • 1800 339 412



Flexi-Seal™ SIGNAL™

Reference: 1. HPIS. Q4 2011 data. Moving Annual Total (MAT) Fecal Control Category. Data on file, ConvaTec.

ConvaTec (Australia) Pty Limited. ABN 70 131 232 570. Unipark Monash, Building 2, Ground Floor, 195 Wellington Road, Clayton VIC 3168 Australia.
PO Box 63, Mulgrave, VIC 3170. Phone: (03) 9239 2700 Facsimile: (03) 9239 2743. Customer Support Freecall: 1800 339 412. www.convatec.com

®/TM indicates trademarks of ConvaTec Inc. ConvaTec, the ConvaTec logo, Flexi-Seal, the Flexi-Seal logo, and Flexi-Seal SIGNAL are registered trademarks in the U.S. © 2012 ConvaTec Inc. AP-012558-MM. July 2012. CCC143