


Nasal Oxygen Cannulas

Sterile

Intended use:

The nasal cannula is a two-pronged device intended for direct low flow administration of oxygen via the nasopharyngeal route. This product is intended for single use only.



List of product range:

REF*	Prongs type	Compress	Length of prongs, mm	Distance between prongs, mm	Total length, mm	Size	Colour	Connector
2900118x	Straight	No	13	14.8	1,500	Adult	Green	Oxygen Connector Male cone connector
2900218x	Straight long	No	17	14.8	1,500	Adult	Clear	Oxygen Connector Male cone connector
2900818x	Straight	No	13	14.8	1,500	Adult	Green	Oxygen Funnel CH 14
2900918x	Straight	Yes	13	14.8	1,500	Adult	Green	Oxygen Funnel CH 14
2901218x	Curved	No	7	16	1,500	Adult	Green	Oxygen Connector Male cone connector
2901318x	Curved	No	7	16	1,500	Adult	Blue	Oxygen Connector Male cone connector
2901418x	Curved	No	7	16	1,500	Adult	Green	Oxygen Funnel CH 14 with cap
2901518x	Curved	No	7	16	2,100	Adult	Clear	Oxygen Funnel CH 10
2901718x	Straight	No	13	14.8	2,100	Adult	Clear	Oxygen Funnel CH 10
2901818x	Curved	No	7	16	2,100	Adult	Green	Oxygen Funnel CH 10
2902018x	Curved	No	7	16	2,100	Adult	Blue	Oxygen Funnel CH 10

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Oxygen connector Polyvinyl chloride
Compress Polyurethane
Tube Piece Polyethylene
Other Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification IIA rule 5 (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN).35201

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life5 years

Nasal Oxygen Cannulas

Sterile

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
2900118x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2900218x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2900818x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2900918x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2901218x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2901318x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2901418x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2901518x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2901718x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2901818x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300
2902018x	50	294 x 125 x 360	300	391 x 304 x 743	10.9	0.092	300

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

