


Drain-Fix™

Securement device

Intended use:

Drain-Fix™ securement devices are intended to be used to secure drains, catheters and tubes ranging in sizes from CH5 to CH22. Drain-Fix™ securement devices are specially designed to stop migration, movement and accidental removal of the drain/catheter. The product is intended for single use only.



List of product range:

	REF*	Size, cm	Applicable drains, CH
Drain-Fix™, Small	680M	9.1 x 10.0	5-14
Drain-Fix™, Large	685M	9.1 x 10.0	12-22

Material specifications:

Contains latex? No
Contains DEHP? No
Carrier Non-woven polypropylene
Adhesive Acrylic adhesive
Adhesive on central part Drain-Fix™ .. Hydrocolloid 0.6 mm
Plastic component Thermoplastic elastomer
Tapes Double sided adhesive tape
Backing paper on lid Paper (siliconized on one side)
Backing paper on skin contact area Paper SRP
Pouch Gas-paper, Polyester/Polyethylene, Self-reinforcing polyethylene film

Certification:

CE-marking? Yes (0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 56630

Sterilisation:

Sterile? Yes
Sterilized by Gamma radiation
Shelf life 3 years

Packaging:

	REF*	Shelf box		Shipping carton			MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Volume, m³	
Drain-Fix™, Small	680M	25	207 x 138 x 135	200	564 x 220 x 290	0.036	2,500
Drain-Fix™, Large	685M	20	207 x 138 x 135	160	564 x 220 x 290	0.036	2,000

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.