


# Silicone Elastomer coated Foley catheter

## Standard, 2-way, hard valve

**Intended use:**

The Silicone Elastomer coated Foley catheter is a urinary catheter intended for short to long term drainage of the bladder. For routine and post-operative drainage.  
The product is intended for single use only.



**List of product range:**

Silicone Elastomer coated Foley catheter, Standard, 2-way, hard valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM51161210	MM51161230
14	MM51161410	MM51161430
16	MM51161610	MM51161630
18	MM51161810	MM51161830
20	MM51162010	MM51162030
22	MM51162210	MM51162230
24	MM51162410	MM51162430
26	MM51162610	MM51162630
28	MM51162810	MM51162830
30	MM51163010	MM51163030

**Material specifications:**

Contains latex? . . . . . Yes  
 Contains DEHP? . . . . . Yes  
 Catheter . . . . . Latex  
 Hard valve . . . . . Polypropylene  
 Sleeve . . . . . Polyvinyl chloride  
 Catheter surface . . . . . Silicone Elastomer coating

**Certification:**

CE-marking? . . . . . Yes (CE 0123)  
 Product classification . . . . . Class IIB (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . .34917

**Sterilisation:**

Sterile? . . . . . Yes  
 Sterilized by . . . . . Gamma radiation  
 Shelf life . . . . . 5 years

# Silicone Elastomer coated Foley catheter

Standard, 2-way, hard valve

## Packaging:

Silicone Elastomer coated Foley catheter, Standard, 2-way, hard valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
MM51161210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51161230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51161410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51161430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51161610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51161630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51161810	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51161830	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162010	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162030	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162810	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51162830	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51163010	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51163030	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Sdn Bhd**  
(15290-P)  
Address: Bakar Arang Industrial Estate  
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**  
Address: Kongevejen 2  
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

