


Feeding Tubes

Intended use:

Feeding tubes are intended for the administration of medicine, liquids and/or nutrition (enteral feeding). A feeding tube is inserted through the nose or mouth via the oesophagus into the stomach using a syringe or tubing. The product is intended for single use only.



List of product range:

	Effective length mm						
		Tube, Outer diameter, mm					
		1.40	1.67	2.00	2.67	3.33	4.00
		CH04 REF*	CH05 REF*	CH06 REF*	CH08 REF*	CH10 REF*	CH12 REF*
	375	1202118x	1202418x	1202718x	1203018x	1203318x	
Metric with x-ray line	475	1202218x	1202518x	1202818x	1203118x	1203418x	
	725	1211518x	1211618x	1211718x	1211818x	1211918x	
Metric with x-ray line marking 400 mm	975		1211018x	1211118x	1211218x	1211318x	1211418x

All REF: Shore A hardness : 78

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Tube Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 14221

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Feeding Tubes

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
1202118x	100	145 x 74 x 576	1,000	391 x 304 x 585	5.3	0.070	2,000
1202218x	100	145 x 74 x 734	1,000	391 x 304 x 743	6.7	0.088	2,000
1202418x	100	145 x 74 x 576	1,000	391 x 304 x 585	5.8	0.070	2,000
1202518x	100	145 x 74 x 734	1,000	391 x 304 x 743	7.0	0.088	2,000
1202718x	100	145 x 74 x 576	1,000	391 x 304 x 585	6.3	0.070	2,000
1202818x	100	145 x 74 x 734	1,000	391 x 304 x 743	7.2	0.088	2,000
1203018x	100	145 x 74 x 576	1,000	391 x 304 x 585	6.7	0.070	2,000
1203118x	100	145 x 74 x 734	1,000	391 x 304 x 743	8.3	0.088	2,000
1203318x	100	145 x 74 x 576	1,000	391 x 304 x 585	6.7	0.070	2,000
1203418x	100	145 x 74 x 734	1,000	391 x 304 x 743	8.3	0.088	2,000
1211018x	100	294 x 125 x 360	600	391 x 304 x 743	5.4	0.088	1,800
1211118x	100	294 x 125 x 360	600	391 x 304 x 743	5.7	0.088	1,800
1211218x	100	294 x 125 x 360	600	391 x 304 x 743	5.9	0.088	1,800
1211318x	100	294 x 125 x 360	600	391 x 304 x 743	6.9	0.088	1,800
1211418x	100	294 x 125 x 360	600	391 x 304 x 743	5.5	0.088	1,200
1211518x	100	294 x 125 x 360	600	391 x 304 x 743	5.2	0.088	1,200
1211618x	100	294 x 125 x 360	600	391 x 304 x 743	5.2	0.088	1,200
1211718x	100	294 x 125 x 360	600	391 x 304 x 743	5.2	0.088	1,200
1211818x	100	294 x 125 x 360	600	391 x 304 x 743	5.5	0.088	4,200
1211918x	100	294 x 125 x 360	600	391 x 304 x 743	6.2	0.088	1,200

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

