

# DuoDERM® *Extra Thin*

Designed to support wound repair and to protect skin...  
**here, there, everywhere!**



DuoDERM® Extra Thin dressing is designed to cover and protect susceptible areas

- Can be used anywhere on the body to gently protect newly-formed tissue or skin at risk of breakdown
- Designed to reduce cost of care

Supports moist wound healing

- Provides a moist wound healing environment that supports epithelialization<sup>1</sup>
- May be used on skin tears, superficial wounds, dry-to-lightly exuding dermal ulcers, and postoperative wounds

# DuoDERM<sup>®</sup>

## Extra Thin



DuoDERM<sup>®</sup> Extra Thin dressings offer so many product benefits... whether you want to support wound repair or to protect skin.

### Designed for patient comfort

- Translucent backing is designed to prevent needless dressing changes that can cause discomfort and disturb healing
- Thin design is flexible and conformable
- Designed to promote a moist wound environment, which supports healing and prevents damage to new epithelium upon removal
- Moisture-resistant for ease of cleaning/showering

### Incredibly versatile

- Can be used as a primary or secondary dressing
- Available in a wide variety of sizes and shapes—including oval, square, and rectangle
- Can be cut to shape

### Creates a moist wound healing environment

- DuoDERM is the only hydrocolloid dressing with both a Triple Hydrocolloid Matrix<sup>®</sup> and a viral/bacterial barrier<sup>\*2-4</sup>


### Helps to reduce cost of care for caregiver

- Preventing breakdown may spare cost of treating wounds that may have developed
- Long wear time means fewer dressings required
- DuoDERM Extra Thin dressing is a low-cost hydrocolloid

DuoDERM Extra Thin dressings come in the following sizes and shapes



Unit Size	Product Number	Quantity Per Box	HCPCS Code
<b>Spot</b>			
1.25" x 1.5"	187932	20	A6234
<b>Rectangle</b>			
2" x 4"	187900	20	A6234
2" x 8"	187961	10	A6234
<b>Square</b>			
3" x 3"	187901	20	A6234
4" x 4"	187955	10	A6234
6" x 6"	187957	10	A6235
<b>Oval</b>			
4" x 6"	187902	10	A6235
<b>Triangle</b>			
6" x 7"	187903	10	A6235



Any wound larger than 1.5" needs a dressing larger than 4"x4"!

References: 1. Greguric S, Budimcic D, Soldo-Belic A, et al. Hydrocolloid dressing versus a conventional dressing using magnesium sulphate paste in the management of venous leg ulcers. *Acta Dermatovenerol Croat.* 1994;2:65-71. 2. Wilson P, Burroughs D, Dunn LJ. Methicillin-resistant *Staphylococcus aureus* and hydrocolloid dressings. *Pharm J.* 1988;243:787-788. 3. Mertz PM, Marshall DA, Eaglstein WH. Occlusive wound dressings to prevent bacterial invasion and wound infection. *J Am Acad Dermatol.* 1985;12:662-668. 4. Bowler PG, Delargy H, Prince D, Fondberg L. The viral barrier properties of some occlusive dressings and their role in infection control. *Wounds.* 1993;5:1-8.

Please see package insert for full product information including instructions for use. \*While intact without leakage.

DuoDERM is a registered trademark of E.R. Squibb & Sons, L.L.C. ©2005 E.R. Squibb & Sons, L.L.C. Printed in U.S.A.

US-05-464

To learn more, call  
**1-800-422-8811**  
 M-Th, 8:30 AM–8:00 PM, EST  
 Fri, 8:30 AM–6:00 PM, EST  
[www.convatec.com](http://www.convatec.com)

 **Convatec**  
 A Bristol-Myers Squibb Company