

# A4

## Sterile urine bags

### Intended use:

The A4 sterile urine bag is intended for the collection of urine. The product is intended for single use only. Fitted with bottom outlet.



### Product range:

REF*	Tube length cm	Standard sample port (for use with needle)	Kombikon™ needle-free sample port	Integrated hanger	Bottom outlet / Cross valve
38 103 1S-xx	10	No	no	No	Yes
38 117 1S-xx	110	Yes	no	No	Yes
38 151 1S-xx	120	Yes	No	No	Yes
38 155 1S-xx	120	No	Yes	No	Yes
38 161 1S-xx	110	Yes	No	Yes	Yes
38 162 1S-xx	120	No	Yes	Yes	Yes

\* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Material specifications:

**Contains latex?** . . . . . No  
**Contains DEHP?** . . . . . No  
**Film** . . . . . Polyvinyl chloride  
**Tube** . . . . . Polyvinyl chloride  
**Standard sample port** . . . . . Polyvinyl chloride  
**Kombikon™ needle-free sample port** . . . Polyvinyl chloride, Polypropylene, High density polyethylene, polyethylene, Poly(styrene-butadiene-styrene)  
**Cap** . . . Low density polyethylene, High density polyethylene  
**Hanger** . . . . . Polypropylene, High density polyethylene

### Certifications:

**CE-marking?** . . . . . Yes (CE 0086)  
**Product classification** . . . . . Class I (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)** . . . . . 58921

### Sterilization:

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

### Packaging:

REF*	Shipping Carton				MOQ, pcs.
	Qty., pcs.	Dimensions (LxWxD), mm	Gross Wt., kg	Volume, m <sup>3</sup>	
38 103 1S-xx	80	375 x 280 x 278	4.8	0.03	3840
38 117 1S-xx	40	375 x 280 x 278	9.6	0.03	1920
38 151 1S-xx	50	375 x 280 x 278	5.7	0.03	2400
38 155 1S-xx	50	375 x 280 x 278	5.7	0.03	2400
38 161 1S-xx	30	388 x 286 x 242	3.8	0.03	1440
38 162 1S-xx	30	388 x 286 x 242	3.8	0.03	1440

\* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

# A4

## Sterile urine bags



### Storage and handling:



**ConvaTec Limited**

First Avenue  
Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK



For further information on these products, please contact your local ConvaTec sales representative.  
*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

All trade marks and logos are the property of a ConvaTec group company.  
© 2016 ConvaTec Inc. AP-xxxxx-MM