

Product Data Sheet

Neria™Soft

Intended use

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency, apomorphine for Parkinson's Disease and morphine for pain management therapy. The Neria™Soft infusion set has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion.



Product specifications

Product code	Description	Cannula length	Tubing length	GTIN (Inner box)	Pcs per inner box	Pcs per outer box
79-030-2522	Neria™Soft infusion set	13 mm	30 cm	05705244000339	10	200
79-060-2522	Neria™Soft infusion set	13 mm	60 cm	05705244000353	10	200
79-080-2522	Neria™Soft infusion set	13 mm	80 cm	05705244000377	10	200
79-110-2522	Neria™Soft infusion set	13 mm	110 cm	05705244000391	10	200
79-060-2622	Neria™Soft infusion set	17 mm	60 cm	05705244000315	10	200
79-080-2622	Neria™Soft infusion set	17 mm	80 cm	05705244001558	10	200
79-110-2622	Neria™Soft infusion set	17 mm	110 cm	05705244001138	10	200

For further information on these products, please contact your local Convatec sales representative.

Neria™ Soft

Product specifications:

Insertion angle	20-45 degrees
Soft cannula length	13 mm and 17 mm
Soft cannula outer diameter	~ 0,68 mm
Introducer needle	27 gauge
Tubing length in cm	30, 60, 80 and 110 cm
Tubing inner diameter	0,41 mm
Tubing outer diameter	1,50 mm

Material specifications:

PVC free	Yes
Made with natural rubber latex	No
Raw materials or derivation of animal origin	No
The product contains phthalates	No

Certifications:

CE-marking	Yes
Product classification	Class IIB (EU)
Quality system certified to	ISO 13485
Tariff code	HS9018390000
GMDN	35833






Sterilization:

Sterile	Yes
Sterilized by	Ethylene oxide
Shelf life	5 years
Usage time (drug depended)	Maximum 72 hours

Drug compatibility/stability tested with¹:

- Immunoglobulins
- Apomorphine
- Morphine
(hydromorphone, morphine sulphate & morphine chloride)

Storage and handling

	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	DO NOT RE-USE
	DO NOT USE IF PACKAGE IS DAMAGED
	CAUTION - See Instruction for Use



Manufacturer:
Unomedical A/S - A Convatec company
 Aaholmvej 1-3, Østød,
 4320 Iejre, Denmark
www.convatec.com
 Made in Mexico

For further information on these products, please contact your local Convatec sales representative.
Please refer to the instructions for use for additional information, warnings, and precautions.
 With reservations for misprints and alterations in product specifications. Printed version not controlled.

Reference: 1. Convatec IC, data on file, 3A02300 (21) Product Specification neria™ soft, 2024.