

Aerosol Masks

Intended use:

Aerosol masks are intended for nebulization via a mask in order to administer medication via respiration. The product is intended for single use only.


Product range:

	REF*	Item	Size	Description
	93-100MM 93-100MMV (Vanilla scented mask)	Aerosol Mask	Adult	Adult Aerosol Mask (Elongated) with aerosol connector, Adjustable nose-clip and elastic strap
	93-220MM 93-220MMV (Vanilla scented mask)	Aerosol Mask	Pediatric	Pediatric Aerosol Mask (Elongated) with aerosol connector, Adjustable nose-clip and elastic strap
	93-260MM	Aerosol Mask	Infant	Infant Aerosol Mask with aerosol connector, Adjustable nose-clip and elastic strap
	93-110MM	Face Tent	Adult	Adult Face Tent with aerosol connector and elastic strap

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Contains bisphenols? No
Aerosol mask Polyvinyl chloride
Nose clip Aluminum
Elastic strap Spandex/Polyamide (Latex-free)
Aerosol connector Polypropylene

Certifications:

CE-marking? Yes (CE 1639)
Product classification Class IIA (EU)
Quality system certified to ISO 13485

Sterilization:

Sterile? No
Shelf life 5 years

Distributed by


ConvaTec

Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
93-100MM	1	10	50	400 x 300 x 240	1.91	1,600
93-100MMV	1	10	50	400 x 300 x 240	2.03	1,600
93-220MM	1	10	50	350 x 300 x 200	1.48	1,600
93-220MMV	1	10	50	350 x 300 x 200	1.28	1,600
93-260MM	1	10	50	350 x 300 x 200	1.40	1,600
93-110MM	1	10	50	350 x 300 x 200	1.48	900

EAN codes:

REF	Primary package	Secondary package	Tertiary package
93-100MM	(01)06945099302052	(01)26945099302056	(01)46945099302050
93-100MMV	(01)06945099302113	(01)26945099302117	(01)46945099302111
93-220MM	(01)06945099302069	(01)26945099302063	(01)46945099302067
93-220MMV	(01)06945099302120	(01)26945099302124	(01)46945099302128
93-260MM	(01)06945099302076	(01)26945099302070	(01)46945099302074
93-110MM	(01)06945099302106	(01)26945099302100	(01)46945099302104

Packaging material:

Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK
critcare.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

© 2021 ConvaTec. ®/™ indicate trademarks of the ConvaTec group of companies.

© 2021 ConvaTec Inc. AP-017302-MM, ver. AUG 2021