












Connecting Tubes

Intended use:

Connecting tubes are intended for various medical procedures, such as suction procedures for the removal of blood and fluids, primarily in hospital departments such as the OR and ICU. The product is intended for single use only.



List of product range:

Connecting tubes		REF*	Length, mm	CH/FG	Shore A
No connector		1603918x	3,000	30	64
One female connector		1605918x	2,100	24	78
		1606418x	3,000	30	78
		1603118x	3,500	24	64
		1606618x	4,800	30	78
One male fingertip connector		1601318x	3,500	30	78
Two female connectors		1604918x	2,100	24	78
		1605018x	3,000	24	78
		1604218x	4,000	24	64
		1605618x	4,000	30	78
One female and one male connector		1617218x	3,000	25	78
		1617318x	5,000	24	78
One female connector and one male fingertip connector		1602118x	2,100	24	78
One female connector (soft funnel) and one male fingertip connector		1602418x	2,100	24	78
One female connector and one nipple connector		1606918x	2,100	24	78
One female connector and one tube piece		1607418x	4,000	24	78
Two Luer connectors (one male and one female)		1617518x	1,460	6	78

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Connecting Tubes

Material specifications:

Contains latex? No
Contains DEHP? Yes
Tube Polyvinyl chloride
Female connectors Polyvinyl chloride
Male connectors Polypropylene
Male fingertip connectors Polyvinyl chloride
Nipple connectors Polyvinyl chloride
Cap for nipple connector Polyethylene
Tube piece Polyvinyl chloride
Luer connectors Methyl methacrylate-acrylonitrile-butadienestyrene (MABS) polymer, Polypropylene

Certification:

CE-marking? Yes (CE 0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 16779

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
1601318x	*	*	50	392 x 306 x 375	7.80	0.045	150
1602118x	*	*	50	392 x 306 x 340	5.50	0.043	150
1602418x	*	*	50	392 x 306 x 340	5.50	0.043	150
1603918x	*	*	50	392 x 306 x 375	8.60	0.045	150
1604218x	*	*	40	392 x 306 x 375	7.90	0.045	120
1604918x	*	*	50	392 x 306 x 340	5.30	0.043	300
1605018x	*	*	50	392 x 306 x 340	6.50	0.045	150
1605618x	*	*	40	392 x 306 x 375	9.00	0.045	120
1605918x	*	*	50	392 x 306 x 340	5.50	0.043	150
1606418x	*	*	50	392 x 306 x 375	8.90	0.045	150
1606618x	*	*	40	391 x 304 x 455	6.50	0.049	120
1606918x	*	*	50	392 x 306 x 340	5.50	0.043	300
1617218x	*	*	50	392 x 306 x 375	5.50	0.045	300
1617318x	*	*	40	392 x 306 x 375	5.50	0.045	280
1607418x	*	*	40	392 x 306 x 375	6.50	0.045	120
1617518x	100	297 x 190 x 366	200	392 x 306 x 375	7.40	0.045	200

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: ConvaTec Limited
 First Avenue
Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

