

# Thyroidea set

## Wound drainage system

### Intended use:

Low vacuum wound drainage system for post operative removal of blood and body fluids after thyroidectomy. The product is intended for single use only.



### Product range:

	REF*	Drains Shore A Hardness 69			
		CH/FG	Length, mm	I.D. / O.D., mm	Perforation Length, mm
Thyroidea Set with 100 ml bellows, 2 Ulmer drains with trocar, 150 mm connecting tube with clamp and connector for optional use of 1 or 2 drains	1300418x	10	500	2.0 / 3.3	50

\* Here noted as x. reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Material specifications:

Contains latex? . . . . . No  
 Contains DEHP? . . . . . No  
 Bellows . . . . . Polyethylene  
 Drains and connecting tube . . . . . Polyvinyl chloride  
 Trocar . . . . . Stainless steel  
 Connector . . . . . Silicone  
 Tube clamp . . . . . Polypropylene

### Certifications:

CE-marking? . . . . . Yes (CE 2797)  
 Product classification . . . . . Class IIA (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN) . . . . . 35824

### Sterilization:

Sterile? . . . . . Yes  
 Sterilized by . . . . . Ethylene oxide  
 Shelf life . . . . . 5 years

### Packaging specifications:

REF	Shelf box		Shipping carton			
	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m <sup>3</sup>
1300418x	10	294 x 125 x 360	60	391 x 304 x 743	5.50	0.088

### Storage and handling:




**ConvaTec Limited**  
 First Avenue  
 Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK  
 Made in Belarus

EC REP  
**Unomedical A/S**  
 Aaholmvej 1-3, Osted, 4320 Lejre  
 Denmark



**ConvaTec**  
 ConvaTec Limited  
 First Avenue, Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK  
[www.convatec.com](http://www.convatec.com)

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled. All trade marks and logos are the property of a ConvaTec group company. © 2020 ConvaTec Inc. AP-030272-EM-NP-SP, ver. AUG 2020