

ClearVue™ Plus Anesthesia Mask

PVC-Free

Intended use:

Anesthesia masks (face masks) are intended for the delivery of anesthetic gases, air and/or oxygen to patients. The products are for single use only.



Product range:

	REF*	Item	Size		Colour (open cushion)	Hook ring (removable)
			Mask	Connector		
	93-681UM	ClearVue™ Plus Anesthesia Mask PVC-Free	Size 1 Neonate	15mm (Male)	Blue	-
	93-682UM	ClearVue™ Plus Anesthesia Mask PVC-Free	Size 2 Infant	15mm (Male)	Grey	-
	93-683UM	ClearVue™ Plus Anesthesia Mask PVC-Free	Size 3 Pediatric	22mm (Female)	White	Yes
	93-684UM	ClearVue™ Plus Anesthesia Mask PVC-Free	Size 4 Pediatric/Adult (Small)	22mm (Female)	Yellow	Yes
	93-685UM	ClearVue™ Plus Anesthesia Mask PVC-Free	Size 5 Adult (Medium)	22mm (Female)	Green	Yes
	93-686UM	ClearVue™ Plus Anesthesia Mask PVC-Free	Size 6 Adult (Large)	22mm (Female)	Orange	Yes

*Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Contains bisphenols? No
Anesthesia mask Polypropylene
Cushion Thermoplastic elastomer
Connector Polypropylene
Hook ring Polycarbonate

Certifications:

CE-marking? Yes (CE 1639)
Product classification Class IIA (EU)
Quality system certified to ISO 13485

Sterilization:

Sterile? No
Shelf life 3 years



Packaging specifications:

REF	Inner package		Outer package/Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	Qty, pcs
93-681UM	1	160 x 170	20	400 x 200 x 150	0.55	0.012	960
93-682UM	1	160 x 170	20	400 x 200 x 150	0.60	0.012	960
93-683UM	1	160 x 170	50	400 x 390 x 260	1.85	0.041	1,200
93-684UM	1	180 x 220	50	400 x 390 x 260	2.10	0.041	1,200
93-685UM	1	180 x 220	50	560 x 400 x 210	2.40	0.047	1,000
93-686UM	1	180 x 220	50	560 x 400 x 210	2.60	0.047	1,000





EAN codes:

REF	Inner package	Outer package
93-681UM	(01)06945099304216	(01)46945099304214
93-682UM	(01)06945099304223	(01)46945099304221
93-683UM	(01)06945099304230	(01)46945099304238
93-684UM	(01)06945099304247	(01)46945099304245
93-685UM	(01)06945099304254	(01)46945099304252
93-686UM	(01)06945099304261	(01)46945099304269

Packaging material:

Inner package	Poly bag - Low density polyethylene / Insert - 80g paper
Outer package	Carton - Corrugated board / Sticker - Coated paper

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2020 ConvaTec Inc. AP-021043-MM, ver. AUG 2020