

Aerosol Masks, Tracheostomy Masks, Face Tent for mask nebulization


Intended use:

Aerosol masks, tracheostomy masks and face tents are intended for nebulization via a mask or face tent in order to administer medication via respiration. The product is intended for single use only.




List of product range:


Aerosol masks, Elongated (under the chin) with adjustable nose-clip and elastic braid*

	REF**	Size	Description	Length of aerosol tubing, m
100MM	100MM	Adult	Aerosol mask (elongated) only	-
	118MM	Adult	Aerosol mask (elongated), with aerosol tubing	1.8
	220MM	Pediatric	Aerosol mask (elongated) only	-
220MM	221MM	Pediatric	Aerosol mask (elongated), with Aerosol tubing	1.5


Aerosol masks, Standard (under the lip) with adjustable nose-clip and elastic braid*

	REF**	Size	Description	Length of aerosol tubing, m
260MM	200MM	Adult	Aerosol mask (Standard) only	-
	260MM	Infant	Aerosol mask (Standard) only	-

Trach masks with adjustable elastic braid*

	REF**	Size	Description	Length of aerosol tubing, m
245MM	245MM	Adult	Trach mask with swivel connector	-
246MM	246MM	Pediatric	Trach mask with swivel connector	-

Face tents

	REF**	Size	Description	Length of aerosol tubing, m
110MM	110MM	Adult	Face tent	-

* Connectors accept standard 22 mm aerosol tubing

** Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Aerosol mask Polyvinyl chloride
 Trach mask Polyvinyl chloride
 Face tent Polyvinyl chloride
 Nose clip Aluminum
 Elastic braid Polyester (Latex-Free)
 Aerosol tubing Ethylene vinyl acetate
 Aerosol connector High density polyethylene

Certification:

CE-marking? Yes (CE 0120)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device
 Nomenclature (GMDN) 35142/35179

Sterilisation:

Sterile? No
 Shelf life 5 years

Aerosol Masks, Tracheostomy Masks, Face Tent for mask nebulization

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
100MM	*	*	50	381 x 254 x 203	2.04	0.019	1,250
110MM	*	*	50	406 x 355 x 254	4.54	0.037	1,250
118MM	*	*	50	406 x 609 x 584	7.26	0.144	1,250
200MM	*	*	50	381 x 254 x 203	1.81	0.019	1,250
220MM	*	*	50	381 x 254 x 203	1.81	0.019	1,250
221MM	*	*	50	406 x 609 x 584	4.99	1.440	1,250
245MM	*	*	50	381 x 254 x 203	2.04	0.019	1,250
246MM	*	*	50	381 x 254 x 203	2.04	0.019	1,250
260MM	*	*	50	381 x 254 x 203	1.36	0.019	1,250

* = Single packed - no shelf box

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Inc.**
5701-1S, Ware Road
Address: McAllen, Texas 78503
USA

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

