


Muco-Safe™











Mucus Extractor

Intended use:

Muco-Safe™ Mucus Extractor is intended to be used for aspiration of secretion from the oro-pharynx in newborn babies to ensure free respiration or to obtain mucus specimen for microbiological examination. The products are intended for single use only.






List of product range:

Muco-Safe™/ Mucus Extractors	REF*	Description	Suction tube					Operator tube
			Shore A hardness	CH	Outer diameter mm	Length mm	Connector	Connector
 	1702818x	Muco-Safe™ with hydrophobic filter and extra lid	78	08	2.7	370	-	Funnel
	1702918x		78	10	3.3	370	-	Funnel
	1703018x		78	12	4.0	370	-	Funnel
	1703118x		78	14	4.7	370	-	Funnel
 	1703618x	Muco-Safe™ with hydrophobic filter	64	14	4.7	370	-	Funnel
 	1703218x	Muco-Safe™ with hydrophobic filter, extra lid and male connector for attachment of funnel suction catheters	78	-	-	-	Male	Funnel
	1700518x	Mucus Extractor (without filter) with extra lid	78	10	3.3	370	-	Funnel
	1700618x		78	12	4.0	370	-	Funnel
	1700818x		78	14	4.7	370	-	Funnel
	1700718x		78	12	4.0	530	-	Funnel
	1700918x		78	14	4.7	530	-	Funnel
	1702318x	Mucus Extractor (without filter)	78	10	3.3	370	-	Funnel
	1702418x		78	12	4.0	370	-	Funnel
	1702518x		78	14	4.7	370	-	Funnel
	1702018x	Mucus Extractor (without filter) with extra lid and male connector for attachment of funnel suction catheters, includes one funnel suction catheter	80	6	2.0	370	Male	Funnel
	1702118x	Mucus Extractor (without filter) with extra lid and male connector for attachment of funnel suction catheters	78	-	-	-	Male	Funnel

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

List of product range:

Muco-Safe™/ Mucus Extractors	REF*	Description	Suction tube					Operator tube
			Shore A hardness	CH	Outer diameter mm	Length mm	Connector	Connector
	1700118x	Double chamber Mucus-Extractor (without filter)	64	10	3.3	370	-	Funnel
	1700218x		64	12	4.0	370	-	Funnel
	1700418x		64	14	4.7	370	-	Funnel
	1702218x		78	10	3.3	370	-	Funnel
	258100210120	Mucus Extractor (without filter) includes one extra lid, with interchangeable Mully catheters	78	14	4.7	430	Male	Funnel
			78	16	5.3			
	258100110120	Mucus Extractor (without filter) includes extra lid, with interchangeable Mully catheters	80	6	2.0	430	Male	Funnel
			78	10	3.3			
			78	14	4.7			

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Chamber. Polypropylene
Filter Polyester, Polyethylene
Lid Polyvinyl chloride
Connectors Polyvinyl chloride
Identification label. Paper

Certification:

CE-marking? Yes (CE 0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN).38472

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume m ³	
1700118x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	200
1700218x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	200
1700418x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	200
1700518x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1700618x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1700718x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1700818x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1700918x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1702018x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1702118x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1702218x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	200
1702318x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1702418x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1702518x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1702818x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1702918x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1703018x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1703118x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1703218x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
1703618x	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	600
258100210120	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	8,000
258100110120	50	297 x 190 x 366	200	391 x 304 x 743	7	0.092	8,000

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

