

To whom it may concern

Subject: Temporary Product Data Sheet for Suction Catheter 84003181

I herewith confirm the below product specifications for the Unomedical reference code 84003181

Product Specifications:

REF	Description
84003181	Suction catheter with fingertip (with strap and cap)

Packaging Information:

REF	Shelf Box			Shipping Carton		
	Qty. pcs	Dimensions, mm	Volume, m ³	Qty. pcs	Dimensions, mm	Volume, m ³
84003181	100	297×92×239	0.0065	1200	391×304×743	0.0883

Material Specifications:

Contains Latex?	No
Contains DEHP?	Yes
Suction catheter with fingertip	Polyvinyl chloride

Certifications:

CE-marking?	Yes (CE 0086)
Product Classification	Ila
Quality system certified to	ISO 13485
Global Medical Device Nomenclature (GMDN)	10749

Sterilization:

Sterile?	No
Sterilization Method	Ethylene Oxide
Shelf life	5 years

Packaging:



Manufacturer:

ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

*TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.*

Schaffhausen, May 12, 2016



Antje Wedler, Product Manager ConvaTec Continence & Critical Care