


neria™ detach

Intended use:

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency, apomorphine for Parkinson's Disease and morphine for pain management therapy. The neria™ detach infusion set has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion.


Product range:

Product code	Description	Needle length	Tubing length	GTIN (Inner box)	Pcs per inner box	Pcs per outer box
neria™ detach						
78-060-2721	neria™ detach infusion set	10 mm	60 cm	05705244000872	10	200
78-060-2728	neria™ detach infusion set	8 mm	60 cm	05705244000858	10	200
78-060-2726	neria™ detach infusion set	6 mm	60 cm	05705244000834	10	200
78-080-2721	neria™ detach infusion set	10 mm	80 cm	05705244000933	10	200
78-080-2728	neria™ detach infusion set	8 mm	80 cm	05705244000919	10	200
78-080-2726	neria™ detach infusion set	6 mm	80 cm	05705244000896	10	200
78-110-2721	neria™ detach infusion set	10 mm	110 cm	05705244000995	10	200
78-110-2728	neria™ detach infusion set	8 mm	110 cm	05705244000971	10	200
78-110-2726	neria™ detach infusion set	6 mm	110 cm	05705244000957	10	200

Product specifications:

Insertion angle 90 degrees
Needle length 6 mm, 8mm and 10mm
Needle gauge 27 gauge
Tubing length in cm 60, 80 and 110
Tubing inner diameter 0,385 mm
Tubing inner outer diameter 1,50 mm

Material specifications:

PVC free Yes
Made with natural rubber latex No
Raw materials or derivation of animal origin No
The product contains phthalates No
Introducer & connector needles may contain . . ≤ 10% nickel

Certifications:

CE-marking Yes
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Tariff code HS9018390000
GMDN 35833

Sterilization:

Sterile Yes
Sterilized by Ethylene oxide
Shelf life 5 years
Usage time (drug depended) Maximum 48 hours

Drug compatibility/stability tested with¹:

Immunoglobulins
Apomorphine
Morphine ((hydromorphone, morphine sulphate & morphine chloride)

Storage and handling:

	AVOID DIRECT SUNLIGHT - Keep the product away from direct sunlight and any humid environment
	KEEP DRY - Store in a dry place at room temperature
	DO NOT RE-USE
	DO NOT USE - If package is damaged
	CAUTION - See Instruction for Use



Unomedical A/S A ConvaTec Company

Aaholmvej 1-3
 Østved, 4320 Lejre, Denmark
www.neria.com