


Epi-Fix™

Securement device

Intended use:

Epi-Fix™ securement device is a sterile device for securing epidural catheters to prevent kinking, migration and infections. The product is intended for single use only.



List of product range:

	REF*	Size, cm
Epi-Fix™	670M	10.4 x 8.7

Material specifications:

Contains latex? No
Contains DEHP? No
Carrier Non-woven polypropylene
Carrier adhesive Skin-compatible medical grade acrylic
Film Polyurethane
Permeability Min. 1500g/m²/24h/37°C
Adhesive in film Medical grade acrylic
White foam Alveolit
Backing paper Paper with Polyethylene

Certification:

CE-marking? Yes (0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 34916

Sterilisation:

Sterile? Yes
Sterilized by Gamma radiation
Shelf life 2 years

Packaging:

	REF*	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
Epi-Fix™	670M	25	102 x 87 x 135	250	400 x 296 x 230	2.52	0.0272	25

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.