

Easl-V™

Securement device

Intended use:

Easl-V™ securement device is designed to secure the peripheral intravenous line firmly in place. The product offers a transparent window for inspection of the puncture site. The product is intended for single use only.



List of product range:

	REF*	Size, cm
Easl-V™, Foam	664M	5.8 x 7.6
Easl-V™, Non-woven	665M	5.8 x 7.6

Material specifications:

Contains latex? No
Contains DEHP? No
Foam border and strips . Polypropylene non-woven (micropore)
 Polyethylene foam (foam)
Adhesive on foam Medical grade acrylic
Film Transparent polyurethane
Permeability Min. 1500g/m²/24h/37°C
Backing liner Siliconized polyester
Backing paper Siliconized paper
Adhesive on film Medical grade acrylic adhesive

Certification:

CE-marking? Yes (0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN)56631

Sterilisation:

Sterile? Yes
Sterilized by Gamma radiation
Shelf life 2 years

Packaging:

	REF*	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
Easl-V™, Foam	664M	50	315 x 240 x 325	800	400 x 296 x 230	4.36	0.0272	800
Easl-V™, Non-woven	665M	50	315 x 240 x 325	800	400 x 296 x 230	4.38	0.0272	800

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.