


neria™ soft 90

Intended use:

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency, apomorphine for Parkinson's Disease and morphine for pain management therapy.

The neria™ soft 90 infusion set has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion.



Product range:

Product code	Description	Cannula length	Tubing length	GTIN (Inner box)	Pcs per inner box	Pcs per outer box
neria™ soft 90						
722030-5226	neria™ soft 90 infusion set	6 mm	30 cm	05705244003637	10	100
722030-5229	neria™ soft 90 infusion set	9 mm	30 cm	05705244003651	10	100
722060-5226	neria™ soft 90 infusion set	6 mm	60 cm	05705244003675	10	100
722060-5229	neria™ soft 90 infusion set	9 mm	60 cm	05705244003699	10	100
722110-5226	neria™ soft 90 infusion set	6 mm	110 cm	05705244003712	10	100
722110-5229	neria™ soft 90 infusion set	9 mm	110 cm	05705244003736	10	100

Product specifications:

Insertion angle 90 degrees
Soft cannula length 6 mm and 9 mm
Soft cannula outer diameter ~ 0,68 mm
Introducer needle 27 gauge
Tubing length in cm 30, 60 and 110
Tubing inner diameter 0,385 mm
Tubing outer diameter 1,50 mm

Material specifications:

PVC free Yes
Made with natural rubber latex No
Raw materials or derivation of animal origin No
The product contains phthalates No
Introducer & connector needles may contain . . ≤ 10% nickel

Certifications:

CE-marking Yes
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Tariff code HS9018390000
GMDN 35833






Sterilization:

Sterile Yes
Sterilized by Ethylene oxide
Shelf life 3 years
Usage time (drug depended) Maximum 72 hours

Drug compatibility/stability tested with¹:

Apomorphine
Immunoglobulins
Morphine (hydromorphone, morphine sulphate & morphine chloride)

Storage and handling:

	AVOID DIRECT SUNLIGHT - Keep the product away from direct sunlight and any humid environment
	KEEP DRY - Store in a dry place at room temperature
	DO NOT RE-USE
	DO NOT USE - If package is damaged
	CAUTION - See Instruction for Use