

















Tracheal Tubes

Standard










Intended use:

Tracheal Tubes are inserted through the mouth or nose into the trachea and connected to ventilation equipment to secure an airway. The products are intended for single use only.


Product range:
Magill shaped Tracheal Tubes with HVLP* cuff

| With Murphy eye REF | Without Murphy eye REF | With stylet and Murphy eye REF | Size | | Length mm* | CRD mm** | Connector colour*** |
|---------------------------|------------------------------|--------------------------------------|------------|------------|---------------|-------------|--|
| | | | I.D. mm | O.D. mm | | | |
| UM61110030 | UM61120030 | | 3.0 | 4.2 | 165 | 12 | Red  |
| UM61110035 | UM61120035 | | 3.5 | 4.9 | 185 | 12 | Red  |
| UM61110040 | UM61120040 | | 4.0 | 5.5 | 210 | 14 | Green  |
| UM61110045 | UM61120045 | | 4.5 | 6.2 | 225 | 14 | Green  |
| UM61110050 | UM61120050 | UM61110150 | 5.0 | 6.9 | 245 | 17 | Blue  |
| UM61110055 | UM61120055 | UM61110155 | 5.5 | 7.5 | 275 | 17 | Blue  |
| UM61110060 | UM61120060 | UM61110160 | 6.0 | 8.2 | 285 | 20 | Blue  |
| UM61110065 | UM61120065 | UM61110165 | 6.5 | 8.8 | 295 | 20 | Black  |
| UM61110070 | UM61120070 | UM61110170 | 7.0 | 9.6 | 310 | 25 | Black  |
| UM61110075 | UM61120075 | UM61110175 | 7.5 | 10.2 | 315 | 25 | White  |
| UM61110080 | UM61120080 | UM61110180 | 8.0 | 10.9 | 330 | 26 | White  |
| UM61110085 | UM61120085 | UM61110185 | 8.5 | 11.5 | 330 | 26 | Green  |
| UM61110090 | UM61120090 | | 9.0 | 12.1 | 330 | 28 | Green  |
| UM61110095 | UM61120095 | | 9.5 | 12.7 | 330 | 28 | Orange  |
| UM61110100 | UM61120100 | | 10.0 | 13.6 | 330 | 28 | Orange  |

Magill shaped Tracheal Tubes without HVLP* cuff

| With Murphy eye REF | Without Murphy eye REF | Size | | Length mm* | Connector colour*** |
|---------------------------|------------------------------|------------|------------|---------------|---|
| | | I.D. mm | O.D. mm | | |
| UM61130020 | UM61140020 | 2.0 | 2.9 | 140 | Clear  |
| UM61130025 | UM61140025 | 2.5 | 3.5 | 145 | Clear  |
| UM61130030 | UM61140030 | 3.0 | 4.2 | 165 | Red  |
| UM61130035 | UM61140035 | 3.5 | 4.9 | 185 | Red  |
| UM61130040 | UM61140040 | 4.0 | 5.5 | 210 | Green  |
| UM61130045 | UM61140045 | 4.5 | 6.2 | 225 | Green  |
| UM61130050 | UM61140050 | 5.0 | 6.9 | 245 | Blue  |
| UM61130055 | UM61140055 | 5.5 | 7.5 | 275 | Blue  |
| UM61130060 | UM61140060 | 6.0 | 8.2 | 285 | Blue  |

* HVLP = High Volume Low Pressure; Total length of tube excluding connector

** CRD = Cuff Resting Diameter

*** 15 mm ISO connector (ISO 8836)

Tracheal Tubes

Standard



Material specifications:

Contains latex? No
Contains phthalates? No
Tube Polyvinyl chloride
Cuff Polyvinyl chloride
X-ray line Polyvinyl chloride
Inflation line Polyvinyl chloride
Connector Polypropylene
Pilot balloon Polyvinyl chloride
Valve Polyvinyl chloride, Polypropylene,
 Acrylonitrile-butadiene-styrene, Silicone, Stainless steel
Stylet Aluminum
Sheath Polyvinyl chloride

Certifications:

CE-marking? Yes (0123)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 46967

Sterilization:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging specifications:

| REF* | Sizes | Primary package/ Single peelpack | | Secondary package/ Shelf box | | Tertiary package/ Shipper | | | | Pallet |
|---|-------------|-------------------------------------|-------------------------|---------------------------------|---------------------------|------------------------------|---------------------------|-----------------|--------------------------|------------|
| | | Qty pcs | Dimensions, mm (LxW) | Qty pcs | Dimensions, mm (LxWxD) | Qty pcs | Dimensions, mm (LxWxD) | Gross Wt. kg | Volume m ³ | Qty pcs |
| Magill shaped Tracheal Tubes with HVLP cuff | | | | | | | | | | |
| UM611100XX UM611200XX | 3.0 to 4.5 | 1 | 271 x 90 | 10 | 280 x 125 x 90 | 100 | 470 x 295 x 270 | 5.0 | 0.037 | 1,800 |
| UM611100XX UM611200XX | 5.0 to 10.0 | 1 | 347 x 90 | 10 | 355 x 160 x 94 | 100 | 495 x 370 x 320 | 6.0 | 0.059 | 1,200 |
| Magill shaped Tracheal Tubes without HVLP cuff | | | | | | | | | | |
| UM611300XX UM611400XX | 2.0 to 4.5 | 1 | 271 x 90 | 10 | 280 x 125 x 90 | 100 | 470 x 295 x 270 | 5.0 | 0.037 | 1,800 |
| UM611300XX UM611400XX | 5.0 to 6.0 | 1 | 347 x 90 | 10 | 355 x 160 x 94 | 100 | 495 x 370 x 320 | 6.0 | 0.059 | 1,600 |
| Magill shaped Tracheal Tubes with preloaded stylet and HVLP cuff | | | | | | | | | | |
| UM611101XX | 5.0 to 8.5 | 1 | 347 x 90 | 10 | 355 x 160 x 94 | 100 | 495 x 370 x 320 | 7.0 | 0.059 | 1,200 |

* XX= size, e.g. 45 = 4.5 mm I.D.



Tracheal Tubes

Standard



EAN codes:

Magill shaped Tracheal Tubes with HVLP cuff and Murphy eye

| REF | Primary package | Secondary package | Tertiary package |
|------------|--------------------|--------------------|--------------------|
| UM61110030 | (01)30768455148173 | (01)50768455148177 | (01)70768455148171 |
| UM61110035 | (01)30768455148180 | (01)50768455148184 | (01)70768455148188 |
| UM61110040 | (01)30768455148197 | (01)50768455148191 | (01)70768455148195 |
| UM61110045 | (01)30768455148203 | (01)50768455148207 | (01)70768455148201 |
| UM61110050 | (01)30768455148210 | (01)50768455148214 | (01)70768455148218 |
| UM61110055 | (01)30768455148227 | (01)50768455148221 | (01)70768455148225 |
| UM61110060 | (01)30768455148234 | (01)50768455148238 | (01)70768455148232 |
| UM61110065 | (01)30768455148241 | (01)50768455148245 | (01)70768455148249 |
| UM61110070 | (01)30768455148258 | (01)50768455148252 | (01)70768455148256 |
| UM61110075 | (01)30768455148265 | (01)50768455148269 | (01)70768455148263 |
| UM61110080 | (01)30768455148272 | (01)50768455148276 | (01)70768455148270 |
| UM61110085 | (01)30768455148289 | (01)50768455148283 | (01)70768455148287 |
| UM61110090 | (01)30768455148296 | (01)50768455148290 | (01)70768455148294 |
| UM61110095 | (01)30768455148302 | (01)50768455148306 | (01)70768455148300 |
| UM61110100 | (01)30768455148319 | (01)50768455148313 | (01)70768455148317 |

Magill shaped Tracheal Tubes with HVLP cuff, without Murphy eye

| REF | Primary package | Secondary package | Tertiary package |
|------------|--------------------|--------------------|--------------------|
| UM61120030 | (01)30768455148326 | (01)50768455148320 | (01)70768455148324 |
| UM61120035 | (01)30768455148333 | (01)50768455148337 | (01)70768455148331 |
| UM61120040 | (01)30768455148340 | (01)50768455148344 | (01)70768455148348 |
| UM61120045 | (01)30768455148357 | (01)50768455148351 | (01)70768455148355 |
| UM61120050 | (01)30768455148364 | (01)50768455148368 | (01)70768455148362 |
| UM61120055 | (01)30768455148371 | (01)50768455148375 | (01)70768455148379 |
| UM61120060 | (01)30768455148388 | (01)50768455148382 | (01)70768455148386 |
| UM61120065 | (01)30768455148395 | (01)50768455148399 | (01)70768455148393 |
| UM61120070 | (01)30768455148401 | (01)50768455148405 | (01)70768455148409 |
| UM61120075 | (01)30768455148418 | (01)50768455148412 | (01)70768455148416 |
| UM61120080 | (01)30768455148425 | (01)50768455148429 | (01)70768455148423 |
| UM61120085 | (01)30768455148432 | (01)50768455148436 | (01)70768455148430 |
| UM61120090 | (01)30768455148449 | (01)50768455148443 | (01)70768455148447 |
| UM61120095 | (01)30768455148456 | (01)50768455148450 | (01)70768455148454 |
| UM61120100 | (01)30768455148463 | (01)50768455148467 | (01)70768455148461 |

Magill shaped Tracheal Tubes without HVLP cuff, with Murphy eye

| REF | Primary package | Secondary package | Tertiary package |
|------------|--------------------|--------------------|--------------------|
| UM61130020 | (01)30768455148500 | (01)50768455148504 | (01)70768455148508 |
| UM61130025 | (01)30768455148517 | (01)50768455148511 | (01)70768455148515 |
| UM61130030 | (01)30768455148524 | (01)50768455148528 | (01)70768455148522 |
| UM61130035 | (01)30768455148531 | (01)50768455148535 | (01)70768455148539 |
| UM61130040 | (01)30768455148548 | (01)50768455148542 | (01)70768455148546 |
| UM61130045 | (01)30768455148555 | (01)50768455148559 | (01)70768455148553 |
| UM61130050 | (01)30768455148562 | (01)50768455148566 | (01)70768455148560 |
| UM61130055 | (01)30768455148579 | (01)50768455148573 | (01)70768455148577 |
| UM61130060 | (01)30768455148586 | (01)50768455148580 | (01)70768455148584 |



Tracheal Tubes

Standard



EAN codes:

Magill shaped Tracheal Tubes without HVLP cuff and Murphy eye

| REF | Primary package | Secondary package | Tertiary package |
|------------|--------------------|--------------------|--------------------|
| UM61140020 | (01)30768455148593 | (01)50768455148597 | (01)70768455148591 |
| UM61140025 | (01)30768455148609 | (01)50768455148603 | (01)70768455148607 |
| UM61140030 | (01)30768455148616 | (01)50768455148610 | (01)70768455148614 |
| UM61140035 | (01)30768455148623 | (01)50768455148627 | (01)70768455148621 |
| UM61140040 | (01)30768455148630 | (01)50768455148634 | (01)70768455148638 |
| UM61140045 | (01)30768455148647 | (01)50768455148641 | (01)70768455148645 |
| UM61140050 | (01)30768455148654 | (01)50768455148658 | (01)70768455148652 |
| UM61140055 | (01)30768455148661 | (01)50768455148665 | (01)70768455148669 |
| UM61140060 | (01)30768455148678 | (01)50768455148672 | (01)70768455148676 |





Magill shaped Tracheal Tubes with preloaded stylet and HVLP cuff

| REF | Primary package | Secondary package | Tertiary package |
|------------|--------------------|--------------------|--------------------|
| UM61110150 | (01)30768455154853 | (01)50768455154857 | (01)70768455154851 |
| UM61110155 | (01)30768455154860 | (01)50768455154864 | (01)70768455154868 |
| UM61110160 | (01)30768455154877 | (01)50768455154871 | (01)70768455154875 |
| UM61110165 | (01)30768455154884 | (01)50768455154888 | (01)70768455154882 |
| UM61110170 | (01)30768455154891 | (01)50768455154895 | (01)70768455154899 |
| UM61110175 | (01)30768455154907 | (01)50768455154901 | (01)70768455154905 |
| UM61110180 | (01)30768455154914 | (01)50768455154918 | (01)70768455154912 |
| UM61110185 | (01)30768455154921 | (01)50768455154925 | (01)70768455154929 |

Packaging material:

| | |
|-------------------|---|
| Primary package | Peelback pouch - Uncoated paper, XPP film |
| Secondary package | Ivory cardboard |
| Tertiary package | Corrugated cardboard |

Storage and handling:

| | |
|---|---|
|  | FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care |
|  | KEEP DRY |
|  | THIS WAY UP - Indicates correct upright position of the transport package |
|  | DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package |


Well Lead Medical Co.,Ltd.
 C-4 Jinhu Industrial Estate
 Hualong, Panyu, Guangzhou
 511434 P.R.China

Made in China



Shanghai Intl. Holding Corp. GmbH
 Eiffestrasse 80, 20537 Hamburg
 Germany

Distributed by



ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2020 ConvaTec Inc. AP-019246-MM, ver. NOV 2020